

## Azlon® Low Form Beakers, Straight



- Straight stable shape, ideal for general laboratory use
- Complies with ISO 7056
- Good chemical resistance
- Autoclavable at 121°C
- Graduations are clearly printed in indelible solvent resistant ink




Code	Alt Ref	Capacity, mL	Graduations, mL	Dims, Ø x h, mm	Pack	Price
<b>BEA5022</b>	BDA212P	25	2	42 x 40	10	£19.70
<b>BEA5024</b>	BDA218P	50	5	51 x 51	10	£17.20
<b>BEA5026</b>	BDA224P	100	10	62 x 62	10	£19.70
<b>BEA5028</b>	BDA230P	250	25	82 x 82	10	£30.30
<b>BEA5030</b>	BDA234P	400	50	94 x 95	10	£35.20
<b>BEA5032</b>	BDA236P	500	50	93 x 101	10	£41.50
<b>BEA5034</b>	BDA238P	600	50	106 x 107	5	£22.60
<b>BEA5036</b>	BDA240P	1000	100	127 x 126	5	£28.40
<b>BEA5038</b>	BDA242P	2000	200	158 x 160	Each	£10.00

## Azlon® Low Form Beakers, Griffin Style



- Traditional, griffin shape
- Complies with ISO 7056
- Manufactured from PMP
- Can withstand temperatures of up to 150°C
- Suitable for food contact applications in accordance with U.S. FDA regulations
- Good chemical resistance
- Autoclavable
- Graduations are clearly printed in indelible solvent resistant ink



Code	Alt Ref	Capacity, mL	Graduations, mL	Dims, Ø x h, mm	Pack	Price
<b>BEA1310</b>	BDD312P	25	5	35 x 50	12	£36.30
<b>BEA1311</b>	BDD318P	50	10	42 x 60	12	£39.80
<b>BEA1312</b>	BDD324P	100	10	51 x 69	12	£55.00
<b>BEA1313</b>	BDD330P	250	25	71 x 94	6	£41.40
<b>BEA1314</b>	BDD334P	400	50	80 x 111	6	£38.20
<b>BEA1315</b>	BDD336P	500	50	87 x 117	6	£54.00
<b>BEA1317</b>	BDD340P	1000	100	107 x 148	6	£82.00