



GHS-CLP Pictogram Labels, Sheets

Self adhesive GHS-CLP hazard warning labels, each 20 x 20mm on a white background with a black pictogram, bordered by a red diamond. Available as single label sheets or a mixed sheet with more of the common labels.

Supplied as a sheet of 120 labels

Code	Description	Price
LAB4574	Explosive	£7.50
LAB4576	Flammable	£7.50
LAB4578	Oxidising	£7.50
LAB4580	Compressed gas	£7.50
LAB4582	Corrosive	£7.99
LAB4584	Acutely toxic	£7.50
LAB4586	Irritant/toxic	£8.90
LAB4588	Hazardous to health	£7.99
LAB4590	Dangerous to the environment	£7.50
LAB4592	Mixed GHS labels	£10.28



SLS Select Mercury Spillage Kit



- All you need to deal quickly and easily with mercury spillages
- Designed to cope with mercury spillages on any surface – flooring, carpet or bedding
- Can be used to deal with up to five spillages before requiring spares
- Long shelf life
- Housed in strong plastic case
- Uses COSHH recommended system for dealing with spillages
- Decontamination procedure helps eliminate effects of 'old' spillages
- Includes waste container, microfine sulphur, calcium hydroxide, pair of gloves, mask, scoop, syringe, brush, and alloy wool

Code	Capacity	Price
SAF7020	5 spillages	£79.00



Economy Chemical Spill Kit

- Portable kits ensure a fast response to small spill
- Easy to keep in small storage areas
- Suitable for use on basic drips and spills
- Absorbs aqueous chemicals and water
- Stored in a bright yellow, lightweight, water resistant bag to ensure sorbents stay dry
- Includes 10 x pads (400 x 500mm), 2 x socks (76mm Ø x 1220mm), pair of gloves, and disposal bag

Code	Absorption Capacity, L	Price
SAF7010	16	£143.00

Accessories

Code	Description	Dims, w x d, mm	Pack	Price
SAF7014	Spare chemical sorbent sheets	280 x 330	50	£340.00