


Heathrow Scientific Reversible Lab Mats



- Durable, chemical inert, FDA approved silicone material
- Noise dampening, anti skid, stain resistant and washable working surface
- Lip around the outer edge keeps benchtops safe and clean by containing spillages
- Designated area with internal metal inserts can be used to keep magnetic items such as stir bars organised
- Protects countertops from hot items up to 200°C
- Autoclavable at 121°C for 15 minutes
- Dims, w x l: 600 x 350mm

Code	Alt Ref	Side 1	Side 2	Price
SAF7090	HS120506	Yellow	Blue design on grey background	£65.00
SAF7092	HS120507	Purple	Grey design on purple background	£65.00
SAF7094	HS120599	Blue	White	£65.00



SLS Select Autoclave Bags, HDPE



- High density polyethylene bags for general disposal
- Autoclavable up to 121°C
- Sizes quoted are approximate

Code	Dims, w x l, mm	Pack	Price
AUT1400	670 x 305	200	£39.40
AUT1402	625 x 406	200	£49.50
AUT1404	840 x 635	100	£48.30



SLS Select Autoclave Bags, PP



- Transparent polypropylene disposal bags
- Autoclavable to 160°C
- Sizes quoted are approximate

Code	Dims, w x l, mm	Pack	Price
AUT1410	634 x 413	500	£138.00
AUT1415	766 x 608	300	£246.00



How prepared are you to deal with a spill?

Contact the team at SLS to arrange a **free site assessment!**

