



SLS Select Autoclave Bags, Paper



- Designed for safe, effective packing of instruments, dressings, and pipettes prior to steam sterilisation
- Manufactured from medical grade paper, adhesive and inks
- Designed with a stepped top for easy filling
- Available with a heat seal band and steam indicator or plain top without indicators
- Serrated top to prevent cuts
- Printed with Class 1 process indicators compliant with ISO 11140-1:2014
- Manufactured under ISO 9001:2015 and ISO 13485:2016 quality system
- Conforms to ISO 11607-1:2019 and BS EN 868-4:2017

Heat Seal Closure

Code	Dims, w x d x l, mm	Pack	Price
AUT1432	90 x 50 x 170	1000	£65.00
AUT1434	90 x 50 x 250	1000	£97.00
AUT1435	110 x 30 x 190	1000	£78.00
AUT1440	125 x 30 x 170	1000	£87.00
AUT1442	125 x 50 x 250	1000	£128.00
AUT1444	140 x 75 x 250	1000	£129.00
AUT1445	140 x 50 x 330	500	£77.00
AUT1448	180 x 95 x 380	500	£112.00
AUT1450	190 x 65 x 330	500	£97.00
AUT1451	250 x 100 x 380	250	£75.00
AUT1454	380 x 125 x 610	250	£150.00

Plain Seal Closure

Code	Dims, w x d x l, mm	Pack	Price
AUT1502	90 x 50 x 170	1000	£71.00
AUT1504	90 x 50 x 250	2000	£95.00
AUT1508	110 x 30 x 190	1000	£76.00
AUT1512	125 x 30 x 170	1000	£81.00
AUT1514	125 x 50 x 250	1000	£123.00
AUT1516	140 x 75 x 250	1000	£124.00
AUT1517	140 x 50 x 330	500	£76.00
AUT1520	180 x 95 x 380	500	£109.00
AUT1521	180 x 75 x 530	250	£79.00
AUT1523	250 x 100 x 380	250	£72.00
AUT1526	380 x 125 x 610	250	£160.00



SLS Select Holders for Autoclave Bags



Code	Description	For use with	Price
AUT14200	Small	300mm wide bags	£124.06
AUT14220	Large	686mm wide bags	£53.80



Waste Bins

Code	Description	Colour	Capacity, L	Price
BIN1002	Bin with swing open, lift off lid	Grey	50	£42.50
BIN1010	Dustbin with lift off lid	Black	95	£33.00