



Endecotts Octagon 200CL Sieve Shaker



Precise, reproducible and error free sieving processes. Several unique features have been developed specifically for this sieve shaker, including the 'Closed Loop' amplitude control for ultimate reproducibility. Designed to work with Endecotts' SieveWare, the new software for easy evaluation and documentation of the sieving process.

- 'Closed Loop' total amplitude control ensures reproducible sieving
- Digital controls for easy and reliable operation
- Easy to use sieve clamping system
- Suitable for dry and wet sieving
- 3D sieving motion allows for high separation efficiency and non blinding sieving action
- Full compatibility with SieveWare evaluation and control software via RS232 port (printed or digital protocols)
- Digital 00:10 – 99:50 minute timer
- No mechanical moving parts
- Compact and portable

Code	Alt Ref	Range, mm	Amplitude, mm	Suitable for Sieve Ø, mm	Max Sieves (full/half height)	Price
SIE3076	OCT200CL/100-240	0.02 - 125	0 - 3 in 0.1 steps	100, 200	8/16	£3053.00



Endecotts EFL 300 Sieve Shaker



Extremely robust but at the same time lightweight and powerful, making it ideal for sieving heavier and larger samples. The sieve shaker EFL 300 is very versatile. Amplitude setting is standard, and the heavy electric motor is replaced by the electromagnetic system found in all modern sieve shakers. Its lighter form means that it can be either floor standing or even bench mounted making it suitable for both, laboratory and industrial environments. Sieving parameters are set by the remote control unit. Its functions are logical and very simple to operate.

- Electromagnetic drive
- Adjustable amplitude
- Floor or table mounted
- Suitable for wet or dry sieving
- Digital controls for easy and reliable operation via external interface
- Economical
- Digital 00:10 – 99:50 minute timer

Code	Alt Ref	Range, mm	Amplitude, mm	Suitable for Sieve Ø, mm	Max Sieves (full/half height)	Price
SIE3078	EFL300/23050	0.02 - 40	0 - 2 in 10 steps	100, 150, 200, 250, 300, 315	6/12	£2655.00



Model SV001 Sieve Shaker

- Incorporates an electric motor to provide an effective sieving action
- Will accept 200, 250, 300 and 315mm, 8 and 12 inch diameter sieves
- Can hold up to 8 x 200mm diameter sieves or 7 x 300mm diameter sieves
- Fitted with a 0 to 60 minute timer

Code	Dims, w x d x h, mm	Weight, kg	Price
SIE1200	400 x 350 x 950	24	£1397.00