



Seward Stomacher® Bags

seward
World leaders in sample preparation


- Designed to fulfil all applications from sample blending, straining and storage, including bags designed for extensive blending without tears or perforations and small samples for life science applications
- High quality food grade virgin LDPE granule
- Extruded to 65µm film thickness (unless otherwise stated)
- No risk of cross contamination between samples
- Irradiated sterile



STANDARD



STRAINER



CLOSURE

Standard Bags

- Simplest and widely used bag, designed for blended homogenisation of diluent samples for most applications

Code	Alt Ref	For use with Stomacher®	Capacity, mL	Inner Pack	Pack	Price
MIX3102	BA6041	400 Classic	80 - 400	50	500	£51.00
MIX3175	BA6041/5/500	400 Classic	80 - 400	5	500	£80.00
MIX3150	BA6141	400 Circulator	80 - 400	50	500	£40.10
MIX3103	BA6141/ARV	400 Circulator	80 - 400	50	1000	£54.00
MIX3177	BA6141/5/500	400 Circulator	80 - 400	5	500	£81.00
MIX3104	BA6042	3500 Series	400 - 3500	50	250	£73.00
MIX3156	BA6040	80 Biomaster	5 - 80	50	2000	£70.00

Strainer Bags

- For samples that produce a large amount of debris or fat or harder samples such as shellfish
- Contain a robust strainer element to make it easier to pipette the sample produced (pore size 0.5mm, prevents pore blockage and non-homogenous sample)
- MIX3164 is a double bag for biosafety

Code	Alt Ref	For use with Stomacher®	Capacity, mL	Pack	Price
MIX3108	BA6041/STR	400 Classic	80 - 400	200	£73.00
MIX3109	BA6141/STR	400 Circulator	80 - 400	200	£89.00
MIX3160	BA6040/STR	80 Biomaster	5 - 80	200	£41.00
MIX3164	BA6040/DBL/STR	80 Biomaster	5 - 80	200	£55.00

Closure Bags

- A robust range of bags with a wire mesh band at the top of the bag, allowing easy opening and tight closure for storage in a bag rack or incubator
- MIX3162 is a closure bag which also has a robust strainer element to make it easier to pipette the sample produced

Code	Alt Ref	For use with Stomacher®	Capacity, mL	Pack	Price
BAG2020	BAG6042/CLR	3500 Series	400 - 3500	200	£149.00
MIX3110	BA6041/CLR	400 Classic	80 - 400	250	£80.00
MIX3111	BA6141/CLR	400 Circulator	80 - 400	250	£58.00
MIX3158	BA6040/CLR	80 Biomaster	5 - 80	200	£43.00
MIX3162	BA6040/CLR/STR	80 Biomaster	5 - 80	200	£55.00