

Hardwood Stands

Single and double hardwood stands for burettes and conical funnels. Plastic clamping screws ensure easy movement when loose and sturdy placement when tightened. Rods are approximately 460mm high. Single base approximately 200 x 110mm, double base 290 x 110mm.

Code	Description	For use with	Price
STA1000	Single	Burettes	£32.80
STA1010	Single	Conical funnels	£51.00
STA1002	Double	Burettes	£35.60
STA1012	Double	Conical funnels	£42.00



Boss Head Clamps

For use with rods up to 16mm diameter.


Code	Price
STA1200	£8.00



Vessel Clamps

A range of natural finish alloy and stainless steel clamps suitable for various applications, as specified.

Code	Description	Price
STA1300	Alloy, cork lined jaws	£18.79
STA1302	Alloy, 3 prong, rubber covered	£18.90
STA1304	Steel, 3 prong, rubber covered	£62.00
STA1306	Alloy, 3 prong, heat resistant cover	£66.00
STA1308	Alloy, 4 prong, heat resistant cover	£69.00
STA1310	Alloy, 2 prong, extendable for reaction vessels	£71.00



Kartell Burette Clamps

Kartell
LABWARE

Polypropylene, spring-loaded burette clamps for one or two burettes.

Code	Alt Ref	Description	Price
STA1400	139	Single clamp	£28.00
STA1402	140	Double clamp	£38.00



Kartell Separating Funnel Holder

Kartell
LABWARE

Polypropylene funnel holder with cutout, for holding separating funnels of diameter 125 to 1000mm. Cannot be used with funnels below 75mm diameter. Supplied with two olives for use with rods from 8 to 14mm.

Code	Alt Ref	Price
STA1094	970	£30.00

