



## SLS Lab Basics Mini Magnetic Stirrer



An economical solution made for ultimate function with optimal magnetic coupling to ensure dependable stirring efficiency. A powerful magnetic coupling of the high-quality stir bar and the unit ensures consistent accurate alignment which significantly reduces spin outs during acceleration or deceleration of spin speeds and ensure maximum stirring performance. Ideal for routine stirring applications in general research, life science, industrial, quality analysis, and clinical lab settings.

- Stirs up to 3L
- Quiet operation with low vibration, noise level, full speed, unloaded: 53dB
- Low profile design with a small footprint to save valuable space on a bench top
- Lightweight for transporting with ease
- Made from an ABS and PC blend that is chemical resistant and easy to clean
- Supplied with 2 stir bars (9 x 35mm) and universal power adapter

Code	No. of Places	Speed Range, RPM	Dims, w x d x h, mm	Weight, kg	Price
SLS9604	1	200 - 2200	132 x 183 x 47	0.45	£137.00




## Cole-Parmer® Analog Magnetic Stirrer

Cole-Parmer

A large, heavy duty, powerful stirrer for use with vessels up to 30L capacity. Designed with safety as well as performance in mind. The compact shape takes up less bench space and makes storage easier.

- Robust construction
- Stainless steel top plate
- Powerful magnets for strong magnetic coupling
- Stirs up to 30L
- Supplied with 2 x 25mm PTFE coated stir bars

Code	Alt Ref	Model	No. of Places	Speed Range, RPM	Plate Dims w x l, mm	Dims, w x d x h, mm	Weight, kg	Price
STI2015	04801-11	SB301	1	100 - 600	300 x 300	328 x 369 x 103	5	£637.00



## Cole-Parmer® Multiposition Magnetic Stirrer

Cole-Parmer

A space saving and economical unit which only requires one power source – ideal for quality control applications where multiple samples require processing under the same conditions.

- Powerful magnets for strong magnetic coupling and stirring capacity
- Easily accommodate a 2L beaker in each position
- Top plate material: aluminium
- Supplied with 3 x 25mm PTFE coated stir bars

Code	Alt Ref	Model	No. of Places	Speed Range, RPM	Plate Dims w x l, mm	Dims, w x d x h, mm	Weight, kg	Price
STI2311	51400-02	SB161-3	3	0 - 1500	210 x 550	270 x 600 x 110	7	£747.00