

Cole-Parmer® Stuart™ Undergrad Analog Stirrers

Cole-Parmer
stuart™

Constructed with powerful magnets and motor, with a square plate design allows you flexibility of vessel sizes and combinations. Up to four 250mL beakers can easily fit on the surface at the same time. For safety, all units are designed to direct spills away from the user controls. The stirrers also feature BioCote® antimicrobial protection which reduces the presence of odour and stain causing microbes such as bacteria, mold, and fungi on the product surface by up to 99.99%. Formulated from silver ion additives, it provides built in antimicrobial protection to support a more hygienic laboratory environment.

- Available with either a chemically resistant glass ceramic plate for added chemical resistance or stainless steel plate for stronger magnet coupling
- Compact, stackable design takes up less benchtop and storage space
- Built with tough aluminum casing for longevity
- Designed to prevent accidental spillage on control panel
- IP32 rated for water and dust resistance



Specifications	
No. of Places	1
Speed Range, RPM	100 – 2000
Max Stirring Capacity, L	15
Plate Dimensions, mm	150 x 150
Dims, w x d x h, mm	180 x 255 x 109
Weight, kg	2.4

Code	Alt Ref	Model	Plate Material	Price
ST12610	04807-54	UC151	Glass ceramic	£257.00
ST12612	04807-56	US151	Stainless steel	£253.00

Accessories

Code	Alt Ref	Description	Dims, h x Ø, mm	Price
ST11062	04805-91	Retort rod	600 x 12	£42.40

Hanna Mini Magnetic Stirrers

HANNA
instruments

Compact and lightweight stirrers making lack of laboratory bench space no longer a concern. These stirrers incorporate electronic controls that allow the user to regulate the speed with greater precision.

- Internal Speedsafe™ mechanism ensures that the maximum speed is never exceeded
- Compact size maximises bench space
- Supplied with micro stirrer bar

Specifications	
No. of Places	1
Speed Range, RPM	100 - 1000
Max Stirring Capacity, L	1
Dims, w x d x h, mm	120 x 120 x 45
Weight, g	610

Code	Alt Ref	Supplied with	Weight, g	Price
ST12500	HI-190M	ABS cover	610	£106.00
ST12510	HI-200M	Stainless steel cover	640	£195.00

