



IKA Topolino Magnetic Mini Stirrers

IKA®

Extremely convenient magnetic mini-stirrers for mixing quantities up to 250mL. Features a durable, brushless motor, polypropylene plate and a continuously adjustable speed range. High magnetic adhesion.

- Durable, brushless motor
- Continuously adjustable speed range
- High magnetic adhesion

Topolino Mobil has the same features as the Topolino, plus:

- Portable unit with long 8 to 12 hour operating time
- Short 2 to 3 hour charging time
- Standard replaceable AA rechargeable batteries
- Optional power modes: a) mains-free with standard batteries, b) with supplied mains adapter (without batteries), c) combined mains/battery operation (with batteries fitted)

Specifications	
No. of Places	1
Speed Range, RPM	500 - 1800
Max Stirring Capacity, mL	250
Max Magnetic Bar Length, mm	30
Plate Ø, mm	80

Code	Alt Ref	Model	Dims, w x d x h, mm	Weight, g	Price
STI2452	0003368000	Topolino	95 x 115 x 37	484	£129.00
STI2450	0003381300	Topolino Mobil	140 x 140 x 42	747	£275.00



IKA Squid Magnetic Stirrers

IKA®

Small magnetic stirrers with digital LED speed display. Outstanding chemical resistance due to glass plate and synthetic bottom made of TPC-ET.

- Electronically controlled motor for more capacity
- Recyclable materials

Specifications	
No. of Places	1
Speed Range, RPM	0 - 2500
Max Magnetic Bar Length, mm	30
Setting Accuracy Speed, RPM	50

Code	Alt Ref	Description	Stirring Capacity, L	Plate Ø, mm	Dims, w x d x h, mm	Weight, g	Price
STI2273	0003671000	Colour squid	1	115	145 x 160 x 45	1074	£360.00
STI2275	0003672000	Big squid	1.5	160	180 x 195 x 40	889	£372.00