



Benchmark Digital Magnetic Stirrer



Provide reliable performance and safe operation for everyday laboratory use. The space efficient stirrer series features a durable, chemical resistant white ceramic top plate and advanced microprocessor controls. In addition, the powerful motor drive system ensures efficient mixing in larger and more viscous samples.

- Mixing speed can be viewed on the large digital display
- Large, ceramic work surface
- 5L stir capacity
- Compatible with borosilicate glass beakers and other vessels

Code	Alt Ref	No. of Places	Speed Range, RPM	Platform Dims, mm	Dims, w x d x h, mm	Weight, kg	Price
MIX1263	H3770-S-E	1	150 - 1500	178 x 178	180 x 279 x 101	2.9	£291.00



SLS Lab Basics Hotplate Stirrer, 280°C



Stainless steel stirrer with a ceramic coated hot plate which provides good chemical resistance. Features a digital control system, ideal for use in physics, chemistry and biotech labs.

- Heating temperature up to 280°C
- LED display shows temperature and speed
- 50°C HOT warning indicates residual hotplate temperature
- ABS chemical resistant casing
- Protection class: IP21

Code	Stirring Capacity, L	Speed Range, RPM	Max Temp, °C	Plate Ø, mm	Price
SL54780	3	100 - 1500	280	135	£365.00



SLS Lab Basics Hotplate Stirrer, 340°C



Stirrer hotplate with ceramic coated stainless steel plate, widely used in chemical synthesis, physical and chemical analysis, bio-tech labs. Featuring PID temperature technology precisely controls the heating process, rapidly reaches target temperature and enhances control accuracy. Automatically identifies temperature sensor PT1000 and transfers to PT1000 operation mode.

- Simple running interface easy to use and program
- Records and displays the last running parameters when switched on
- LCD display shows set and actual temperatures, and speed
- Brushless DC motor - maintenance free and explosion-proof
- Fully closed assembly guarantees safety and long service life, effectively prevents the liquid damage
- HOT warning indicates residual hotplate temperature
- Free software StirPC which enables PC control via RS232 interface
- Protection class: IP42

Code	Stirring Capacity, L	Speed Range, RPM	Max Temp, °C	Plate Ø, mm	Price
SLS4800	20	100 - 1500	340	135	£638.00