

Cole-Parmer® Stuart™ Undergrad Hotplate Stirrers Continued

Cole-Parmer
stuart™

Model	US152	UC152
Plate Material	Coated aluminium/silicon	Glass ceramic
Plate Dimensions, mm	150 x 150	
Heated Area, mm	150 x 150	120 x 120
Heater Power, W	700	500
Max Plate Temp, °C	325	450
Speed Range, RPM	100 - 2000	
Maximum Stirring Capacity, L	15	
Dims, w x d x h, mm	172 x 248 x 120	172 x 248 x 122
Weight, kg	2.9	

Code	Alt Ref	Model	Supplied with	Price
STI2602	04807-60	US152 analogue	2 x 25mm PTFE coated stir bars	£353.00
STI2615	US152/KIT	US152 analogue	SCT1 controller, probe holder, retort rod, 2 x 25mm PTFE coated stir bars	£576.00
STI2604	04807-64	US152 digital	2 x 25mm PTFE coated stir bars	£419.00
STI2617	US152D/KIT	US152 digital	SCT1 controller, probe holder, retort rod, 2 x 25mm PTFE coated stir bars	£619.00
STI2600	04807-58	UC152 analogue	2 x 25mm PTFE coated stir bars	£367.00
STI2614	UC152/KIT	UC152 analogue	SCT1 controller, probe holder, retort rod, 2 x 25mm PTFE coated stir bars	£587.00
STI2606	04807-62	UC152 digital	2 x 25mm PTFE coated stir bars	£419.00

Accessories

Code	Alt Ref	Description	Dims, h x Ø, mm	Price
STI5500	SCT1	Digital contact thermometer	—	£220.00
STI1062	04805-91	Retort rod	600 x 12	£42.40



Cole-Parmer® Analog Hotplate Stirrer

Cole-Parmer

General purpose stirring hotplate designed with safety as well as performance in mind. The compact shape takes up less bench space and makes storage easier. Powerful magnets and motor give stirring speeds up to 1500RPM and is capable of mixing large volumes (up to 15L). The chemically resistant glass ceramic top plate stays cool around the edges to help prevent burns.

- “Hot” warning light will flash whenever the plate temperature is above 50°C even when the hotplate is turned off (lasts up to 30 minutes without power)
- Compact, space saving design suits even cluttered laboratories
- Powerful magnets provide strong stirring action
- Integral fitting for attaching a support rod
- Supplied with 2 x 25mm PTFE coated stir bars

Code	Alt Ref	Model	Plate Dims, mm	Heated Area, mm	Heater Power, W	Max Temp, °C	Speed Range, RPM	Dims, w x d x h, mm	Price
STI2014	04800-01	CB302	300 x 300	200 x 200	1200	450	100 - 1500	300 x 365 x 105	£752.00