

## IKA EUROSTAR Control Overhead Stirrers

**IKA®**

Extremely powerful laboratory stirrer designed with a removable wireless controller and a digital TFT display. It automatically adjusts the speed through microprocessor controlled technology. The stirrer comes equipped with a RS 232 and a USB interface to control and document all parameters. An integrated torque trend display is provided for the measurement of viscosity changes. Safety circuits installed ensures automatic cut off in an anti-stall or overload conditions. Continuous comparison of shaft speed to desired speed is maintained and variations are adjusted automatically. This guarantees a constant speed even with changes in viscosities of the sample.

- Multilingual TFT display
- Programmable functions
- Integrated temperature measurement
- Interval operation
- Timer function
- Adjustable safety circuit
- Locked function
- Infinitely adjustable speed
- Push through agitator shafts
- Overload protection
- Short term overload operation
- Slim casing
- Quiet operation
- Error code display
- Supplied with H 67.60 temperature sensor, WH11, OS 1.0 WiCo power supply unit, 2 x USB cables, screw driver, protective cover for WiCo contact and chuck key



Code	Alt Ref	Description	Capacity, L	Max Viscosity, mPas	Speed Range, RPM	Max Torque at Stirrer Shaft, Ncm	Dims, w x d x h, mm	Price
<b>STI2876</b>	0003992002	EUROSTAR 200	100	100000	6 - 400, 30 - 2000	200	91 x 231 x 297	<b>£3694.00</b>
<b>STI2878</b>	0004028502	EUROSTAR 100	100	70000	30 - 1300	100	86 x 230 x 267	<b>£2819.00</b>
<b>STI2886</b>	0004440002	EUROSTAR 60	40	50000	30 - 2000	60	86 x 230 x 267	<b>£2215.00</b>

## IKA RW Digital Overhead Stirrers

**IKA®**

Overhead stirrers with digital LED display, and a robust, slimline, ergonomic design.

- Constant power drive
- Two speed ranges for universal use
- Push through agitator shafts (only when stationary)
- For stirring quantities of up to 200L (H<sub>2</sub>O)
- Supplied with chuck key



Code	Alt Ref	Capacity, L	Max Viscosity, mPas	Speed Range, RPM	Max Torque at Stirrer Shaft, Ncm	Dims, w x d x h, mm	Price
<b>STI2836</b>	0003593002	20	10000	60 - 2000	150	88 x 212 x 294	<b>£1044.00</b>
<b>STI2896</b>	0005040000	80	50000	60 - 1400	900	123 x 252 x 364	<b>£3962.00</b>
<b>STI2898</b>	0004050000	200	100000	57 - 1300	3000	145 x 358 x 478	<b>£7143.00</b>