


IKA Overhead Stirrer Accessories



General Accessories

Code	Alt Ref	Description	Price
STI6328	0002336000	FK 1 flexible coupling, clamping range 6 - 10mm	£324.00
STI6334	0002603000	R 301 stirring shaft protection, length adjustment 190 - 310mm	£192.00
STI2822	0002604000	R 301.1 support holder, length: 310 mm	£228.00



4 Bladed Propeller Stirrers


Code	Alt Ref	Stirrer Ø, mm	Shaft Dims, Ø x l, mm	Max Speed, RPM	Price
STI6324	0000741000	50	8 x 350	2000	£160.00
STI2808	0000741300	100	8 x 350	800	£197.00
STI2802	0000739000	150	13 x 800	600	£689.00



3 Bladed Propeller Stirrers


Code	Alt Ref	Stirrer Ø, mm	Shaft Dims, Ø x l, mm	Max Speed, RPM	Price
STI6382	0001296000	45	8 x 350	2000	£320.00
STI6390	0002343600	45	8 x 350	2000	£579.00
STI6384	0001295900	55	8 x 350	2000	£333.00
STI6386	0000477700	140	10 x 550	800	£320.00
STI6388	0000477800	140	10 x 800	400	£382.00
STI2810	0001242900	55*	—	—	£275.00
STI2814	0001289800	45*	—	—	£221.00

* Propeller only, no shaft




Turbine Stirrers

Code	Alt Ref	Stirrer Ø, mm	Shaft Dims, Ø x l, mm	Max Speed, RPM	Price
STI6326	0002332900	30	8 x 350	2000	£216.00
STI6330	0002333000	50	8 x 350	2000	£274.00
STI6392	0002333100	70	10 x 400	800	£336.00



Dissolver Stirrers

Code	Alt Ref	Stirrer Ø, mm	Shaft Dims, Ø x l, mm	Max Speed, RPM	Price
STI6325	0000513500	80	8 x 350	2000	£283.00
STI6396	0002746700	100	10 x 350	1000	£443.00
STI6394	0002387900	40	8 x 350	2000	£599.00



Centrifugal Stirrers

Code	Alt Ref	Stirrer Ø, mm	Shaft Dims, Ø x l, mm	Max Speed, RPM	Price
STI6398	0000756900	60/15	8 x 350	2000	£178.00
STI6400	0001132700	100/24	8 x 550	800	£239.00