

PTFE Covered Rare Earth Stirrer Magnets Continued

Oval

- For use in round bottomed flasks
- The magnet core is also oval for maximum stirring effect

Code	Dims, Ø x h, mm	Price
STI4240	5 x 10	£0.91
STI4242	6 x 15	£0.84
STI4244	10 x 20	£1.29
STI4246	10 x 25	£1.50
STI4248	12 x 25	£1.50
STI4250	10 x 30	£1.96
STI4252	16 x 30	£2.41
STI4254	13 x 35	£2.64
STI4256	16 x 35	£2.71
STI4258	13 x 40	£2.85
STI4260	20 x 40	£3.59
STI4262	17 x 50	£4.41
STI4264	20 x 50	£5.00
STI4266	20 x 64	£6.14
STI4268	20 x 70	£8.12
STI4270	25 x 70	£9.79



Elliptical

Code	Dims, w x l, mm	Price
STI4610	6 x 10	£1.97
STI4612	10 x 15	£3.07
STI4614	14 x 25	£4.07
STI4616	24 x 50	£11.27
STI4618	28 x 70	£20.65



Triangular

- Particularly effective for dissolving solids and mixing sediments because of the scraper-like action on the base of the container

Code	Dims, Ø x h, mm	Price
STI5100	6 x 12	£1.19
STI5102	8 x 20	£1.45
STI5104	8 x 25	£1.29
STI5105	14 x 25	£2.23
STI5106	9 x 35	£2.01
STI5108	14 x 40	£2.41
STI5110	12 x 50	£2.75
STI5112	14 x 55	£3.29
STI5114	18 x 80	£4.38

