

PTFE Covered Rare Earth Stirrer Magnets Continued



Cross

- A very stable general purpose stirrer


Code	Ø, mm	Price
STI4110	10	£1.28
STI4112	20	£1.48
STI4114	25	£1.91
STI4116	30	£2.33
STI4118	38	£2.89



Glass Covered Stirrer Magnets

- Ideal for use with very abrasive media which may erode PTFE
- Alnico V permanent magnet core

Code	Dims, Ø x h, mm	Price
STI4130	5 x 6	£3.72
STI4132	5 x 12	£3.85
STI4134	6 x 25	£4.11
STI4136	8 x 45	£4.47
STI4138	8 x 60	£4.83



A

Magnet Retrievers

- Retrievers for the removal or recovery of magnetic stirring bars from containers or liquids
- Inert and leak free
- Moulded in pure, FDA approved PTFE or polypropylene
- Encapsulated Alnico V magnet
- Autoclavable

Code	Material	Dims, Ø x h, mm	Price
STI5200	PTFE	10 x 150	£8.48
STI5204	PTFE	10 x 250	£11.00
STI5206	PTFE	10 x 350	£12.00
STI5208	Polypropylene	10 x 300	£4.95
STI5210	Polypropylene	10 x 350	£5.24
STI5212	Polypropylene	10 x 450	£5.59



A

Thermo Scientific™ Nalgene™ Stirring Rods

Thermo
SCIENTIFIC

- Flat paddle at one end and triangular paddle at the other
- Polypropylene
- Chemically compatible with a wide range of lab reagents
- Usable over a temperature range of 0 to 135°C
- Autoclavable

Code	Alt Ref	Length, mm	Pack	Price
ROD1220	6169-0010	248	12	£47.90