



Turnover End Rubber Stoppers

- Natural red rubber
- With a turnover flange and external serrations to grip the inside of the container

Supplied in packs of 10

Code	Ø, mm	Price
ST06300	6.5	£8.66
ST06302	8	£7.61
ST06304	9.5	£10.10
ST06306	11	£9.97
ST06308	12.5	£8.66
ST06310	14	£9.07
ST06312	16	£10.76
ST06314	17.5	£12.85
ST06316	19	£17.00
ST06318	20.5	£19.39
ST06320	22	£24.03
ST06322	24	£24.90
ST06324	25.5	£25.57



A

Bacti-Caps

- Made from chemically resistant plastic
- Provide a firm, positive culture tube closure whilst allowing a controlled gas interchange
- An excellent and more practical alternative to cotton wool plugs
- Three flexible fins moulded on the inner wall of each cap ensure a positive fit even when minor variations in tube sizes occur
- Chemically resistant plastic
- Exclusive use of virgin material eliminates any possibility of leachable components
- Autoclavable

Supplied in packs of 100

Code	Size, mm	Colour	Price
CAP1000	13	Red	£42.00
CAP1002	13	White	£42.00
CAP1004	13	Blue	£42.00
CAP1006	13	Yellow	£42.00
CAP1008	13	Green	£42.00
CAP1010	13	Black	£42.00
CAP1012	16	Red	£44.00
CAP1014	16	White	£44.00
CAP1016	16	Blue	£44.00
CAP1018	16	Yellow	£44.00
CAP1020	16	Green	£44.00
CAP1022	16	Black	£44.00

5
YEARS

Gratnells F1 Shallow Trays



- Dimensions, w x d x h: 312 x 427 x 75mm

Code	Alt Ref	Colour	Price
SED11020	F0120	Translucent	£8.42
SED11022	F0102	Sunshine yellow	£8.42
SED11024	F0101	Tropical orange	£8.42
SED11026	F0109	Flame red	£8.42
SED11028	F0107	Fuschia pink	£8.42
SED11030	F0106	Royal blue	£8.42
SED11032	F0110	Grass green	£8.42