

3M™ Clean-Trace™ Hygiene Monitoring and Management System Continued

3M Science.
Applied to Life.™

ATP Test Swabs

Supplied in packs of 100

Code	Alt Ref	Description	Price
FSA1094	UXL100	Surface ATP test	£335.00
FSA1090	AQT200	Water plus - total ATP test	£426.00
FSA1092	AQF100	Water-free ATP test	£478.00



Surface Protein Test Swabs

Code	Alt Ref	Description	Pack	Price
FSA1108	PRO50	Surface protein plus test swab	50	£155.00
PRO100	PRO100	Surface protein plus test swab	100	£221.00
FSA1110	CIT50	Surface protein instant	50	£73.00
FSA1104	ALLTEC60	Surface protein (allergen) test swab	60	£200.00



3M™ Sponges and Swabs

3M Science.
Applied to Life.™

Sample collection and preparation solutions provide quality products for every step of your process. Convenient, ready-to-use products with consistent, reliable results.

- Easy to use for product sampling and environmental testing
- Exceptional convenience and quality
- Increased productivity and efficiency
- Reliable and accurate - time after time
- Improved consistency from technician to technician and plant to plant

Sponge Sticks

- 1.5 x 3 inch cellulose sponges mounted at one end of a stick and pre-hydrated
- Allows sample collection without touching the inside of the bag
- Stick helps access hard to reach areas and to collect the sample without directly handling the sponge
- Easily snap off the sponge from the handle after swabbing
- Pre-moistened sponge is durable and biocide-free and does not affect organism viability
- Cellulose material is absorbent and will effectively hold and release the sample without excess dripping
- 6 x 11" sample bag with a 30oz capacity for safe, convenient transportation

Supplied in packs of 100



Code	Alt Ref	Description	Price
FSA1166	SSL100	Without media	£194.00
FSA1180	SSL10BPW	With 10mL buffered peptone water broth	£184.00
FSA1168	SSL10NB	With 10mL neutralising buffer	£189.00