



Thermo Scientific™ Nunc™ EasyFill™ Cell Factory™ Systems

Thermo
SCIENTIFIC

- Each 10 layer, single use system holds the equivalent of 36 T175 flasks
- Linear multi layer cell culture format is easy to use and promotes easy scalability, ideal for both research and commercial scale applications
- "Out of the box" system
- Large, vented screw closure
- Versatile port design for pouring and aseptic filling
- Fill and empty five times faster than T175 flasks
- No accessories required

Code	Tray Layers	Culture Area, cm ²	Working Volume, mL	Pack	Price
140000	1	632	200	6	£290.00
140250	2	1264	400	6	£370.00
140360	4	2528	800	4	£408.00
140400	10	6320	2000	2	£361.00



Thermo Scientific™ Nunc™ Immuno MicroWell™ Plates

Thermo
SCIENTIFIC

- Plates intended for quantitative and qualitative solid phase immunoassays (e.g. ELISA) and binding assays
- Choice of flat, round bottom and flat with curved edges for optimal washing and liquid retrieval
- Raised rims on the wells
- Certified reproducibility of binding
- Alphanumeric orientation system
- Without lid
- Non sterile

Code	Surface Treatment	Colour	Description	Well Volume, µL	Bottom	Inner Pack	Pack	Price
475094	PolySorp	Clear	Pinchbar	400	Flat	5	60	£172.00
456529	PolySorp	Clear	High flange	400	Flat	10	180	£472.00
446140	PolySorp	Clear	Pinchbar	350	C	5	60	£162.00
475434	PolySorp	Clear	High flange	300	U	5	60	£165.00
467320	MediSorp	Clear	High flange	400	Flat	5	60	£237.00
467340	Multisorp	Clear	High flange	400	Flat	5	60	£237.00
456537	MaxiSorp	Clear	High flange	400	Flat	10	180	£615.00
439454	MaxiSorp	Clear	Pinchbar	400	Flat	5	60	£229.00
442404	MaxiSorp	Clear	Pinchbar	400	Flat	5	60	£168.00
430341	MaxiSorp	Clear	Pinchbar	350	C	5	60	£162.00
446612	MaxiSorp	Clear	Pinchbar	350	C	5	60	£236.00
449824	MaxiSorp	Clear	High flange	300	U	5	60	£158.00
436110	MaxiSorp	White	Pinchbar	400	Flat	10	80	£373.00
437111	MaxiSorp	Black	Pinchbar	400	Flat	10	80	£371.00