

GETINGE Lancer ULTIMA Series Labware Washing Systems

GETINGE ✱



GETINGE Lancer ULTIMA series washers meet space and budget requirements in laboratories performing a range of critical applications where clean labware is essential.

The exclusive fully programmable Prolux microprocessor controller commands a full range of prewash, wash, rinse and drying functions through simple touchscreen menus on all ULTIMA series washers. The convenience of five preset programs for light to heavy loads simplifies operation, and now comes with the inclusion of a new pre-programmed "eco-cycle" for those loads which need cleaning with the least demand on valuable resources. Full details and specifications for individual models can be found below.

As standard with the GETINGE Lancer range, all units are delivered to point of use, installed by qualified OEM engineers, commissioned to GETINGE Lancer standards, and supplied with FREE chemical starter pack, FREE first year service contract and 12 months warranty as standard.

ULTIMA Undercounter Washers

The ULTIMA series undercounter washers offer 3 high-performance models that offer a range of efficient labware loading options using baskets and injector racks. The compact footprint maximises valuable research space.

	810LX	815LX	820LX
Dims, w x d x h, mm	609 x 739 x 850	609 x 739 x 850	909 x 739 x 850
Effective Chamber Volume, L	140		
Wash Temperature, °C	Up to 95		
Drying	No	Yes - convection	Yes - HEPA filtered
View in Process (VIP)	No	No	Yes
Touchscreen Controller	3.5" colour stores up to 40 programs		
Low Pressure Valve & Pump Kit	Yes - for trouble free DI water connection		
Steam Condenser	Yes	Yes	Yes
Water Softener	No	No	No
Chemical Storage	No	No	No
Chemical Dosing Pumps	2 pumps to dispense detergent & acid rinse		
Chemical Levels Sensors	Option	Option	2
Number of Wash Levels	2	2	2
Number of Runner Positions	2	2	2
Noise Level, dB(A)	<63		



Code	Alt Ref	Model	Description	Price
WAS2010	1260365	810LX	No drying	£8320.00
WAS2015	1260366	815LX	Forced-air drying	£9093.00
WAS2020	1260367	820LX	Forced-air drying, injectors, HEPA filtered with on board chemical storage and additive level sensor	£13862.00