

## Thermo Scientific™ X4/X4R Multifuge Centrifuges Continued

**Thermo**  
SCIENTIFIC

## Rotors and Accessories

Code	Description	Capacity	Max Speed, RPM	Max RCF, xg	Pack	Price
75003017	TX-1000 rotor (without buckets)	—	3800*	3374**	Each	£1679.00
75003001	Buckets for TX-1000 rotor	—	—	—	4	£1679.00
75007309	ClickSeal biocontainment lids for TX-1000 rotor buckets	—	—	—	4	£513.00
75007303	Adapters for TX-1000 rotor	4 x 6 microplate or 4 x 2 deep well plates	—	—	4	£482.00
75007306	Adapters for TX-1000 rotor	24 x 15mL conical tubes	—	—	4	£482.00
75003671	Adapters for TX-1000 rotor	49 x 5/7mL Vacutainers™	—	—	4	£483.00
75003674	Adapters for TX-1000 rotor	10 x 50mL tubes	—	—	4	£612.00
75003180	TX-750 rotor (without buckets)	—	4700	4816	Each	£1521.00
75003608	750mL round buckets for TX-750 rotor	—	—	—	4	£1358.00
75003609	ClickSeal biocontainment lids for TX-750 rotor buckets	—	—	—	4	£349.00
75003795	Microplate carriers and T-75 flask for TX-750 rotor	4 x microplates or 1 x deepwell	—	—	2	£943.00
75003617	Microplate carriers and T-75 flask for TX-750 rotor	4 x microplates or 1 x deepwell	—	—	4	£1838.00
75003638	Adapters for TX-750 rotor	7 x 50mL conical tubes	—	—	4	£616.00
75003723	Adapters for TX-750 rotor	27 x 5/7mL Vacutainers™	—	—	4	£616.00
75003639	Adapters for TX-750 rotor	14 x 15mL conical tubes	—	—	4	£616.00
75003661	Fiberlite™ F13-14x50cy fixed angle rotor	14 x 50mL	10000	17105	Each	£4318.00
75003662	Fiberlite™ F14-6x250LE fixed angle rotor	6 x 250mL	11000	18533	Each	£4288.00
75003663	Fiberlite™ F15-8x50cy fixed angle rotor	8 x 50mL	14500	24446	Each	£2755.00

\* Ventilated centrifuge only, 4200RPM in refrigerated centrifuge

\*\* Ventilated centrifuge only, 4122xg in refrigerated centrifuge

## OHAUS Frontier™ 5706 Multi Centrifuge



A small benchtop multi function centrifuge ideal for a multitude of spin down applications. Characterised by reliability, various safety features, the ability to accommodate test tubes of various sizes and intuitive operation, Frontier 5706 was built to meet the basic need for liquid separation in the life sciences.

- Simple and intuitive operation
- Easy to turn control knob
- Large LCD display and the capabilities to set and view both speed in RPM and g-force
- Precise results and enhanced functionality with extremely simple operation
- Can accommodate all standard tube sizes from 1.5 to 50mL, providing centrifugation for a wide variety of applications
- Sturdy mechanical lid with locking system and a rotor imbalance sensor that stops the centrifuge if the rotors become imbalanced or sudden shock occurs
- Rotor not included and must be ordered separately



Code	Alt Ref	Model	Speed Range, RPM	Max RCF, xg	Dims, w x d x h, mm	Price
CEN2628	30130875	FC5706	200 - 6000	4427	291 x 392 x 266	£1285.00