

Kinesis Vials and Caps Continued



Screw Neck Vials

- 2mL screw neck vials
- Dims, h x Ø: 32 x 11.6mm
- Neck Ø: 10mm
- Thread type: 10-425 thread
- Flat bottom

Supplied in packs of 100

Code	Alt Ref	Description	Price
CHR2660	SCV12-02L	Clear	£146.00
CHR2662	SCV12-02LA	Amber	£156.00



Screw Neck Caps

Supplied in packs of 1000

Code	Alt Ref	Description	Price
CHR3418	SCCSEP10-02	10mm combination seal PP screw cap black centre hole 10-425 thread natural rubber red-orange/TEF transparent 60 shore A 1.3mm	£194.00
CHR3428	SCCSEP10-04S	10mm combination seal PP screw cap black centre hole 10-425 thread silicone white/PTFE blue 55 shore A 1.5mm slitted	£483.00



Headspace Vials

- Clear glass

Supplied in packs of 1000

Code	Alt Ref	Volume, mL	Bottom	Neck Style	Dims, h x Ø, mm	Price
CHR2516	CRV10-01	10	Flat	Crimp	46 x 23	£250.00
CHR2584	CRV20-01C	20	Round	Crimp	76 x 23	£321.00
CHR2726	SNV22-25	25	Flat	Snap	65 x 26	£224.00



Glass Vial Inserts

Supplied in packs of 1000

Code	Alt Ref	Description	Price
CHR2974	INWC-01	Insert wide opening vial 0.1mL [31 x 6mm] conical	£182.43
CHR2984	INWF-01	Insert wide opening vial 0.3mL [31 x 6mm] flat	£253.00