



Epredia™ Wax Dispensers

epredia | Authorized Distributor

Offered in two capacities: a 7.5L unit and a 15L, dual-well unit. Both units offer digital temperature control and an ultra-fast heating system ensures rapid melting of pelletised wax, with melted wax readily available thirty minutes after the wax dispenser is filled.

- Non-drip, three-position lever tap-push dispense, lock open for rapid delivery, or lock closed for safety
- High tap position (153mm from the bench) accommodates large containers for transferring molten wax
- Integral tap filter screen (0.5mm mesh) prevents coarse particles from blocking the delivery tap
- Anti-microbial coating (Ag⁺) inhibits growth of bacteria
- Ultra-fast heating system minimises wait times
- Digital temperature control—high-capacity model features individual controls for each of the dual reservoirs
- Aluminium inner tank provides good thermal conductivity and improved temp control
- Insulated paraffin tank prevents heat loss from the tank and ensures that outer surfaces of the tank are safe to touch
- Dedicated tap-heating system prevents blockages and ensures even flow of wax

Code	Alt Ref	Description	Capacity, L	Price
HIS3040	A84600001	Slimline wax dispenser	7.5 (6kg)	£1928.00
HIS3042	A84600010	High capacity wax dispenser	15 (12kg)	£2597.00



Epredia™ HM 325 Rotary Microtomes

epredia | Authorized Distributor

Cut high-quality paraffin sections without a power supply. The HM 325 features retraction at the return travel to protect the specimen and manual coarse advance for fast working.

- One-handed zero orientation head
- Quick-trim button
- Section counter
- Guide-rail guidance for the knife holder
- Large, easy-to-clean waste tray extends below the blade holder
- All supplied with standard tools, brushed aluminium cover plate, dust cover and instruction manual
- Note: blade/knife holder must be ordered separately unless purchased as a package

Specifications	
Total Specimen Advance, mm	28
Sectioning Thickness Range, µm	0.5 to 60
Mechanical Trim Stages, µm	10/30
Retraction, µm	40
Dims, w x d x h, mm	490 x 420 x 280

Code	Alt Ref	Supplied with	Price
HIS3002	902100	—	£6896.00
HIS3003	902100A	Universal cassette clamp (715020S), disposable blade carrier E (705800)	£10545.00
HIS3005	902100ER	Universal cassette clamp (715020S), disposable blade carrier ER (705830)	£10456.00

Accessories

Code	Alt Ref	Description	Price
HIS3024	715010	Standard clamp, silver	£486.00
HIS3026	715020	Universal cassette, quick release clamp, orienting head	£565.00
HIS3030	716120	Universal adjustable macro cassette clamp	£741.00
HIS3162	705830	Disposable blade carrier (ER)	£1416.00