

Corning® Ice Buckets and Pans



Non-toxic, recyclable ethyl-vinyl acetate (EVA) foam containers for use with ice, dry ice, liquid nitrogen, alcohol slurries. A great non-toxic alternative to conventional ice pans. EVA has little or no odour and is inert. Containers will not “sweat” or deform and have a textured finish on the bottom to prevent slipping on the benchtop.

- Highly insulative, lightweight and leak-proof
- Stackable
- Ideal for use with ice, dry ice (-78°C), liquid nitrogen (-196°C), alcohol or saline solutions, or warm solutions (+93°C)
- Can be cleaned with soap and water, 10% bleach or 70% alcohol
- Supplied without lid unless otherwise stated



Rectangular Pans

Code	Description	Capacity, L	Colour	Price
432115	Mini	1	Green	£70.00
432116	Mini	1	Blue	£70.00
432117	Mini	1	Red	£70.00
432118	Mini	1	Orange	£70.00
432119	Mini	1	Lime green	£70.00
432120	Mini	1	Pink	£70.00
432121	Mini	1	Purple	£70.00
432103	Midi	4	Green	£109.00
432104	Midi	4	Blue	£109.00
432105	Midi	4	Red	£109.00
432106	Midi	4	Orange	£109.00
432107	Midi	4	Lime green	£109.00
432108	Midi	4	Pink	£109.00
432109	Midi	4	Purple	£109.00
432110	Midi with lid	4	Blue	£133.00
432111	Midi with lid	4	Orange	£133.00
432112	Midi with lid	4	Lime green	£133.00
432113	Midi with lid	4	Pink	£133.00
432114	Midi with lid	4	Purple	£133.00
432093	Maxi	9	Green	£146.00
432094	Maxi	9	Blue	£146.00
432095	Maxi	9	Red	£146.00
432096	Maxi	9	Orange	£146.00
432097	Maxi	9	Lime green	£146.00
432098	Maxi	9	Pink	£146.00
432099	Maxi	9	Purple	£146.00
432100	Maxi with lid	9	Blue	£169.00
432101	Maxi with lid	9	Lime green	£169.00
432102	Maxi with lid	9	Purple	£169.00