



Corning® Falcon® Bacteriological Petri Dishes

FALCON

- Durable construction for stable dish manipulation
- Crystal-grade, virgin polystyrene for optical clarity
- Sterile

Code	Description	Dims, Ø x h, mm	Inner Pack	Pack	Price
351008	Easy-Grip	35 x 10	20	500	£145.00
351006	Tight-fit lid	50 x 9	20	500	£205.00
351007	Standard	60 x 15	20	500	£139.00
351029	Standard	100 x 15	20	500	£178.00
351058	Standard	150 x 15	10	100	£120.00



A

Glass Petri Dishes

- Special hybrid glass which can be autoclaved up to 160°C

Code	Dims, Ø x h, mm	Pack	Price
PET1002	60 x 12	18	£23.20
PET1004	80 x 15	18	£26.60
PET1006	90 x 15	18	£30.70
PET1008	100 x 15	18	£27.40
PET1010	120 x 20	18	£46.20
PET10120	150 x 25	18	£74.00
PET10142	200 x 30	6	£47.60




A

Nickel Electro Petri Dish Sterilising Boxes



- Capacity to hold 10 Petri dishes in a removable carrier for transportation and sterilisation
- Can be steam sterilised at 120°C and withstand dry heat temperatures of up to 180°C
- Dims, h (external) x Ø: 254 x 115mm

Code	Alt Ref	Material	Price
PET4050	6301	Aluminium	£77.00
PET4052	6501	Stainless steel	£111.00



Petri Dish Holder

- Comprising of three compartments, each holding up to 10 Petri dishes, maximum 100mm diameter
- Fold down stainless steel handle for transport
- Inverted base

Code	Compartments	Dims, w x d x h, mm	Price
PET4000	3	110 x 325 x 170	£36.50