

MBL® Class B Bulb Pipettes, One Mark

MBL®
VOLUMETRICS

- One mark
- Calibrated to Class B tolerances
- Complies with ISO 648 and DIN 12690
- Calibrated for delivery
- Manufactured from soda-lime glass with tooled jets to reduce breakage
- Markings in permanent amber stain diffused into the surface of the glass for increased chemical resistance


Supplied in packs of 5



Code	Alt Ref	Capacity, mL	Delivery Time, secs	Pack	Price
PIP11000	PRA700B	1	5 - 20	5	£14.60
PIP11020	PRA702B	2	5 - 25	5	£14.60
PIP11080	PRA704B	5	7 - 30	5	£15.40
PIP11100	PRA706B	10	8 - 40	5	£18.50
PIP11120	PRA708B	15	9 - 50	5	£21.80
PIP11140	PRA710B	20	9 - 50	5	£23.10
PIP11160	PRA712B	25	10 - 50	5	£23.60
PIP11180	PRA714B	50	13 - 60	5	£28.20
PIP11200	PRA716B	100	25 - 60	2	£21.10

Volac Class AS Bulb Pipettes, One Mark

- One mark
- Class AS
- Clear glass
- Conformity approved
- With batch certificate
- Blue graduations



Code	Capacity, mL	Tolerance, ± mL	Works Certificate	Pack	Price
PIP12000	1	0.008	No	5	£77.00
PIP12300	1	0.008	Yes	5	£55.00
PIP12020	2	0.01	No	5	£54.00
PIP12320	2	0.01	Yes	5	£55.00
PIP12040	3	0.015	No	5	£31.20
PIP12340	3	0.015	Yes	5	£58.00
PIP12060	4	0.015	No	5	£46.60
PIP12360	4	0.015	Yes	5	£66.00
PIP12080	5	0.015	No	5	£61.00
PIP12380	5	0.015	Yes	5	£60.00
PIP12100	10	0.02	No	5	£49.10
PIP12400	10	0.02	Yes	5	£69.00
PIP12120	15	0.025	No	5	£55.00
PIP12140	20	0.03	No	5	£58.00
PIP12440	20	0.03	Yes	5	£76.00
PIP12160	25	0.03	No	5	£59.00
PIP12460	25	0.05	Yes	5	£77.00
PIP12180	50	0.05	No	5	£64.00
PIP12480	50	0.08	Yes	5	£97.00
PIP1220	100	0.08	No	Each	£57.00
PIP1250	100	0.08	Yes	5	£76.00