

## Corning® Stripette® Serological Pipettes



- Certified non pyrogenic and DNase/RNase free
- Sterile, individually wrapped
- Bi-directional graduations provide choice of ascending and descending scales
- Negative graduations allow additional working volume
- Colour coded magnifier stripes make volume reading easier
- Volumetric accuracy of  $\pm 2\%$  at full capacity



### Plastic Wrapped

- Plastic wrapped for reduced static cling and particulate-free opening
- Single bin recycling

Code	Volume, mL	Graduations, mL	Negative Graduations, mL	Colour	Inner Pack	Pack	Price
4012	1	0.01	0.2	Yellow	100	200	£59.00
4011	1	0.01	0.2	Yellow	100	1000	£215.00
4021	2	0.01	0.2	Green	100	1000	£221.00
4051	5	0.01	2.5	Blue	50	200	£68.00
4101	10	0.1	3	Orange	50	200	£69.00
4251	25	0.2	10	Red	50	200	£137.00
4501	50	0.5	10	Purple	25	100	£150.00
4484	100	1	—	Aqua	10	100	£638.00

### Paper/Plastic Wrapped

Code	Volume, mL	Graduations, mL	Negative Graduations, mL	Colour	Inner Pack	Pack	Price
4485	1	0.01	0.2	Yellow	50	1000	£215.00
4486	2	0.01	0.2	Green	50	1000	£239.00
4487	5	0.1	2.5	Blue	50	200	£74.00
4488	10	0.1	3	Orange	50	200	£74.00
4489	25	0.2	10	Red	25	200	£137.00
4490	50	0.5	10	Purple	25	100	£152.00
4491	100	1	—	Aqua	10	100	£638.00

### Triple Paper/Plastic Wrapped

- With an inner "dust cover" liner for the carton as well as a second sealed sleeve over each inner package, eliminating the need for alcohol wipe downs and similar procedures required in clean room production facilities
- Ideal for cell-based vaccine and drug production facilities, or other laboratory clean rooms
- Sterility Assurance Level (SAL) of  $10^{-6}$

Code	Volume, mL	Graduations, mL	Negative Graduations, mL	Colour	Inner Pack	Pack	Price
7041	1	0.01	0.2	Yellow	50	1000	£294.00
7042	2	0.01	0.2	Green	50	1000	£333.00
7045	5	0.1	2.5	Blue	50	200	£123.00
7015	10	0.1	3	Orange	50	200	£107.00
7016	25	0.2	10	Red	25	200	£224.00
7017	50	0.2	10	Purple	25	100	£235.00
7000	100	1	—	Aqua	25	100	£454.00