



Thermo Scientific™ S1 Pipette Fillers

Thermo
SCIENTIFIC

Lightweight and cordless, it features a heavy-duty rechargeable lithium-ion battery that has a rapid charge rate of 80% in one hour and delivers up to 15 hours of continuous operation. For use with 1 to 100mL serological pipettes.

- Filter type: 0.45µm standard
- Large backlit LCD display
- Eight speeds and gravity dispense to choose from
- Supplied with universal charger, table stand, wall-mounted holder, 1mL pipette support, instruction manual and warranty card

Code	Alt Ref	Colour	Price
PIP7490	9501	White	£591.00
PIP7492	9511	Clear	£591.00
PIP7494	9521	Blue	£591.00
PIP7496	9531	Red	£591.00
PIP7498	9541	Green	£591.00



Pi-Pump Pipette Fillers

A high quality acid resistant plastic instrument which is easily dismantled for cleaning. Simple and easy to use with finely adjustable thumb wheel and release valve lever. Available in four different sizes with each size being colour coded for easy identification.

- Innovative, easy to use and cost-efficient solution for handling liquids
- Designed specifically to allow for pipetting to be done easily and safely by hand
- Manufactured, and assembled and tested in Germany
- Individually packaged and supplied with usage instructions


Code	Capacity, mL	Colour Code	Price
PIP6150	Up to 0.2	Yellow	£18.14
PIP6152	Up to 2	Blue	£7.00
PIP6154	Up to 10	Green	£11.14
PIP6156	Up to 25	Red	£11.10



Rubber Pipette Fillers

- Simple pipetting aids for one-mark and graduated pipettes
- Control the functions by squeezing the appropriate valves between thumb and forefinger

Code	Description	Price
PIP6000	Standard model up to 10mL	£12.38
PIP6002	Universal up to 100mL	£16.00



PVC Pipette Filler

- Saffron pattern (bulb on side) is a low-cost design to simplify the problems of pipette filling

Code	Material	Price
PIP6010	PVC	£38.00