


Thermo Scientific™ Finntip™ Tips

Designed to provide perfect sealing with optimal pipetting accuracy and precision.

- Manufactured from highest quality virgin polypropylene raw materials
- Covers 0.2µL to 10mL volumes with extended length options
- Colour coded racks to match the compatible Finnpiquette™
- The products detailed below represent a proportion of the full product range available, please see website for full details, including Finntip™ Flex and wide orifice options



Loose Tips

Code	Alt Ref	Description	Volume Range, µL	Sterility	Pack	Price
PIP3600	9400310	Finntip 10	0.2 - 10	Non sterile	1000	£74.00
PIP3604	9400620	Finntip 20	0.2 - 20	Non sterile	1000	£96.00
PIP3606	9400360	Finntip 50	0.2 - 50	Non sterile	1000	£96.00
PIP3612	9400230	Finntip 250 universal	0.5 - 250	Non sterile	1000	£72.00
PIP3618	9401260	Finntip 300	5 - 300	Non sterile	1000	£75.00
PIP3622	9401030	Finntip 1000	100 - 1000	Non sterile	1000	£73.00
PIP3626	9401410	Finntip 1000 extended length	100 - 1000	Non sterile	1000	£172.00
PIP3636	9402030	Finntip 5mL	500 - 5000	Non sterile	500	£98.00
PIP3640	9402151	Finntip 10mL	1000 - 10000	Non sterile	100	£35.00

Non Filtered, Racked Tips

Code	Alt Ref	Description	Volume Range, µL	Sterility	Pack	Price
PIP3644	9400300	Finntip 10	0.2 - 10	Non sterile	10 x 96	£98.00
PIP3712	9400303	Finntip 10	0.2 - 10	Sterile	10 x 96	£124.00
PIP3656	9400610	Finntip 20	0.2 - 20	Non sterile	10 x 384	£600.00
PIP3716	9400613	Finntip 20	0.2 - 20	Sterile	10 x 384	£769.00
PIP3658	9400370	Finntip 50	0.2 - 50	Non sterile	10 x 384	£585.00
PIP3718	9400373	Finntip 50	0.2 - 50	Sterile	10 x 384	£716.00
PIP3660	9400130	Finntip 200 extended length	5 - 200	Non sterile	10 x 96	£194.00
PIP3720	9400133	Finntip 200 extended length	5 - 200	Sterile	10 x 96	£244.00
PIP3668	9400260	Finntip 250 universal	0.5 - 250	Non sterile	10 x 96	£93.00
PIP3724	9400263	Finntip 250 universal	0.5 - 250	Sterile	10 x 96	£113.00
PIP3676	9401250	Finntip 300	5 - 300	Non sterile	10 x 96	£96.00
PIP3728	9401253	Finntip 300	5 - 300	Sterile	10 x 96	£118.00
PIP3686	9401110	Finntip 1000	100 - 1000	Non sterile	10 x 96	£100.00
PIP3732	9401113	Finntip 1000	100 - 1000	Sterile	10 x 96	£125.00
PIP3690	9401420	Finntip 1000 extended length	100 - 1000	Non sterile	5 x 96	£93.00
PIP3734	9401423	Finntip 1000 extended length	100 - 1000	Sterile	5 x 96	£117.00
PIP3708	9402070	Finntip 5mL	500 - 5000	Non sterile	5 x 54	£96.00
PIP3742	9402073	Finntip 5mL	500 - 5000	Sterile	5 x 54	£121.00
PIP3710	9402160	Finntip 10mL	1000 - 10000	Non sterile	5 x 24	£64.00
PIP3744	9402163	Finntip 10mL	1000 - 10000	Sterile	5 x 24	£77.00