




## ecoSHIELD™ Eco Nitrile PF 250 Gloves



- Green nitrile
- Powder and latex free
- Textured fingertips for grip
- Beaded cuff
- Ambidextrous
- Recycled cardboard packaging with less printing to reduce the polluting de-inking recycling process
- Enhanced design to package 50% more gloves in the same volume
- twinSHIELD™ double walled protection
- Length: 250mm (S - L), 260mm (XL)
- Category III
- AQL 0.25
- EN ISO 374-1:2016 + A1:2018 Type B JKPT (chemical risks)
- EN 16523-1:2015 + A1:2018 (chemical permeation)
- ISO 374-4:2019 (resistance to degradation by chemicals)
- EN ISO 374-5:2016 VIRUS (protection from viruses, bacteria and fungi)
- ISO 16604:2004 Procedure B (viral penetration test)
- EN 455-1/2/3/4:2000 (medical gloves for single use)

Supplied in packs of 150

Code	Alt Ref	Size	Price
SAF4320	625121	X small (6)	£55.00
SAF4322	625122	Small (7)	£55.00
SAF4324	625123	Medium (8)	£55.00
SAF4326	625124	Large (9)	£55.00
SAF4328	625125	X large (10)	£55.00



## duoSHIELD™ ICE NITRILE™ 240 Gloves



- Blue nitrile
- Powder and latex free
- Textured fingertips for grip
- Beaded cuff
- Ambidextrous
- Length: 240mm
- Category III
- AQL 1.5
- EN ISO 374-1:2016 + A1:2018 Type B KPT (chemical risks)
- EN 16523-1:2015 + A1:2018 (chemical permeation)
- ISO 374-4:2019 (resistance to degradation by chemicals)
- EN ISO 374-5:2016 VIRUS (protection from viruses, bacteria and fungi)
- ISO 16604:2004 Procedure B (viral penetration test)
- EN 455-1/2/3/4:2000 (medical gloves for single use)

Supplied in packs of 100

Code	Alt Ref	Size	Price
SAF1058	655121	X small (6)	£66.00
SAF1060	655122	Small (7)	£66.00
SAF1062	655123	Medium (8)	£66.00
SAF1064	655124	Large (9)	£66.00
SAF1066	655125	X large (10)	£66.00