



PYREX®

Safety in the Lab: Glassware

PYREX® Plastic Coated Safety Bottles, Clear Glass

These laboratory bottles have an external plastic coating that provides protection from spills and splashes if the bottle is broken. The plastic coating will contain the bottle contents until it can be transferred to another bottle.

- Bottles only. Supplied without cap or pouring ring
- Complies with ISO 4796
- Borosilicate glass with clear PU plastic coating to provide added protection against mechanical breakage
- Reduces the hazard of flying glass and liquid splashes in the event of breakage
- Printed with Trace Code for downloadable batch certificate
- Coating withstands steam sterilisation up to 135°C



SLS Code	Supplier Part No	Capacity	Pack	Offer Price,€
WEB2300	1518/04	100mL	10	129.66
WEB2302	1518/06	250mL	10	136.46
WEB2304	1518/08	500mL	10	165.18
WEB2306	1518/10	1L	10	237.81
WEB2308	1518/12	2L	10	487.72
WEB2310	1518/16	5L	Each	250.39
WEB2312	1518/18	10L	Each	413.17

PYREX® Accessories

SLS Code	Supplier Part No	Description	Pack	Offer Price,€
WEB2314	4506/45D	Screwcap, GL45, Blue PP	10	22.23
WEB2316	4516/45D	Clear Pouring Ring GL45, PP	10	9.75
WEB2318	4526/45D	Screwcap, GL45, PTFE Seal, Red	10	69.15
WEB2320	4536/45D	ETFE Pouring Ring, GL45, Red	10	17.68

Offers valid until 30th June 2023.

United Kingdom:
slsinfo@scientific-labs.com

International:
slsexport@scientific-labs.com

Northern Ireland:
NI.sales@scientific-labs.com

Republic of Ireland:
ROI.sales@scientific-labs.com





Safety in the Lab: Glassware

DURAN® Protect Plastic Coated (PU) Amber Glass Bottles

With easy-to-read scale. Fired-on, highly durable white ceramic. Service temperature limit of the PU plastic coating: -30°C to $+135^{\circ}\text{C}$. The coating provides scratch, leak and splinter protection and is ideally suited to both the transport and storage of hazardous media or valuable samples. UV protection up to approx. 500nm wavelength. Unchanged DURAN® properties within the bottle, as colouration is only on the outer surface. Very uniform, durable and chemically resistant amber colour due to use of innovative technology. Suitable for microwaving. Typical applications: storage, transport and safe handling of hazardous or valuable substances.



SLS Code	Supplier Part No	Capacity	Pack	Offer Price,€
WEB2326	218062438	100	10	350.52
WEB2328	218063631	250mL	10	456.69
WEB2330	218064439	500mL	10	566.67
WEB2332	218065435	1L	10	736.13
WEB2334	218066337	2L	10	1133.54
WEB2336	218067333	5L	Each	365.79

DURAN® Accessories

For DURAN® laboratory glass bottles with DIN thread

SLS Code	Supplier Part No	Description	Pack	Offer Price,€
WEB2346	292392809	Screw-caps with lip seal, GL45, PP, blue	10	16.50
WEB2348	292422809	Pouring ring, GL45, PP, blue	10	5.38
WEB2350	292402807	Screw-caps, GL45, PBT, red, with PTFE protected seal	10	66.00
WEB2352	292442802	Pouring ring, GL45, ETFE, red	10	13.86

Offers valid until 30th June 2023.

United Kingdom:
slsinfo@scientific-labs.com

International:
slsexport@scientific-labs.com

Northern Ireland:
NI.sales@scientific-labs.com

Republic of Ireland:
ROI.sales@scientific-labs.com





Safety in the Lab: Glassware

DURAN® Pressure Plus+ Pressure-Resistant Clear Glass Bottles

The DURAN® pressure plus+ GL 45 protect coated bottles offer a safer option for laboratory applications involving a vacuum or pressure. Standard glass laboratory bottles are not designed to safely withstand pressure differentials. As a result of a modified geometry, the pressure plus+ bottles have a guaranteed pressure-resistance, and consequently offer optimal safety for both the user and contents. The special geometry still conforms to ISO 4796-1:2016 to maintain instrument compatibility. Certified in accordance with DIN EN 1595. The bottles are moulded from high-purity, clear, USP/EP Type I borosilicate glass 3.3, renowned for its excellent thermal performance, and high resistance to chemical attack. The glass is coated with a clear PU plastic coating that offers protection for the user and contents in the event of breakage due to impact.



SLS Code	Supplier Part No	Capacity	Pack	Offer Price,€
WEB2338	218152402	100mL	10	159.46
WEB2340	1175925	250mL	10	205.31
WEB2342	1175926	500mL	10	256.74
WEB2344	218155408	1L	10	338.60

DURAN® Accessories

For DURAN® laboratory glass bottles with DIN thread

SLS Code	Supplier Part No	Description	Pack	Offer Price,€
WEB2346	292392809	Screw-caps with lip seal, GL45, PP, blue	10	16.50
WEB2348	292422809	Pouring ring, GL45, PP, blue	10	5.38
WEB2350	292402807	Screw-caps, GL45, PBT, red, with PTFE protected seal	10	66.00
WEB2352	292442802	Pouring ring, GL45, ETFE, red	10	13.86

Offers valid until 30th June 2023.