


Special Instruction

Filtration · Rapid Tests · Water Analysis · Chromatography · Bioanalysis
 Filtration · Schnellteste · Wasseranalytik · Chromatographie · Bioanalytik

Edition: 04.08.2010

Parameter	Peroxide
Application	Determination of Peroxide in solvents
Product group	<i>QUANTOFIX</i> [®]
Reagents and equipment	- <i>QUANTOFIX</i> [®] Peroxide
Procedure	<ol style="list-style-type: none"> 1. Dip the test strip into the test solution for 1 second. 2. Shake off excess liquid. 3. Allow solvent to dry completely. 4. Once dry, moisten the test field with a drop of water. 5. Wait 15 seconds. 6. Compare with the colour scale. If hydrogen peroxide is present, the test field turns blue.



MACHERY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6-8 · D-52355 Düren · Germany
Germany and international:
 Tel.: +49 (0) 24 21 96 90
 Fax: +49 (0) 24 21 96 91 99
 e-mail: sales-de@mn-net.com

Switzerland:
MACHERY-NAGEL AG
 Tel.: +41 (0) 62 388 55 00
 Fax: +41 (0) 62 388 55 05
 e-mail: sales-ch@mn-net.com

France:
MACHERY-NAGEL EURL
 Tel.: +33 (0) 3 88 68 22 68
 Fax: +33 (0) 3 88 51 76 88
 e-mail: sales-fr@mn-net.com

USA:
MACHERY-NAGEL Inc.
 Tel.: +1 484 821 09 84
 Fax: +1 484 821 12 72
 e-mail: sales-us@mn-net.com