

# hybex™

## Media Storage Bottles



**Borosilicate  
Glass**

— Borosilicate 3.3 hybex glass

— Autoclavable caps (GL45) & drip rings

— Permanent enamel graduation markings

— Large area for bottle marking/identification

— Autoclavable, reuseable & breakage resistant

— Five cap colors

### Technical Data

Material: hybex™ borosilicate glass (3.3)  
Cap Thread Size: GL45\*  
Temp. Range: -70°C to 500°C  
Sterilization: Autoclave or dry heat  
Cap Material: Polypropylene (140°C)  
Cap Colors: Blue, green, purple  
yellow, white or assorted

\*50ml bottles are only supplied with blue GL32 sealing caps

### Also Available

hybex™  
Erlenmeyer Flasks



### Ordering Information

B3000-50 Media storage bottle, 50ml with standard (GL32) blue cap, 10/pk.  
B3000-100\* Media storage bottle, 100ml with standard (GL45) cap, 10/pk.  
B3000-250\* Media storage bottle, 250ml with standard (GL45) cap, 10/pk.  
B3000-500\* Media storage bottle, 500ml with standard (GL45) cap, 10/pk.  
B3000-1000\* Media storage bottle, 1L with standard (GL45) cap, 10/pk.  
B3000-2000\* Media storage bottle, 2L with standard (GL45) blue cap, 5/pk.  
B3000-5L Media storage bottle, 5L with standard (GL45) blue cap, 1/ea.  
B3000-10L Media storage bottle, 10L with standard (GL45) blue cap, 1/ea.  
B3000-SP\* Media storage bottle, *Starter Pack* with standard (GL45) cap  
(*Starter Pack* of ten bottles includes 2x100ml, 3x250ml, 3x500ml & 2x1000ml)

B3000-HTC High temperature red caps, up to 180°C (GL45), 10/pk  
B3000-CAP\* Replacement cap, standard blue (GL45), 10/pk.  
B3000-RIN Replacement sealing ring, (GL45), 10/pk.  
B3000-CAP2 Replacement cap, standard blue (GL32), 10/pk.  
B3000-RIN2 Replacement sealing ring, (GL32), 10/pk.

\*Please add (-B) for blue caps, (-G) for green caps, (-P) for purple caps, (-Y) for yellow caps, (-W) for white caps or (-AS) for assorted caps (2 of ea. color)

**Benchmark**  
Scientific



**Benchmark Scientific Inc.**  
PO Box 709 Edison, NJ 08818  
Phone: 908-769-5555 Fax: 908-222-1864  
Web: [www.BenchmarkScientific.com](http://www.BenchmarkScientific.com)  
Email: [Info@BenchmarkScientific.com](mailto:Info@BenchmarkScientific.com)