

# **UNIVERSAL 320**

### BLOOD / CELL CULTURE / URINE / CYTOLOGY / PLATES

# Universal Appeal

The UNIVERSAL 320 is a compact, versatile and indispensable general-purpose centrifuge. Excellent performance and a comprehensive range of accessories enable the UNIVERSAL 320 to carry out virtually any centrifuging task.

MAX. CAPACITY	4 x 200 mL / 6 x 85 mL
MAX. RPM/RCF	15,000/21,382
TEMP. CONTROL	from -20 to +40°C (320 R)





#### **QUIET OPERATION**

Quiet as a normal conversation. Will not increase noise levels or contribute to distractions.



#### **FIVE YEAR WARRANTY**

All Hettich centrifuges are protected by our 5-year warranty and full technical support.



#### **FUNCTIONAL DESIGN**

Easy to lift lid with front-facing control panel for effortless interaction and programming.



#### **FULLY CERTIFIED**

Engineered and manufactured in accordance with all relevant safety and production standards.



#### **ENHANCED SAFETY**

Self-closing, auto-locking lid with steel latches. All-metal chamber provides reinforced security.



#### **MADE IN GERMANY**

German engineered and manufactured to be quiet, reliable and safe since 1904.

## HETTICH UNIVERSAL 320 APPLICATION PACKAGES



#### UNIVERSAL 320 CONICAL TUBE PACKAGE - 5 mL - BIO-CONTAINMENT

|--|

Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF
16	5 Conical	17 x 59	≤ 4,500	≤ 3,328
Includes: (1) Centrifuge. (1) Horizontal Rotor. (4) Buckets. (4) Adapters 5 mL conical. (4) Bio-containment Lids				

PKG. NO.	PKG. NO. (TEMP. CONTROL)
320CONICAL5-BC	320RCONICAL5-BC

#### UNIVERSAL 320 BLOOD TUBE PACKAGE



Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF
28	1.6 - 7	13 x 75/100	≤ 5,000	≤ 3,913
28	4 - 10	16 x 75/100	≤ 5,000	≤ 3,913
Includes: (1) Centrifuge, (1) Horizontal Rotor, (4) Buckets, (4) Adapters 13 x 75/100 mm, (4) Adapters 16 x 75/100 mm				

PKG. NO.	PKG. NO. (TEMP. CONTROL)	
320BL00D	320RBL00D	

#### UNIVERSAL 320 BLOOD TUBE PACKAGE - 4 - BIO-CONTAINMENT



Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF
28	1.6 - 5	13 x 75/100	≤ 4,500	≤ 3,215
20	4 - 10	16 x 75/100	≤ 4,500	≤ 3,215
Includes (I) Outside (A) Universal Data (A) Publish (I) Advance 10 or 75 (100 or re)				

PKG. NO.	PKG. NO. (TEMP. CONTROL)
320BL00D4-BC	320RBL00D4-BC

#### UNIVERSAL 320 URINALYSIS PACKAGE



Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF
12	12 KOVA	17 x 102	≤ 4,500	≤ 3,102
Includes: (1) Centrifuge,	(1) Horizontal Rotor, (4) Bu	ckets, (4) Adapters KOVA		

PKG. NO.	PKG. NO. (TEMP. CONTROL)
320URINALYSIS	320RURINALYSIS

#### UNIVERSAL 320 CELL CULTURE PACKAGE - 4



Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF
16	15 Conical	17 x 120	≤ 4,500	≤ 3,305
4 50 Conical 29 x 115 ≤ 4,500 ≤ 3,2				
Includes: (1) Centrifuge, (1) Horizontal Rotor, (4) Buckets, (4) Adapters 15 mL, (4) Adapters 50 mL				

PKG. NO.	PKG. NO. (TEMP. CONTROL)
320CELLCULTURE4	320RCELLCULTURE4

#### UNIVERSAL 320 CELL CULTURE PACKAGE - 4 - BIO-CONTAINMENT



Tube Qty	Volume (mL)	Size (mm)	RPM	RCF
8	15 Conical	17 x 120	≤ 4,500	≤ 3,328
4	50 Conical	29 x 115	≤ 4,500	≤ 3,328
Includes: (1) Centrifuge, (1) Horizontal Rotor, (4) Buckets, (4) Adapters 15 mL, (4) Adapters 50 mL, (4) Bio-containment Lids				

PKG. NO.	PKG. NO. (TEMP. CONTROL)
320CELLCULTURE4-BC	320RCELLCULTURE4-BC

#### UNIVERSAL 320 CELL CULTURE PACKAGE - 8



Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF
8	15 Conical	17 x 120	≤ 4,500	≤ 3,260
8	50 Conical	29 x 115	≤ 4,500	≤ 3,260
Includes: (1) Centrifuge, (1) Horizontal Rotor, (4) Buckets, (8) Adapters 15 mL, (8) Adapters 50 mL				

PKG. NO.	PKG. NO. (TEMP. CONTROL)
320CELLCULTURE8	320RCELLCULTURE8

#### UNIVERSAL 320 PLATE PACKAGE



Τι	ıbe Qty.	Volume (mL)	Size (mm)	RPM	RCF
10	)	Microtiter		≤ 4,000	≤ 2,218
Includes: (1) Centrifuge, (1) Horizontal Rotor, (2) Buckets					

PKG. NO.	PKG. NO. (TEMP. CONTROL)
320PLATE	320RPLATE

#### UNIVERSAL 320 CUSTOM PACKAGE



Supports over 413 application configurations, please inquire for more information

Tube types include: Microliter up to 5 mL, Blood, Pediatric, Urine, Conical and Skirted up to 50 mL, Standard up to 100 mL, Plates, Bottles up to 200 mL, Cytology

J

#### hettweb.com/universal320

CONTACT YOUR LOCAL HETTICH REPRESENTATIVE TO LEARN MORE:



