


THE PATH FROM VEGETABLE COOKING OIL TO EPPENDORF TUBES® BIOBASED



> The food oil wastes and residues used to produce the renewable feedstock based polypropylene are traceable to the first collection points.

> The renewable feedstock based polypropylene and all manufacturing steps are sustainability certified by ISCC PLUS – a reliable, world-leading certification system. The renewable feedstock content is determined via the ISCC mass balance approach.

> The resulting materials for the polyolefins used meet the same strict specifications, purity and quality, regardless of whether they are fossil-based or renewable feedstock based recycled raw materials.

