

dispensman™



BOTTLE-TOP DISPENSER



Dispense  
Liquids Simply  
and Accurately

| Precise

| Safe

| No Loss of Liquid

| Simple & Quick Operation

 **GILSON®**



## Accurate & Reproducible

Dispense just the liquid needed every time with a quick, easy, and precise volume setting.

## Safe

Keep your bench and other devices safe with the anti-drip system and the return of liquid into the bottle after use.

## No Loss of Liquid

Prime your operation by purging the air in RECYCLE position keeping the liquid inside the bottle.

## Easy & Comfortable Use

Simple one hand operation and volume adjustment with the sliding volume adjuster.

## Flexible

- Adaptable to every type of liquid with an easy calibration by the user.
- Numerous bottle adaptors.
- Universal chemical compatibility for many applications.
- Fully autoclavable at 121°C.



Convenient 3-position nozzle:

- DISPENSE
- RECYCLE
- EMPTY

MODEL	PART NUMBER	VOLUME RANGE	SCALING (mL)	MAXIMUM PERMISSIBLE ERRORS			
				GILSON		ISO 8655-5	
				SYSTEMATIC ERROR (mL)	RANDOM ERROR (mL)	SYSTEMATIC ERROR (mL)	RANDOM ERROR (mL)
DISPENSAMAN 2.5 mL	<b>F110101</b>	0.25-2.5 mL	0.05	±0.012	≤ 0.002	±0.030	≤ 0.010
DISPENSAMAN 5.0 mL	<b>F110102</b>	0.5-5 mL	0.10	±0.030	≤ 0.005	±0.030	≤ 0.010
DISPENSAMAN 10 mL	<b>F110103</b>	1-10 mL	0.20	±0.060	≤ 0.010	±0.060	≤ 0.020
DISPENSAMAN 25 mL	<b>F110104</b>	2.5-25 mL	0.50	±0.150	≤ 0.025	±0.150	≤ 0.050
DISPENSAMAN 50 mL	<b>F110105</b>	5-50 mL	1.00	±0.300	≤ 0.050	±0.300	≤ 0.100

### Gilson S.A.S.

19, avenue des Entrepreneurs BP 145 • F-95400 Villiers-le-Bel, France  
T +33 (0) 1 34 29 50 00 • F +33 (0) 1 34 29 50 20

[www.gilson.com/ContactUs](http://www.gilson.com/ContactUs)

### Gilson, Inc.

3000 Parmenter Street • PO Box 620027  
Middleton, WI 53562-0027, USA  
T 608-836-1551 or 800-445-7661 • F 608-831-4451



LT800563 © 06/2017 Gilson SAS.  
All rights reserved. All trademarks are the property of Gilson Inc. and its subsidiaries.

