



Petrifilm™

3M Food Safety Lab Automation

Supports a wide assortment of food types

Automated reporting helps users turn results into insights with easy to use software*



Enumeration in
6
seconds
or less

Productivity. Advanced.

3M™ Petrifilm™ Plate Reader Advanced

There's no slowing down with automated plate enumeration and advanced reporting capabilities of the 3M™ Petrifilm™ Plate Reader Advanced. It can help you save time, improve productivity and capacity.

It's time—upgrade to advanced testing at 3M.com/LabAutomation

*Computer not included.

Advancing the efficiency of food safety.

Automate the imaging, enumeration and storage of results of 3M™ Petrifilm™ Plate tests.



- LED bar indicates reader status
- Designed to be easy to clean and maintain
- Fixed artificial intelligence (AI) networks to enumerate 3M™ Petrifilm™ Plates
- Reads 10 barcode symbologies



Productive

Improve laboratory efficiency in a space saving design

**Plate Enumeration in 6 seconds or less—
up to 900 plates per hour**



Proficient

Standardized technology and easy to use software

Accurate colony counting




Proactive

Software stores results and offers data analytics

**Automated Reports and Data Trending for
Actionable Insights**

For more details on colony counting accuracy, read the comparison technical bulletin.



3M™ Petrifilm™ Rapid Aerobic Count Plates

3M™ Petrifilm™ Rapid Coliform Count Plates

3M™ Petrifilm™ Rapid E. coli/Coliform Count Plates

3M™ Petrifilm™ Rapid Yeast and Mold Count Plates

3M™ Petrifilm™ Aerobic Count Plates

3M™ Petrifilm™ Coliform Count Plates

3M™ Petrifilm™ Enterobacteriaceae Count Plates

3M™ Petrifilm™ E. coli/Coliform Count Plates

3M™ Petrifilm™ Select E. coli Count Plates

3M™ Petrifilm™ Lactic Acid Bacteria Count Plates

3M™ Petrifilm™ Staph Express Count Plates

3M™ Petrifilm™ Staph Express Disk

3M

Food Safety Division

3M Center

Bldg. 275-5W-05

St. Paul, MN 55144-1000 USA

Phone 1-800-328-6553

Web 3M.com/foodsafety/Petrifilm

3M and Petrifilm are trademarks of 3M.

Please recycle. Printed in USA.

© 3M 2022. All rights reserved. 70-2011-5232-2