ProSieve[™] ProTrack[™] Dual Color Protein Loading Buffer

Protect and Track Protein Samples


Benefits

- Protects proteins from degradation during sample preparation
- Two colors for tracking electrophoresis progress (blue) and monitoring Western transfer (pink)
- Contains SDS and DTT for complete protein denaturing





ProSieve[™] ProTrack[™] Blue Dye monitors protein separation on the gel, while the pink dye confirms transfer of the proteins onto the blot.



Lanes 1 and 3 are proteins protected by ProSieve[™] ProTrack[™] Dual Color Loading Buffer, Lanes 2 and 4 are proteins prepared and run in a standard loading buffer.

Ordering Information – ProSieve[™] ProTrack[™] Dual Color Protein Loading Buffer

Cat. No. NA	Cat. No. EU	Product Name	Application	Size
193861	193861	ProSieve™ ProTrack™ Dual Color Protein Loading Buffer (4X)	(4X), 5 μL	5 mL

9