

# The World's Most Trusted Agarose Brands



Lonza is the leading innovator and world's most trusted supplier of agarose. Our well known product brands set the standard in quality, purity and performance for electrophoresis.

- We are the inventors of agarose and have more than 40 years manufacturing experience
- Lonza offers supply chain value options such as vendor managed inventory and safety stock
- We are the only direct source of agarose in the world
- Other agarose suppliers are distributors of repackaged agarose
- Lonza quality standards are the highest in the industry
- We are the market leader with brands recognized worldwide
- We have the ability to customize agarose characteristics for specific applications and requirements
- Lonza Agarose is finely tuned by application for optimal performance





Lonza offers a range of agarose types designed by fragment size, sample type and application. A few of the well known Lonza Agarose products are included here. See www.lonza.com for a complete listing.

# **Ordering Information**

## SeaKem® LE Agarose

- The standard for routine analysis
- Wide resolution for most any application 100 bp 23 kb

Cat. No.	Size
50001	25 g
50002	100 g
50000	125 g
50004	500 g
50005	1 kg

#### SeaPlaque® Agarose

- The original low melting agarose for recovery
- Ideal for DNA and RNA recovery and cloning of tissue culture cells

Cat. No.	Size
50101	25 g
50100	125 g
50105	1 kg

## NuSieve® 3:1 Agarose

- Easy to use for small fragment PCR analysis
- Exceptional resolution of fragments 50 bp 1 kb
- Superior gel strength for easy handling and blotting

Cat. No.	Size
50091	25 g
50090	125 g
50094	500 g
50095	1 kg

#### MetaPhor® Agarose

- Resolution that rivals polyacrylamide
- Resolves 2% size differences in fragments 20 bp 800 bp

Cat. No.	Size
50181	25 g
50180	125 g
50184	500 g
50185	1 kg

Larger package sizes of Lonza Agarose are available upon request.

## **Contact Information**

## Europe

Customer Service: +32 87 321 611

Research Solutions

Scientific Support: +49 221 99199 400 scientific.support.eu@lonza.com

Research Solutions,

Therapeutic Cell Solutions, Testing Solutions
Scientific Support: +32 87 321 611
scientific.support.eu@lonza.com
Online Ordering: www.lonza.com

## International

Contact your local Lonza Distributor
Customer Service: +1 301 898 7025, ext. 2400
scientific.support@lonza.com

## International Offices

For additional international offices , please visit www.lonza.com/contactus

#### Lonza Cologne AG 50829 Cologne, Germany

For Research Use Only. Not for use in diagnostic procedures.

SeaKem® is a registered trademark of FMC.

Unless otherwise noted, all trademarks herein are marks of the Lonza Group or its affiliates.

The information contained herein is believed to be correct and corresponds to the latest state of scientific and technical knowledge. However, no warranty is made, either expressed or implied, regarding its accuracy or the results to be obtained from the use of such information and no warranty is expressed or implied concerning the use of these products. The buyer assumes all risks of use and/or handling. No statement is intended or should be construed as a recommendation to infringe any existing patent.

©Copyright 2010, Lonza Cologne AG. All rights reserved. FL-AgaroseBrands 08/10 MB-AN013

