## AccuGENE™ Protein Electrophoresis Buffers

Optimum Performance

AccuGENE™ Electrophoresis Buffers are formulated to match PAGEr™ Precast Gels. AccuGENE™ Buffers for protein electrophoresis are prepared with high quality reagents and use 18 meg0hm water. Products are filtered using a 0.2-micron filter.

- 18°C to 24°C
- www.lonza.com/sourcebook

## Benefits

- Reliable Manufactured according to strict quality control standards to ensure lot-to-lot consistency
- Efficient Ready-to-use solutions eliminate preparation time
- Flexible Customized solutions are available to meet individual needs

## Ordering Information - AccuGENE™ 10X Tris-Glycine Buffer

Cat. No. NA	Cat. No. EU	Product Name	Product Description	Storage Conditions	Size
50879	BE50879	AccuGENE™ 10X Tris-Glycine Buffer	0.25 M Tris base, 1.92 M Glycine	18°C to 24°C	1 L
50880	BE50880	AccuGENE™ 10X Tris-Glycine SDS Buffer	0.25 M Tris base, 1.92 M Glycine, 1% SDS	18°C to 24°C	1 L
50882	50882	AccuGENE™ 10X Tris-Glycine SDS Buffer	0.25 M Tris base, 1.92 M Glycine, 1% SDS	18°C to 24°C	4 L