

GelBond® Film

Agarose Support Film

GelBond® Film is a transparent, flexible polyester film designed to support agarose gels. Gels cast on GelBond® Film remain permanently attached to the film through electrophoresis or immunodiffusion and all subsequent fixing, staining, destaining, and drying procedures (gels remain flexible after drying). GelBond® Film is available either as precut sheets or rolls.

■ Benefits

- **Reliable** – Agarose gels cast on GelBond® Film retain their original dimensions during staining and after drying
- **Durable** – Gels, particularly thin ones, are easier to handle during staining, destaining, and drying when supported
- **Convenient** – Gel orientation can be recorded directly on the GelBond® Film prior to casting

NOTE: Polyester films will not transmit light of less than 310 nm, and will fluoresce at higher wavelengths.

■ Applications

- Drying and support of agarose gels

■ 18°C to 26°C

Ordering Information – GelBond® Film Sheets and Rolls

Cat. No. NA	Cat. No. EU	Product Name	Product Description	Sheet Size (mm)	Chamber Compatibility
53734	53734	GelBond® Film Sheets	For agarose gels, 0.2 mm thick	85 mm × 100 mm (100 sheets)	
53745	53745	GelBond® Film Sheets	For agarose gels, 0.2 mm thick	110 mm × 125 mm (100 sheets)	
53746	53746	GelBond® Film Sheets	For agarose gels, 0.2 mm thick	100 mm × 150 mm (100 sheets)	Bio-Rad® Wide Mini-Sub® Cell, Bio-Rad® Sub-Cell® (H)
53748	53748	GelBond® Film Sheets	For agarose gels, 0.2 mm thick	110 mm × 205 mm (100 sheets)	
53749	53749	GelBond® Film Sheets	For agarose gels, 0.2 mm thick	160 mm × 180 mm (100 sheets)	Hofer® SE400, Hofer® SE600 (V), Bio-Rad PROTEAN® II xi (V)
53759	53759	GelBond® Film Sheets	For agarose gels, 0.2 mm thick	125 mm × 245 mm (100 sheets)	
53761	53761	GelBond® Film Sheets	For agarose gels, 0.2 mm thick	124 mm × 258 mm (100 sheets)	GE Multiphor® (H)
53740	53740	GelBond® Film Rolls	For agarose gels, 0.2 mm thick	102 mm × 16.5 m (roll)	
53750	53750	GelBond® Film Rolls	For agarose gels, 0.2 mm thick	102 mm × 16.5 m (roll)	
53780	53780	GelBond® Film Rolls	For agarose gels, 0.2 mm thick	203 mm × 16.5 m (roll)	

Custom-cut GelBond® Film is available upon special request. Please inquire for pricing and availability.

(H) = Horizontal; (V) = Vertical