

TheraPEAK™ X-Vivo™ Defined, serum-free hematopoietic cell media




TheraPEAK™ X-Vivo™ chemically defined, serum-free hematopoietic cell media, now exists with recombinant transferrin.

- More regulatory friendly as a human component is replaced by a recombinant one
- Used by pharmaceutical and cell therapy customers in authorized clinical trials
- Performance is not impacted by the change from human to recombinant transferrin

TheraPEAK™ X-Vivo™ formulations provide nutritionally adequate and balanced environments for a variety of cells including:

- Lymphokine activated killer (LAK)
- Peripheral blood lymphocytes (PBL)
- Tumor infiltrating lymphocytes (TIL)
- Dendritic cells (DC)



X-Vivo™ made with recombinant transferrin is giving same level of yield than X-Vivo™ with human source

All TheraPEAK™ X-Vivo™ medium products are manufactured according to current GMP regulations. These media do not contain any exogenous growth factors, artificial stimulators of cellular proliferation, or undefined supplements. The products are devoid of any protein-kinase C stimulators and are suitable for the investigation of second messenger systems in the activation of human and murine lymphocytes. The formulations contain pharmaceutical grade human albumin, recombinant human insulin and recombinant human transferrin.

Ordering information

Cat. no.	Description	Size
TheraPEA	√™ X-Vivo™ with rec transferrin	
BE02-053Q	X-Vivo™ 15 with rec transferrin	1 L
BE02-054Q	X-Vivo™ 15 with rec transferrin w/o phenol red or gentamycin	1 L
BE02-0550	X-Vivo™ 10 with rec transferrin w/o phenol red or gentamycin	1 L
TheraPEA	√ X-Vivo™ with human transferrin	
04-3800	X-Vivo™ 10 with human transferrin	1 L
04-418F	X-Vivo™ 15 with human transferrin	0.5 L
04-4180	X-Vivo™ 15 with human transferrin	1 L
04-4480	X-Vivo™ 20 with human transferrin	1 L
04-7430	X-Vivo™ 10 with human transferrin w/o phenol red or gentamycin	1 L
04-7440	X-Vivo™ 15 with human transferrin w/o phenol red or gentamycin	1 L
Platinum l	JitraPAK™ bioprocess containers for co	ell therapy
ASCTC-1067A	Platinum UltraPAK™ bag	0.5 L
08-E1418	Platinum UltraPAK™ bag	0.25 L
08-E1769	Platinum UltraPAK™ bag	0.5 L
08-E1427	Platinum UltraPAK™ bag	2 L
08-E1322	Platinum UltraPAK™ bag	1.0 L



Contact customer service to place an order

North America

Tel 800 638 8174 order.us@lonza.com

Europe

Tel +32 87 321 611 order.europe@lonza.com

Contact information

Online ordering

https://shop.lonza.com

North America

Customer service: 800 638 8174 (toll free) order.us@lonza.com Scientific support: 800 521 0390 (toll free) scientific.support@lonza.com

Europe

Customer service: +32 87 321 611 order.europe@lonza.com Scientific support: +49 221 99199 400 scientific.support.eu@lonza.com

International

Contact your local Lonza distributor Customer service: +1 301 898 7025 scientific.support@lonza.com

International offices

+61 3 9550 0883 Δustralia 0800 201 538 (toll free) Austria +32 87 321 611 Belgium +55 11 5641 3325 Brazil Denmark 808 83 159 (toll free) France 0800 91 19 81 (toll free) Germany 0800 182 52 87 (toll free) India +91 22 4342 4000 Ireland 1 800 654 253 (toll free) Italy 800 789 888 (toll free) Japan +81 3 5566 0622 +32 87 321 611 Luxemburg +48 781 120 300 Poland +65 6521 4379 Singapore Spain 900 963 298 (toll free) Sweden 020 790 220 (toll free) Switzerland 0800 83 86 20 (toll free) The Netherlands 0800 022 4525 (toll free) United Kingdom 0808 234 97 88 (toll free)

Lonza Verviers Sprl - 4800 Verviers, Belgium

For research use only. Not for use in diagnostic procedures. Unless otherwise noted, all trademarks herein are marks of the Lonza Group or its affiliates. The information contained herein is believed to be correct and corresponds to the latest state of scientific and technical knowledge. However, no warranty is made, either expressed or implied, regarding its accuracy or the results to be obtained from the use of such information and no warranty is expressed or implied concerning the use of these products. The buyer assumes all risks of use and/or handling. No statement is intended or should be construed as a recommendation to infringe any existing patent.

©Copyright 2011, Lonza Verviers Sprl. All rights reserved. FL-ThPkXVivo 08/11 CT-TS001