


FOUNDED 1978



PRODUCT  
CATALOGUE



**EO Labs**  
Microbiological Culture Media



# EO Labs

Microbiological Culture Media

E&O Laboratories Ltd have successfully grown our Culture Media business over the last 27 years to become a leading manufacturer servicing microbiology laboratories worldwide. Our comprehensive media range is continually expanding into specialised industry sectors as well as aiding the rapid diagnosis of bacterial species exhibiting antibiotic resistance. Our professional team encompasses a wide range of skills, delivering compliant, high quality products with an exceptional level of personal service. We have strong relationships with our raw material suppliers and a proven track record for continuity of supply. E&O's two manufacturing sites are based in central Scotland with excellent road/rail and airport accessibility. Prospective customers and distributors are welcome to visit/audit and meet our loyal and dedicated team.

## Industry sector:

- Clinical, Veterinary & Aquaculture
- Food, Water & Environmental
- Pharmaceutical, Industrial & Cosmetics



## ORDER OF CONTENTS:

4 Dehydrated Culture Media

12 Supplements and Reagents

14 Donor Animal Blood and Serum

16 Ready-to-use Petri-Dishes

31 Ready-to-use Bottled and Bagged Media

46 Containers

**Manufacture:**

E&O's extensive product portfolio is manufactured across two fully automated cGMP facilities at Westfield and Burnhouse (100,000sq ft). All operations are validated to the highest pharmaceutical standards with all environmental and sterilisation processes guaranteeing minimal batch to batch variability. Pumps, balances and dispensing equipment are calibrated by highly trained technicians using industry standard procedures. Disaster strategy implementation ensures customers' continuity of supply not only with product manufacture but also with raw material supply.

**Compliance:**

Company compliance ensures all manufacturing processes and services meet the required Environmental, Health and Safety and Quality Standards regulated within the different industry sectors using Microbiological Culture Media both nationally and internationally.

- ISO 9001:2015 & ISO 13485:2003
- Animals (Scientific Procedure) Act 1986 Project Licence
- Quality control test methods accredited by UKAS to ISO 17025:2005
  - Methodology based on ISO 11133:2014 & ISO 2859-1:1999
- IVD CE Marking – European *In-Vitro* Diagnostics Medical Devices Directive 98/79/EC
- Harmonised USP/EP/JP Pharmacopoeia Methods

**Bespoke – OEM – Design:**

E&O employ a very highly qualified team of scientists to formulate new media, improve on existing formulations and work with clients to design specialised diluents/reagents and new culture media. A dedicated laboratory area with calibrated modern equipment allows a wide variety of media assays to be performed and designed before scaling up for final customer verification. E&O can source raw materials worldwide and has the in-house expertise to mill, blend, sterilise, aliquot and quality control to any specification in volumes, weights and containers to satisfy clients' individual requirements.

**Logistics & Distribution:**

E&O's logistics department dispatch a wide range of Culture Media and related products with accompanying documentation worldwide. All consignments are packed in secure, insulated (if required), strong, well designed cardboard boxes to ensure contents are not damaged during the transportation process. The quality assurance department have validated (temperature data logged) information ensuring no degradation of product occurs. All UK locations are serviced on a standard *next day* service. Timed deliveries can be booked to suit emergency requirements as well as Saturday morning deliveries.

Export consignments can be delivered via air, road or sea depending on product type and customer requirement. Having built strong supportive relationships with many distributors worldwide, our product range can be manufactured and shipped on an OEM basis or under the E&O brand.

**Account Management:**

E&O's finance department control all quotations, discounts, payment terms and any return of goods credit. Our sales department review all tender notices, distribution costs and account turnovers on an annual basis. E&O strive to deliver competitive packages and are proud of our extremely loyal and growing global customer account network. Detailed up-to-date *Terms & Conditions* can be requested through our Customer Service Department.

**Service:**

E&O's front line Customer Service team are always available and happy to direct you to the right department for all routine enquiries. Our technical team support and advise on all microbiological and scientific product dialogue. E&O's product portfolio also includes innovative service advantages from integrated stock management systems, to special delivery schedules using returnable packaging and enhanced Quality Control reporting online 24/7.



| Product Code | Description                        | Shelf Life (days) | Storage Temp (°C) | Related Supplement |  | Related Supplement Shelf Life (days) | Related Supplement Storage Temp (°C) |
|--------------|------------------------------------|-------------------|-------------------|--------------------|--|--------------------------------------|--------------------------------------|
| KM8010       | Acid Casein                        | 1095              | 10-30°C           |                    |  |                                      |                                      |
| KM0075       | Alkaline Peptone Water             | 1095              | 10-30°C           |                    |  |                                      |                                      |
| KM0196       | Artificial Sea Water               | 1095              | 10-30°C           |                    |  |                                      |                                      |
| KM0186       | Bacillus cereus Agar (PEMBA)       | 1095              | 10-30°C           | LS1051             | Bacillus cereus Selective Supplement (PEMBA) | 730                                  | 2-8°C                                |
|              |                                    |                   |                   | BM0140             | Egg Yolk Emulsion                            | 182                                  | 2-8°C                                |
| KM0084       | Bacillus cereus Agar (MYP)         | 1095              | 10-30°C           | LS0020             | Bacillus cereus Selective Supplement         | 730                                  | 2-8°C                                |
|              |                                    |                   |                   | BM0140             | Egg Yolk Emulsion                            | 182                                  | 2-8°C                                |
| KM9001       | Bacteriological Agar               | 1095              | 10-30°C           |                    |  |                                      |                                      |
| KM8001       | Bacteriological Peptone            | 1095              | 10-30°C           |                    |  |                                      |                                      |
| KM0037       | Baird Parker Agar                  | 1095              | 10-30°C           | BM0530             | Egg Yolk Tellurite                           | 182                                  | 2-8°C                                |
| KM8006       | Beef Extract                       | 1095              | 10-30°C           |                    |  |                                      |                                      |
| KM0085       | Bile Aesculin Agar                 | 1095              | 10-30°C           |                    |  |                                      |                                      |
| KM0034       | Blood Agar Base No.2               | 1095              | 10-30°C           | DHB –              | Defibrinated Horse Blood                     | 56                                   | 2-8°C                                |
|              |                                    |                   |                   | DSC –              | Defibrinated Sheep Blood                     | 56                                   | 2-8°C                                |
| KM0091       | Brain Heart Infusion (BHI) Agar    | 1095              | 10-30°C           | DHB –              | Defibrinated Horse Blood                     | 56                                   | 2-8°C                                |
| KM0092       | Brain Heart Infusion (BHI) Broth   | 1095              | 10-30°C           |                    |  |                                      |                                      |
| KM8009       | Brain Heart Infusion Solids (BHI)  | 1095              | 10-30°C           |                    |  |                                      |                                      |
| KM0038       | Brilliant Green Agar               | 1095              | 10-30°C           |                    |  |                                      |                                      |
| KM0042       | Buffered Listeria Enrichment Broth | 1095              | 10-30°C           | LS5006             | Buffered Listeria Enrichment Supplement      | 730                                  | 2-8°C                                |
| KM0039       | Buffered Peptone Water (ISO)       | 1095              | 10-30°C           |                    |  |                                      |                                      |

| Product Code | Description                          | Shelf Life (days) | Storage Temp (°C) | Related Supplement                                 | Related Supplement Shelf Life (days) | Related Supplement Storage Temp (°C) |
|--------------|--------------------------------------|-------------------|-------------------|--|--------------------------------------|--------------------------------------|
| KM0076       | Burkholderia cepacia Agar Base       | 1095              | 10-30°C           | LS0125 B.cepacia Selective Supplement              | 730                                  | 2-8°C                                |
|              |                                      |                   |                   | LS0026 Pseudomonas CFC Selective Supplement        | 730                                  | 2-8°C                                |
| KM0005       | CLED Double Indicator Agar           | 1095              | 10-30°C           |  |                                      |                                      |
| KM0004       | CLED Single Indicator Agar           | 1095              | 10-30°C           |  |                                      |                                      |
| KM0044       | Campylobacter Agar Base              | 1095              | 10-30°C           | LS0010 Campylobacter (CCDA) Selective Supplement   | 730                                  | 2-8°C                                |
|              |                                      |                   |                   | HBL – Horse Blood Lysed                            | 182                                  | 2-8°C                                |
| KM0057       | Campylobacter Blood Free Agar (CCDA) | 1095              | 10-30°C           | LS0010 Campylobacter (CCDA) Selective Supplement   | 730                                  | 2-8°C                                |
| KM0045       | Campylobacter Enrichment Broth       | 1095              | 10-30°C           | LS1054 Campylobacter (Bolton) Selective Supplement | 730                                  | 2-8°C                                |
|              |                                      |                   |                   | HBL – Horse Blood Lysed                            | 182                                  | 2-8°C                                |
| KM0052       | CEMO Agar                            | 1095              | 10-30°C           | LS0041 CEMO Selective Supplement                   | 730                                  | 2-8°C                                |
|              |                                      |                   |                   | DHB – Defibrinated Horse Blood                     | 56                                   | 2-8°C                                |
| KM0059       | Cetrimide Agar                       | 1095              | 10-30°C           |  |                                      |                                      |
| KM0185       | Charcoal Agar                        | 1095              | 10-30°C           | LS0018 Bordetella Selective Supplement             | 730                                  | 2-8°C                                |
|              |                                      |                   |                   | DHB – Defibrinated Horse Blood                     | 56                                   | 2-8°C                                |
| KM0001       | Columbia Agar                        | 1095              | 10-30°C           | DHB – Defibrinated Horse Blood                     | 56                                   | 2-8°C                                |
|              |                                      |                   |                   | DSC – Defibrinated Sheep Blood                     | 56                                   | 2-8°C                                |
|              |                                      |                   |                   | HBL – Horse Blood Lysed                            | 182                                  | 2-8°C                                |
|              |                                      |                   |                   | LS0008 Staph/Strep Selective Supplement            | 730                                  | 2-8°C                                |
|              |                                      |                   |                   | LS0017 Neomycin Selective Supplement               | 730                                  | 2-8°C                                |
|              |                                      |                   |                   | LS0022 Clostridium difficile Selective Supplement  | 730                                  | 2-8°C                                |

| Product Code | Description                                       | Shelf Life (days) | Storage Temp (°C) | Related Supplement | Related Supplement Shelf Life (days)           | Related Supplement Storage Temp (°C) |
|--------------|---|-------------------|-------------------|--------------------|--|--------------------------------------|
|              |   |                   |                   | LS0012             | Bacitracin Selective Supplement                | 730<br>2-8°C                         |
|              |   |                   |                   | LS0011             | Campylobacter Growth Supplement                | 730<br>2-8°C                         |
|              |   |                   |                   | LS0009             | Campylobacter (Skirrow) Selective Supplement   | 730<br>2-8°C                         |
| KM0050       | Desoxycholate Citrate Agar (Hynes)                | 1095              | 10-30°C           |                    |  |                                      |
| KM0184       | Dermatophyte Test Medium                          | 1095              | 10-30°C           | LS0050             | Chloramphenicol Selective Supplement (100mg/L) | 730<br>2-8°C                         |
| KM0078       | Dichloran Rose-Bengal Chloramphenicol (DRBC) Agar | 1095              | 10-30°C           |                    |  |                                      |
| KM0081       | Dichloran-Glycerol (DG18) Agar                    | 1095              | 10-30°C           |                    |  |                                      |
| KM0053       | DNase Agar  | 1095              | 10-30°C           |                    |  |                                      |
| KM0066       | DNase Test Agar with Methyl Green                 | 1095              | 10-30°C           |                    |  |                                      |
| KM0026       | Edwardsiella ictaluri Medium                      | 1095              | 10-30°C           |                    |  |                                      |
| KM0104       | EE Broth (Mossel)                                 | 1095              | 10-30°C           |                    |  |                                      |
| KM0105       | Endo Agar Base                                    | 1095              | 10-30°C           |                    |  |                                      |
| KM0028       | Eosin Methylene Blue (EMB) Agar                   | 1095              | 10-30°C           |                    |  |                                      |
| KM0006       | Fastidious Anaerobe Agar                          | 1095              | 10-30°C           | LS0015             | Actinomycete Selective Supplement              | 730<br>2-8°C                         |
|              |   |                   |                   | LS0017             | Neomycin Selective Supplement                  | 730<br>2-8°C                         |
|              |   |                   |                   | LS0022             | Clostridium difficile Selective Supplement     | 730<br>2-8°C                         |
|              |   |                   |                   | DHB –              | Defibrinated Horse Blood                       | 56<br>2-8°C                          |
|              |   |                   |                   | DSC –              | Defibrinated Sheep Blood                       | 56<br>2-8°C                          |
|              |   |                   |                   | BM0140             | Egg Yolk Emulsion                              | 182<br>2-8°C                         |
| KM0188       | Fastidious Anaerobe Broth                         | 1095              | 10-30°C           |                    |  |                                      |

| Product Code | Description                        | Shelf Life (days) | Storage Temp (°C) | Related Supplement | Related Supplement Shelf Life (days)     | Related Supplement Storage Temp (°C) |       |
|--------------|------------------------------------|-------------------|-------------------|--------------------|--|--------------------------------------|-------|
| KM0165       | Fluid Thioglycollate Medium (USP)  | 1095              | 10-30°C           |                    |  |                                      |       |
| KM0046       | Fraser Broth – Half Strength       | 1095              | 10-30°C           | LS5003             | Half Fraser Broth Selective Supplement   | 730                                  | 2-8°C |
|              |                                    |                   |                   | LS5004             | Ferric Ammonium Citrate (FAC) Supplement | 730                                  | 2-8°C |
| KM0016       | GC Agar                            | 1095              | 10-30°C           | LS0001             | GC LCAT Selective Supplement             | 730                                  | 2-8°C |
|              |                                    |                   |                   | LS0002             | GC VCAT Selective Supplement             | 730                                  | 2-8°C |
|              |                                    |                   |                   | BM0478             | Suplex Solution                          | 182                                  | 2-8°C |
|              |                                    |                   |                   | DHB –              | Defibrinated Horse Blood                 | 56                                   | 2-8°C |
| KM8013       | Gelatin Peptone                    | 1095              | 10-30°C           |                    |  |                                      |       |
| KM0032       | Hektoen Enteric Agar               | 1095              | 10-30°C           |                    |  |                                      |       |
| KM0009       | Helicobacter pylori Agar           | 1095              | 10-30°C           | LS0031             | Helicobacter pylori Selective Supplement | 730                                  | 2-8°C |
|              |                                    |                   |                   | SHS –              | Sterile Horse Serum                      | 1095                                 | -20°C |
| KM0073       | Hoyle's Agar                       | 1095              | 10-30°C           | BM2230             | 3.5% Potassium Tellurite Solution        | 365                                  | 2-8°C |
|              |                                    |                   |                   | HBL –              | Horse Blood Lysed                        | 182                                  | 2-8°C |
| KM0111       | Kligler Iron Agar                  | 1095              | 10-30°C           |                    |  |                                      |       |
| KM8008       | Lactalbumin hydrolysate            | 1095              | 10-30°C           |                    |  |                                      |       |
| KM0116       | Lactose Broth                      | 1095              | 10-30°C           |                    |  |                                      |       |
| KM0015       | Legionella Agar                    | 1095              | 10-30°C           | LS0052             | Legionella GVPC Selective Supplement     | 730                                  | 2-8°C |
|              |                                    |                   |                   | LS0053             | Legionella BCYE Growth Supplement        | 730                                  | 2-8°C |
|              |                                    |                   |                   | LS0054             | Legionella BCYE without Cysteine         | 730                                  | 2-8°C |
| KM0121       | Listeria Isolation Medium (Oxford) | 1095              | 10-30°C           | LS0030             | Listeria Oxford Selective Supplement     | 730                                  | 2-8°C |

| Product Code | Description                                  | Shelf Life (days) | Storage Temp (°C) | Related Supplement                               | Related Supplement Shelf Life (days) | Related Supplement Storage Temp (°C) |
|--------------|--|-------------------|-------------------|--|--------------------------------------|--------------------------------------|
| KM0194       | Listeria Enhanced Enrichment Broth           | 1095              | 10-30°C           | LS0260 Listeria Enhanced Selective Supplement    | 730                                  | 2-8°C                                |
| KM8015       | Liver Peptone                                | 1095              | 10-30°C           |  |                                      |                                      |
| KM0060       | Lowenstein-Jensen Medium                     | 1095              | 10-30°C           |  |                                      |                                      |
| KM0115       | LB Agar (Lennox)                             | 1095              | 10-30°C           | LS0035 Ampicillin Selective Supplement (100mg/L) | 730                                  | 2-8°C                                |
|              |  |                   |                   | LS0037 Kanamycin Selective Supplement (50mg/L)   | 730                                  | 2-8°C                                |
| KM0114       | LB Agar (Miller)                             | 1095              | 10-30°C           | LS0035 Ampicillin Selective Supplement (100mg/L) | 730                                  | 2-8°C                                |
|              |  |                   |                   | LS0037 Kanamycin Selective Supplement (50mg/L)   | 730                                  | 2-8°C                                |
| KM0113       | LB Broth (Lennox)                            | 1095              | 10-30°C           |  |                                      |                                      |
| KM0112       | LB Broth (Miller)                            | 1095              | 10-30°C           |  |                                      |                                      |
| KM0122       | Lysine Iron Agar                             | 1095              | 10-30°C           |  |                                      |                                      |
| KM0080       | MRS Agar (de Man, Rogosa & Sharpe Agar)      | 1095              | 10-30°C           |  |                                      |                                      |
| KM0027       | MacConkey Agar No.3                          | 1095              | 10-30°C           | LS0189 Cefotaxime Supplement (1mg/L)             | 730                                  | 2-8°C                                |
| KM0124       | MacConkey Agar with Salt                     | 1095              | 10-30°C           |  |                                      |                                      |
| KM0011       | MacConkey Agar without Salt & Crystal Violet | 1095              | 10-30°C           |  |                                      |                                      |
| KM0125       | MacConkey Broth (USP)                        | 1095              | 10-30°C           |  |                                      |                                      |
| KM8005       | Malt Extract                                 | 1095              | 10-30°C           |  |                                      |                                      |
| KM0048       | Malt Extract Agar                            | 1095              | 10-30°C           | LS0050 Chloramphenicol Selective Supplement      | 730                                  | 2-8°C                                |
| KM0128       | Malt Extract Broth                           | 1095              | 10-30°C           |  |                                      |                                      |
| KM0031       | Mannitol Salt Agar                           | 1095              | 10-30°C           |  |                                      |                                      |
| KM0062       | Marine Agar                                  | 1095              | 10-30°C           |  |                                      |                                      |



| Product Code | Description   | Shelf Life (days) | Storage Temp (°C) | Related Supplement  | Related Supplement Shelf Life (days) | Related Supplement Storage Temp (°C) |
|--------------|---|-------------------|-------------------|---|--------------------------------------|--------------------------------------|
| KM0195       | Marine Salt Agar  | 1095              | 10-30°C           |   |                                      |                                      |
| KM0040       | Maximum Recovery Diluent (MRD)                          | 1095              | 10-30°C           |   |                                      |                                      |
| KM8014       | Meat Extract  | 1095              | 10-30°C           |   |                                      |                                      |
| KM8004       | Meat Peptone  | 1095              | 10-30°C           |   |                                      |                                      |
| KM0183       | Mueller-Hinton Agar                                     | 1095              | 10-30°C           | LS0005 β-Nicotinamide-Adenine Dinucleotide (NAD) (20mg/L) | 730                                  | 2-8°C                                |
|              |   |                   |                   | DHB – Defibrinated Horse Blood                            | 56                                   | 2-8°C                                |
|              |   |                   |                   | DSC – Defibrinated Sheep Blood                            | 56                                   | 2-8°C                                |
| KM0063       | Muller-Kauffmann Tetrathionate-Novobiocin Broth (MKTTn) | 1095              | 10-30°C           | LS0024 Novobiocin Supplement (20mg/L)                     | 730                                  | 2-8°C                                |
|              |   |                   |                   | BM0946 Iodine & Potassium Iodide Complex                  | 182                                  | 15-30°C                              |
| KM0055       | Mycological Agar  | 1095              | 10-30°C           |   |                                      |                                      |
| KM0141       | Nutrient Agar   | 1095              | 10-30°C           |   |                                      |                                      |
| KM0142       | Nutrient Broth  | 1095              | 10-30°C           |   |                                      |                                      |
| KM0041       | Nutrient Broth No.2                                     | 1095              | 10-30°C           |   |                                      |                                      |
| KM0146       | PALCAM Agar Base  | 1095              | 10-30°C           | LS0038 PALCAM Selective Supplement                        | 730                                  | 2-8°C                                |
| KM0149       | Peptone Water   | 1095              | 10-30°C           |   |                                      |                                      |
| KM0051       | Phosphate Buffered Saline (PBS)                         | 1095              | 10-30°C           |   |                                      |                                      |
| KM0064       | Plate Count Agar  | 1095              | 10-30°C           |   |                                      |                                      |
| KM0056       | Primary Chromogenic Coliform Agar (CCA)                 | 1095              | 10-30°C           |   |                                      |                                      |
| KM0058       | Primary Listeria Agar                                   | 1095              | 10-30°C           | LS0210 Primary Listeria Selective Supplement              | 730                                  | 2-8°C                                |
|              |   |                   |                   | BM8050 Primary Listeria Diagnostic Supplement             | 182                                  | 15-30°C                              |

| Product Code | Description                                      | Shelf Life (days) | Storage Temp (°C) | Related Supplement |                                      | Related Supplement Shelf Life (days) | Related Supplement Storage Temp (°C) |
|--------------|--|-------------------|-------------------|--------------------|--------------------------------------|--------------------------------------|--------------------------------------|
| KM0187       | Primary membrane Lactose Glucuronide Agar (mLGA) | 1095              | 10-30°C           |                    |                                      |                                      |                                      |
| KM0162       | Primary Tryptone Bile X-Glucuronide (TBX) Agar   | 1095              | 10-30°C           |                    |                                      |                                      |                                      |
| KM0007       | Primary UTI Agar                                 | 1095              | 10-30°C           | SHS –              | Sterile Horse Serum                  | 1095                                 | -20°C                                |
|              |  |                   |                   |                    | KM8020 Kaolin                        | 730                                  | 10-30°C                              |
| KM0079       | Pseudomonas Agar Base                            | 1095              | 10-30°C           | LS0006             | Pseudomonas CN Supplement            | 730                                  | 2-8°C                                |
|              |  |                   |                   | LS0026             | Pseudomonas CFC Selective Supplement | 730                                  | 2-8°C                                |
| KM0082       | R2A Agar   | 1095              | 10-30°C           |                    |                                      |                                      |                                      |
| KM0047       | Rappaport-Vassiliadis Broth                      | 1095              | 10-30°C           |                    |                                      |                                      |                                      |
| KM0193       | Ringers ¼ Strength                               | 1095              | 10-30°C           |                    |                                      |                                      |                                      |
| KM0003       | Sabouraud Dextrose Agar                          | 1095              | 10-30°C           |                    |                                      |                                      |                                      |
| KM0068       | Sabouraud Dextrose Agar with Chloramphenicol     | 1095              | 10-30°C           |                    |                                      |                                      |                                      |
| KM0154       | Sabouraud Liquid Medium (USP)                    | 1095              | 10-30°C           |                    |                                      |                                      |                                      |
| KM0033       | Salmonella Shigella (SS) Agar Modified           | 1095              | 10-30°C           |                    |                                      |                                      |                                      |
| KM0158       | Selenite Cystine Broth Base                      | 1095              | 10-30°C           |                    |                                      |                                      |                                      |
| KM0074       | Selenite F Broth Base                            | 1095              | 10-30°C           |                    |                                      |                                      |                                      |
| KM0022       | Selenite Mannitol Broth Base                     | 1095              | 10-30°C           |                    |                                      |                                      |                                      |
| KM0189       | Sheep Blood Agar Base                            | 1095              | 10-30°C           | LS0008             | Staph/Strep Selective Supplement     | 730                                  | 2-8°C                                |
|              |  |                   |                   |                    | DSC – Defibrinated Sheep Blood       | 56                                   | 2-8°C                                |
| KM0160       | Simmons Citrate Agar                             | 1095              | 10-30°C           |                    |                                      |                                      |                                      |
| KM0161       | Slanetz & Bartley Medium                         | 1095              | 10-30°C           |                    |                                      |                                      |                                      |

| Product Code | Description   | Shelf Life (days) | Storage Temp (°C) | Related Supplement                                    | Related Supplement Shelf Life (days) | Related Supplement Storage Temp (°C) |
|--------------|---|-------------------|-------------------|---|--------------------------------------|--------------------------------------|
| KM0010       | Sorbitol MacConkey Agar                             | 1095              | 10-30°C           | LS0013 Escherichia coli O157 Selective Supplement     | 730                                  | 2-8°C                                |
| KM8007       | Soy Peptone   | 1095              | 10-30°C           |   |                                      |                                      |
| KM0008       | Thiosulphate Citrate Bile Salts Sucrose (TCBS) Agar | 1095              | 10-30°C           |   |                                      |                                      |
| KM0070       | Tryptone Yeast Extract Salt (TYES) Agar             | 1095              | 10-30°C           |   |                                      |                                      |
| KM0168       | Triple Sugar Iron Agar                              | 1095              | 10-30°C           |   |                                      |                                      |
| KM8002       | Tryptone  | 1095              | 10-30°C           |   |                                      |                                      |
| KM0024       | Tryptone Soya Agar                                  | 1095              | 10-30°C           | DHB – Defibrinated Horse Blood                        | 56                                   | 2-8°C                                |
|              |   |                   |                   | DSC – Defibrinated Sheep Blood                        | 56                                   | 2-8°C                                |
| KM0025       | Tryptone Soya Broth                                 | 1095              | 10-30°C           |   |                                      |                                      |
| KM0029       | Tryptone Water                                      | 1095              | 10-30°C           |   |                                      |                                      |
| KM8011       | Tryptose  | 1095              | 10-30°C           |   |                                      |                                      |
| KM0172       | Urea Agar Base                                      | 1095              | 10-30°C           | BM3000 40% Urea Solution                              | 90                                   | 2-8°C                                |
| KM0173       | Urea Broth Base                                     | 1095              | 10-30°C           | BM3000 40% Urea Solution                              | 90                                   | 2-8°C                                |
| KM0065       | Violet Red Bile Agar (VRBA)                         | 1095              | 10-30°C           |   |                                      |                                      |
| KM0067       | Violet Red Bile Glucose Agar (VRBGA)                | 1095              | 10-30°C           |   |                                      |                                      |
| KM0013       | Xylose Lysine Desoxycholate Agar (XLD)              | 1095              | 10-30°C           |   |                                      |                                      |
| KM8003       | Yeast Extract                                       | 1095              | 10-30°C           |   |                                      |                                      |
| KM0177       | Yeast Extract Agar                                  | 1095              | 10-30°C           |   |                                      |                                      |
| KM0019       | Yeast Extract Glucose Chloramphenicol Agar          | 1095              | 10-30°C           | LS0019 Oxytetracycline Selective Supplement (100mg/L) | 730                                  | 2-8°C                                |

| Product Code | Container | Volume          | Description  | Pack Size       | Shelf Life (days) | Storage Temp (°C) |
|--------------|-----------|-----------------|--|-----------------|-------------------|-------------------|
| LS0015       | V001      | 500ml / 1 Litre | Actinomycete Selective Supplement                    | Box of 10 Vials | 730               | 2-8°C             |
| LS0035       | V001      | 500ml / 1 Litre | Ampicillin Selective Supplement (100mg/L)            | Box of 10 Vials | 730               | 2-8°C             |
| LS0125       | V001      | 500ml / 1 Litre | B.cepacia Selective Supplement                       | Box of 10 Vials | 730               | 2-8°C             |
| LS1051       | V001      | 500ml / 1 Litre | Bacillus cereus PEMBA Supplement                     | Box of 10 Vials | 730               | 2-8°C             |
| LS0020       | V001      | 500ml / 1 Litre | Bacillus cereus Selective Supplement                 | Box of 10 Vials | 730               | 2-8°C             |
| LS0012       | V001      | 500ml / 1 Litre | Bacitracin Selective Supplement                      | Box of 10 Vials | 730               | 2-8°C             |
| LS0005       | V001      | 500ml / 1 Litre | β-Nicotinamide-Adenine Dinucleotide (NAD) Supplement | Box of 10 Vials | 730               | 2-8°C             |
| LS0018       | V001      | 500ml / 1 Litre | Bordetella Selective Supplement                      | Box of 10 Vials | 730               | 2-8°C             |
| LS5006       | V001      | 500ml / 1 Litre | Buffered Listeria Enrichment Supplement              | Box of 10 Vials | 730               | 2-8°C             |
| LS0041       | V001      | 500ml / 1 Litre | CEMO Selective Supplement                            | Box of 10 Vials | 730               | 2-8°C             |
| LS0010       | V001      | 500ml / 1 Litre | Campylobacter (CCDA) Selective Supplement            | Box of 10 Vials | 730               | 2-8°C             |
| LS1054       | V001      | 500ml / 1 Litre | Campylobacter (Bolton) Selective Supplement          | Box of 10 Vials | 730               | 2-8°C             |
| LS0009       | V001      | 500ml / 1 Litre | Campylobacter (Skirrow) Selective Supplement         | Box of 10 Vials | 730               | 2-8°C             |
| LS0011       | V001      | 500ml / 1 Litre | Campylobacter Growth Supplement                      | Box of 10 Vials | 730               | 2-8°C             |
| LS0189       | V001      | 500ml / 1 Litre | Cefotaxime Selective Supplement (1mg/L)              | Box of 10 Vials | 730               | 2-8°C             |
| LS0050       | V001      | 500ml / 1 Litre | Chloramphenicol Selective Supplement (100mg/L)       | Box of 10 Vials | 730               | 2-8°C             |
| LS0022       | V001      | 500ml / 1 Litre | Clostridium difficile Selective Supplement           | Box of 10 Vials | 730               | 2-8°C             |
| LS0025       | V001      | 500ml           | Clostridium perfringens (SFP) Selective Supplement   | Box of 10 Vials | 730               | 2-8°C             |
| LS0023       | V001      | 500ml / 1 Litre | Clostridium perfringens (TSC) Selective Supplement   | Box of 10 Vials | 730               | 2-8°C             |
| LS0013       | V001      | 500ml / 1 Litre | Escherichia coli O157 Selective Supplement           | Box of 10 Vials | 730               | 2-8°C             |
| LS5004       | V001      | 500ml / 1 Litre | Ferric Ammonium Citrate (FAC) Supplement             | Box of 10 Vials | 730               | 2-8°C             |
| LS5002       | V001      | 500ml / 1 Litre | Fraser Broth Selective Supplement                    | Box of 10 Vials | 730               | 2-8°C             |
| LS0002       | V001      | 500ml / 1 Litre | GC VCAT Selective Supplement                         | Box of 10 Vials | 730               | 2-8°C             |
| LS0001       | V001      | 500ml / 1 Litre | GC LCAT Selective Supplement                         | Box of 10 Vials | 730               | 2-8°C             |

| Product Code | Container     | Volume                  | Description                                    | Pack Size   | Shelf Life (days) | Storage Temp (°C) |
|--------------|---------------|-------------------------|--|---|-------------------|-------------------|
| LS5003       | V001          | 225ml/500ml/<br>1 Litre | Half Fraser Broth Selective Supplement         | Box of 10 Vials   | 730               | 2-8°C             |
| LS0031       | V001          | 500ml/1 Litre           | Helicobacter pylori Selective Supplement       | Box of 10 Vials   | 730               | 2-8°C             |
| LS0037       | V001          | 500ml/1 Litre           | Kanamycin Selective Supplement (50mg/L)        | Box of 10 Vials   | 730               | 2-8°C             |
| LS0052       | V001          | 500ml                   | Legionella GVPC Selective Supplement           | Box of 10 Vials   | 730               | 2-8°C             |
| LS0053       | V001          | 500ml                   | Legionella BCYE Growth Supplement              | Box of 10 Vials   | 730               | 2-8°C             |
| LS0054       | V001          | 500ml                   | Legionella BCYE without Cysteine Supplement    | Box of 10 Vials   | 730               | 2-8°C             |
| LS0030       | V001          | 500ml/1 Litre           | Listeria Oxford Selective Supplement           | Box of 10 Vials   | 730               | 2-8°C             |
| LS0260       | V001          | 500ml/1 Litre           | Listeria Enhanced Selective Supplement         | Box of 10 Vials   | 730               | 2-8°C             |
| LS0017       | V001          | 500ml/1 Litre           | Neomycin Selective Supplement                  | Box of 10 Vials   | 730               | 2-8°C             |
| LS0024       | V001          | 500ml/1 Litre           | Novobiocin Supplement (20mg/L)                 | Box of 10 Vials   | 730               | 2-8°C             |
| LS0019       | V001          | 500ml/1 Litre           | Oxytetracycline Selective Supplement (100mg/L) | Box of 10 Vials   | 730               | 2-8°C             |
| LS0038       | V001          | 500ml/1 Litre           | PALCAM Selective Supplement                    | Box of 10 Vials   | 730               | 2-8°C             |
| LS0210       | V001          | 500ml/1 Litre           | Primary Listeria Selective Supplement          | Box of 10 Vials   | 730               | 2-8°C             |
| LS0006       | V001          | 500ml/1 Litre           | Pseudomonas CN Supplement                      | Box of 10 Vials   | 730               | 2-8°C             |
| LS0026       | V001          | 500ml/1 Litre           | Pseudomonas CFC Selective Supplement           | Box of 10 Vials   | 730               | 2-8°C             |
| LS0008       | V001          | 500ml/1 Litre           | Staph/Strep Selective Supplement               | Box of 10 Vials   | 730               | 2-8°C             |
| BM0140       | N100;<br>N500 | 100ml; 500ml            | Egg Yolk Emulsion                              | 1 x 100ml Square PETG Bottle;<br>1 x 500ml Square PETG Bottle | 182               | 2-8°C             |
| BM0530       | N100;<br>N500 | 100ml; 500ml            | Egg Yolk Tellurite                             | 1 x 100ml Square PETG Bottle;<br>1 x 500ml Square PETG Bottle | 182               | 2-8°C             |
| BM0946       | AM07          | 10ml                    | Iodine & Potassium Iodide Complex              | 50 x Moulded Wall Universal                                   | 182               | 15-30°C           |
| BM2230       | R125          | 100ml                   | Potassium Tellurite Solution (3.5% Solution)   | 1 x 125ml Syrup Bottle  | 365               | 2-8°C             |
| BM8050       | R125          | 40ml                    | Primary Listeria Diagnostic Supplement         | 10 x 125ml Syrup Bottle                                       | 182               | 15-30°C           |
| BM0478       | R125;<br>R500 | 100ml; 500ml            | Suplex Solution                                | 10 x 125ml Syrup Bottle; 6 x 500ml Syrup Bottle               | 182               | 2-8°C             |
| BM3000       | R125          | 100ml                   | 40% Urea Solution                              | 1 x 125ml Syrup Bottle  | 90                | 2-8 °C            |
| LRP020       | X001          | 20ml                    | Lyophilised Rabbit Plasma 20ml                 | 1 x Glass Vial  | 365               | 2-8 °C            |

| Product Code | Container  | Volume              | Description                        | Pack Size  | Shelf Life (days) | Storage Temp (°C) |
|--------------|------------|---------------------|------------------------------------|--|-------------------|-------------------|
| DHB025       | P030       | 25ml                | Defibrinated Horse Blood 25ml      | 1 x Polystyrene Universal                            | 56                | 2-8°C             |
| DHB050       | N100; Y125 | 50ml; 50ml          | Defibrinated Horse Blood 50ml      | 1 x 100ml Square PETG Bottle; 1 x 125ml Syrup Bottle | 56                | 2-8°C             |
| DHB100       | N100; Y125 | 100ml; 100ml        | Defibrinated Horse Blood 100ml     | 1 x 100ml Square PETG Bottle; 1 x 125ml Syrup Bottle | 56                | 2-8°C             |
| DHB300       | Y500       | 300ml               | Defibrinated Horse Blood 300ml     | 1 x 500ml Syrup Bottle                               | 56                | 2-8°C             |
| DHB500       | N500; Y500 | 500ml; 500ml        | Defibrinated Horse Blood 500ml     | 1 x 500ml Square PETG Bottle; 1 x 500ml Syrup Bottle | 56                | 2-8°C             |
| DHB1L        | 1LB; Y1L   | 1,000ml;<br>1,000ml | Defibrinated Horse Blood 1Lt       | 1 x 1 Litre Blood Bag; 1 x 1000ml Syrup Bottle       | 56                | 2-8°C             |
| DHB2L        | 4LB        | 2,000ml             | Defibrinated Horse Blood 2Lt       | 1 x 4 Litre Blood Bag                                | 56                | 2-8°C             |
| DHB4L        | 4LB        | 4,000ml             | Defibrinated Horse Blood 4Lt       | 1 x 4 Litre Blood Bag                                | 56                | 2-8°C             |
| DHB10L       | 10LB       | 10,000ml            | Defibrinated Horse Blood 10Lt      | 1 x 10 Litre Blood Bag                               | 56                | 2-8°C             |
| HIS100       | N100       | 100ml               | Heat Inactivated Horse Serum 100ml | 1 x 100ml Square PETG Bottle                         | 1095              | -20°C             |
| HIS500       | N500       | 500ml               | Heat Inactivated Horse Serum 500ml | 1 x 500ml Square PETG Bottle                         | 1095              | -20°C             |
| HBL025       | P030       | 25ml                | Horse Blood Lysed 25ml             | 1 x Polystyrene Universal                            | 182               | 2-8°C             |
| HBL100       | N100; Y125 | 100ml; 100ml        | Horse Blood Lysed 100ml            | 1 x 100ml Square PETG Bottle; 1 x 125ml Syrup Bottle | 182               | 2-8°C             |
| HBL500       | N500       | 500ml               | Horse Blood Lysed 500ml            | 1 x 500ml Square PETG Bottle                         | 182               | 2-8°C             |
| SHS100       | N100       | 100ml               | Sterile Horse Serum 100ml          | 1 x 100ml Square PETG Bottle                         | 1095              | -20°C             |

| Product Code | Container  | Volume       | Description                    | Pack Size   | Shelf Life (days) | Storage Temp (°C) |
|--------------|------------|--------------|--------------------------------|---|-------------------|-------------------|
| SHS500       | N500       | 500ml        | Sterile Horse Serum 500ml      | 1 x 500ml Square PETG Bottle                            | 1095              | -20°C             |
| DSC010       | P030       | 10ml         | Defibrinated Sheep Blood 10ml  | 1 x Polystyrene Universal                               | 56                | 2-8°C             |
| DSC025       | P030       | 25ml         | Defibrinated Sheep Blood 25ml  | 1 x Polystyrene Universal                               | 56                | 2-8°C             |
| DSC050       | N100; Y125 | 50ml; 50ml   | Defibrinated Sheep Blood 50ml  | 1 x 100ml Square PETG Bottle;<br>1 x 125ml Syrup Bottle | 56                | 2-8°C             |
| DSC100       | N100; Y125 | 100ml; 100ml | Defibrinated Sheep Blood 100ml | 1 x 100ml Square PETG Bottle; 1 x 125ml Syrup Bottle    | 56                | 2-8°C             |
| DSC200       | Y300       | 200ml        | Defibrinated Sheep Blood 200ml | 1 x 300ml Syrup Bottle                                  | 56                | 2-8°C             |
| DSC500       | N500       | 500ml        | Defibrinated Sheep Blood 500ml | 1 x 500ml Square PETG Bottle                            | 56                | 2-8°C             |
| DSC1L        | 1LB        | 1,000ml      | Defibrinated Sheep Blood 1Lt   | 1 x 1 Litre Blood Bag                                   | 56                | 2-8°C             |
| DSC3L        | 4LB        | 3,000ml      | Defibrinated Sheep Blood 3Lt   | 1 x 4 Litre Blood Bag                                   | 56                | 2-8°C             |
| DSC10L       | 10LB       | 10,000ml     | Defibrinated Sheep Blood 10Lt  | 1 x 10 Litre Blood Bag                                  | 56                | 2-8°C             |
| SBC025       | P030       | 25ml         | Sheep Blood Citrate 25ml       | 1 x Polystyrene Universal                               | 28                | 2-8°C             |
| SBC100       | Y125       | 100ml        | Sheep Blood Citrate 100ml      | 1 x 125ml Syrup Bottle                                  | 28                | 2-8°C             |
| SCA025       | P030       | 25ml         | Sheep Cells in Alsevers 25ml   | 1 x Polystyrene Universal                               | 28                | 2-8°C             |
| SCA050       | Y125       | 50ml         | Sheep Cells in Alsevers 50ml   | 1 x 125ml Syrup Bottle                                  | 28                | 2-8°C             |
| SCA100       | N100       | 100ml        | Sheep Cells in Alsevers 100ml  | 1 x 100ml Square PETG Bottle                            | 28                | 2-8°C             |

| Product Code                                  | Container | Volume  | Description   | Pack Size        | Packaging | Shelf Life (days) | Storage Temp (°C) |
|---|-----------|---------|---|------------------|-----------|-------------------|-------------------|
| <b>Clinical, Veterinary &amp; Aquaculture</b> |           |         |   |                  |           |                   |                   |
| PP2210  | P090      | 18-22ml | Actinomyces Agar  | 10 x 90mm Plates | Cellulose | 49                | 2-8°C             |
| PP2211  | P090      | 24-26ml | Actinomyces Agar (25ml)   | 10 x 90mm Plates | Cellulose | 49                | 2-8°C             |
| PP2151  | P090      | 18-22ml | Aeromonas Agar  | 10 x 90mm Plates | Cellulose | 77                | 2-8°C             |
| PP0051  | P090      | 18-22ml | Bacillus cereus Selective Agar (PEMBA)  | 10 x 90mm Plates | Cellulose | 70                | 2-8°C             |
| PP0050  | P090      | 18-22ml | Bacillus cereus (Mannitol Egg Yolk Polymixin – MYP) Agar                                    | 10 x 90mm Plates | Cellulose | 56                | 2-8°C             |
| PP0030  | P090      | 18-22ml | Bile Aesculin Agar  | 10 x 90mm Plates | Cellulose | 56                | 2-8°C             |
| PP1200  | P090      | 18-22ml | Blood Agar Base   | 10 x 90mm Plates | Cellulose | 56                | 2-8°C             |
| PP1220  | P090      | 18-22ml | Blood Agar No.2 with 7% Defibrinated Horse Blood  | 10 x 90mm Plates | Cellulose | 42                | 2-8°C             |
| PP1590  | P090      | 18-22ml | Brain Heart Infusion Agar   | 10 x 90mm Plates | Cellulose | 70                | 2-8°C             |
| PP1980  | P090      | 18-22ml | Brain Heart Infusion Agar with 10% Defibrinated Horse Blood                                 | 10 x 90mm Plates | Cellulose | 42                | 2-8°C             |
| PP1597  | P090      | 24-26ml | Brain Heart Infusion Agar with Vancomycin (6mg/L) (25ml)                                    | 10 x 90mm Plates | Cellulose | 42                | 2-8°C             |
| PP4070  | P090      | 18-22ml | Brazier's CCEY Agar with 1% Horse Blood Lysed, Cycloserine, Cefoxitin, 4% Egg Yolk Emulsion | 10 x 90mm Plates | Cellulose | 42                | 2-8°C             |
| PP1664  | P090      | 24-26ml | Brucella Agar with 5% Defibrinated Sheep Blood, Vitamin K & Haemin (25ml)                   | 10 x 90mm Plates | Cellulose | 42                | 2-8°C             |
| PP0160  | P090      | 18-22ml | Burkholderia cepacia Agar   | 10 x 90mm Plates | Cellulose | 56                | 2-8°C             |
| PP0080  | P090      | 18-22ml | CLED Agar   | 10 x 90mm Plates | Cellulose | 84                | 2-8°C             |
| PP3145  | P090      | 18-22ml | CLED Agar with Aztreonam (8mg/L)  | 10 x 90mm Plates | Cellulose | 56                | 2-8°C             |
| PP0540  | P090      | 18-22ml | CLED Agar (Bevis)   | 10 x 90mm Plates | Cellulose | 84                | 2-8°C             |
| PP0130  | P090      | 18-22ml | Campylobacter Blood Free CCDA   | 10 x 90mm Plates | Barrier   | 56                | 2-8°C             |
| PP0110  | P090      | 18-22ml | Campylobacter Selective Agar (Skirrow)  | 10 x 90mm Plates | Cellulose | 42                | 2-8°C             |



| Product Code | Container  | Volume           | Description  | Pack Size                                    | Packaging | Shelf Life (days) | Storage Temp (°C) |
|--------------|------------|------------------|--|--|-----------|-------------------|-------------------|
| PP0170       | P090       | 18-22ml          | Charcoal Agar with 10% Defibrinated Horse Blood  | 10 x 90mm Plates                             | Cellulose | 42                | 2-8°C             |
| PP0520       | P090       | 18-22ml          | Charcoal Agar with 10% Defibrinated Horse Blood & Cephalaxin                             | 10 x 90mm Plates                             | Cellulose | 42                | 2-8°C             |
| PP0100       | P090       | 18-22ml          | Chocolate Agar with 5% Defibrinated Horse Blood  | 10 x 90mm Plates                             | Cellulose | 84                | 2-8°C             |
| PP0530       | P090; PLSQ | 18-22ml; 55-57ml | Chocolate Agar with 7% Defibrinated Horse Blood & Bacitracin                             | 10 x 90mm Plates;<br>5 x 120mm Square Plates | Cellulose | 84                | 2-8°C             |
| PP0534       | P090       | 18-22ml          | Chocolate Agar with 1% Polyvitex   | 10 x 90mm Plates                             | Cellulose | 84                | 2-8°C             |
| PP1770       | PB90       | 10-12ml per side | Columbia Agar & Defibrinated Horse Blood / Chocolate Bacitracin Agar                     | 10 x 90mm Bi-plates                          | Cellulose | 42                | 2-8°C             |
| PP0090       | P090       | 18-22ml          | Columbia Agar Base   | 10 x 90mm Plates                             | Cellulose | 56                | 2-8°C             |
| PP0092       | P090       | 24-26ml          | Columbia Agar Base (25ml)  | 10 x 90mm Plates                             | Cellulose | 56                | 2-8°C             |
| PP2182       | P090       | 24-26ml          | Columbia Agar Base with 2% Sodium Chloride (25ml)  | 10 x 90mm Plates                             | Cellulose | 56                | 2-8°C             |
| PP0120       | P090       | 18-22ml          | Columbia Agar Base with 5% Defibrinated Horse Blood                                      | 10 x 90mm Plates                             | Cellulose | 42                | 2-8°C             |
| PP1084       | P090       | 18-22ml          | Columbia Agar Base with 5% Defibrinated Horse Blood & CAP Selective Supplement           | 10 x 90mm Plates                             | Cellulose | 42                | 2-8°C             |
| PP0040       | P090       | 18-22ml          | Columbia Agar Base with 5% Defibrinated Horse Blood & NAD (40mg/L)                       | 10 x 90mm Plates                             | Cellulose | 42                | 2-8°C             |
| PP1090       | P090       | 18-22ml          | Columbia Agar Base with 5% Defibrinated Horse Blood & Streptococcal Selective Supplement | 10 x 90mm Plates                             | Cellulose | 42                | 2-8°C             |
| PP1510       | P090       | 18-22ml          | Columbia Agar Base with 7% Defibrinated Horse Blood & Gardnerella Selective Supplement   | 10 x 90mm Plates                             | Cellulose | 42                | 2-8°C             |
| PP1730       | P090       | 18-22ml          | Columbia Agar Base with 7% Defibrinated Horse Blood & Neomycin (50mg/L)                  | 10 x 90mm Plates                             | Cellulose | 42                | 2-8°C             |
| PP1080       | P090       | 18-22ml          | Columbia Agar Base with 7% Defibrinated Horse Blood & CNA Selective Supplement           | 10 x 90mm Plates                             | Cellulose | 42                | 2-8°C             |
| PP1610       | P090       | 18-22ml          | Columbia Agar Base with 7% Defibrinated Sheep Blood & CNA Selective Supplement           | 10 x 90mm Plates                             | Cellulose | 42                | 2-8°C             |
| PP1747       | P090       | 18-22ml          | Columbia Agar Base with 5% Defibrinated Sheep Blood                                      | 10 x 90mm Plates                             | Cellulose | 42                | 2-8°C             |

| Product Code | Container  | Volume           | Description   | Pack Size                                     | Packaging | Shelf Life (days) | Storage Temp (°C) |
|--------------|------------|------------------|---|---|-----------|-------------------|-------------------|
| PP1750       | P090       | 18-22ml          | Columbia Agar Base with 7% Defibrinated Sheep Blood   | 10 x 90mm Plates                              | Cellulose | 42                | 2-8°C             |
| PP1785       | PB90       | 10-12ml per side | Columbia CNA & Defibrinated Sheep Blood/ Sabouraud Dextrose Agar with Chloramphenicol                             | 10 x 90mm Bi-plates                           | Cellulose | 42                | 2-8°C             |
| PP4120       | P090       | 18-22ml          | Cornmeal Agar   | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP1240       | P090; PCON | 18-22ml; 15-17ml | Czapek Dox Agar (Modified)  | 10 x 90mm Plates;<br>10 x 55mm Contact Plates | Cellulose | 56                | 2-8°C             |
| PP0390       | P090       | 18-22ml          | Dermasel Agar Base with Chloramphenicol & Cyclohexamide (Actidione)   | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP2220       | P090       | 18-22ml          | Dermatophyte Test Medium (DTM) with Chloramphenicol (100mg/L)   | 10 x 90mm Plates                              | Cellulose | 70                | 2-8°C             |
| PP0450       | P090       | 18-22ml          | Dermatophyte Test Medium (DTM) with Chloramphenicol (500mg/L)   | 10 x 90mm Plates                              | Cellulose | 70                | 2-8°C             |
| PP3150       | P090       | 18-22ml          | Dermatophyte Test Medium (DTM) with Chloramphenicol & Cyclohexamide (Actidione)                                   | 10 x 90mm Plates                              | Cellulose | 70                | 2-8°C             |
| PP0370       | P090       | 18-22ml          | Desoxycholate Citrate Agar (Hynes)  | 10 x 90mm Plates                              | Cellulose | 70                | 2-8°C             |
| PP0400       | P090; PLSQ | 18-22ml; 55-57ml | DNase Agar  | 10 x 90mm Plates;<br>5 x 120mm Square Plates  | Cellulose | 56                | 2-8°C             |
| PP0560       | P090       | 18-22ml          | DNase Agar with Methyl Green  | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP0710       | P090       | 18-22ml          | Edwards Medium (Modified) with 7% Defibrinated Sheep Blood  | 10 x 90mm Plates                              | Cellulose | 42                | 2-8°C             |
| PP1630       | P090       | 18-22ml          | Enterococcosel Agar   | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP1563       | PB90       | 10-12ml per side | Fastidious Anaerobe Agar (FAA) /Fastidious Anaerobe Agar (FAA) with Neomycin (75mg/L)                             | 10 x 90mm Bi-Plates                           | Cellulose | 84                | 2-8°C             |
| PP0141       | PB90       | 10-12ml per side | Fastidious Anaerobe Agar (FAA) with Neomycin (75mg/L)/Fastidious Anaerobe Agar (FAA) with Nalidixic Acid (20mg/L) | 10 x 90mm Bi-plates                           | Cellulose | 77                | 2-8°C             |
| PP1560       | P090       | 18-22ml          | Fastidious Anaerobe Agar (FAA) with 7% Defibrinated Horse Blood   | 10 x 90mm Plates                              | Cellulose | 84                | 2-8°C             |
| PP3180       | P090       | 18-22ml          | Fastidious Anaerobe Agar (FAA) with 7% Defibrinated Horse Blood & NAT Selective Supplement                        | 10 x 90mm Plates                              | Cellulose | 77                | 2-8°C             |

| Product Code | Container | Volume            | Description   | Pack Size            | Packaging | Shelf Life (days) | Storage Temp (°C) |
|--------------|-----------|-------------------|---|----------------------|-----------|-------------------|-------------------|
| PP0140       | P090      | 18-22ml           | Fastidious Anaerobe Agar (FAA) with 7% Defibrinated Horse Blood & Neomycin (75mg/L)                           | 10 x 90mm Plates     | Cellulose | 84                | 2-8°C             |
| PP0144       | P090      | 18-22ml           | Fastidious Anaerobe Agar (FAA) with 7% Defibrinated Horse Blood, Neomycin (75mg/L) & Aztreonam (2mg/L)        | 10 x 90mm Plates     | Cellulose | 84                | 2-8°C             |
| PP1500       | P090      | 18-22ml           | Fastidious Anaerobe Agar (FAA) with 7% Defibrinated Horse Blood, Vancomycin (6mg/L) & Nalidixic Acid (10mg/L) | 10 x 90mm Plates     | Cellulose | 42                | 2-8°C             |
| PP1561       | P090      | 18-22ml           | Fastidious Anaerobe Agar (FAA) with 5% Defibrinated Sheep Blood   | 10 x 90mm Plates     | Cellulose | 42                | 2-8°C             |
| PP1782       | PB90      | 10-12ml per side  | GC Selective Agar with LCAT / GC Selective Agar with VCAT   | 10 x 90mm Bi-plates  | Cellulose | 56                | 2-8°C             |
| PP1570       | P090      | 18-22ml           | GC Chocolate Agar with VCAT   | 10 x 90mm Plates     | Cellulose | 70                | 2-8°C             |
| PP0782       | P090      | 18-22ml           | GC Chocolate Agar with LCAT   | 10 x 90mm Plates     | Cellulose | 70                | 2-8°C             |
| PP0780       | P090      | 18-22ml           | GC Selective Agar with LCAT   | 10 x 90mm Plates     | Cellulose | 56                | 2-8°C             |
| PP0790       | P090      | 18-22ml           | GC Selective Agar with VCAT   | 10 x 90mm Plates     | Cellulose | 56                | 2-8°C             |
| PP0800       | P090      | 18-22ml           | GC Selective Agar with VCNT   | 10 x 90mm Plates     | Cellulose | 56                | 2-8°C             |
| PP2071       | P090      | 18-22ml           | GC Agar Non-Selective   | 10 x 90mm Plates     | Cellulose | 56                | 2-8°C             |
| PP1822       | P090      | 24-26ml           | Haemophilus Test Medium (25ml)  | 10 x 90mm Plates     | Cellulose | 56                | 2-8°C             |
| PP1678       | PT90      | 7-9ml per section | Hektoen Enteric / Salmonella Shigella (SS) / XLD Agars  | 10 x 90mm Tri-plates | Cellulose | 70                | 2-8°C             |
| PP1520       | P090      | 18-22ml           | Hektoen Enteric Agar  | 10 x 90mm Plates     | Cellulose | 70                | 2-8°C             |
| PP2110       | P090      | 18-22ml           | Helicobacter Pylori Medium  | 10 x 90mm Plates     | Cellulose | 42                | 2-8°C             |
| PP0610       | P090      | 18-22ml           | Hoyles Medium   | 10 x 90mm Plates     | Cellulose | 42                | 2-8°C             |
| PP2138       | P090      | 24-26ml           | Iso-Sensitest Agar (25ml)   | 10 x 90mm Plates     | Cellulose | 70                | 2-8°C             |

| Product Code | Container | Volume  | Description  | Pack Size        | Packaging | Shelf Life (days) | Storage Temp (°C) |
|--------------|-----------|---------|--|------------------|-----------|-------------------|-------------------|
| PP2146       | P090      | 24-26ml | Iso-Sensitest Agar with 5% Defibrinated Horse Blood (25ml)                                   | 10 x 90mm Plates | Cellulose | 42                | 2-8°C             |
| PP2148       | P090      | 24-26ml | Iso-Sensitest with 5% Defibrinated Horse Blood & NAD (20mg/L) (25ml)                         | 10 x 90mm Plates | Cellulose | 42                | 2-8°C             |
| PP3160       | P090      | 18-22ml | Kanamycin Aesculin Agar with Vancomycin (6mg/L)  | 10 x 90mm Plates | Cellulose | 56                | 2-8°C             |
| PP0850       | P090      | 18-22ml | Kanamycin Aesculin Azide Agar  | 10 x 90mm Plates | Cellulose | 84                | 2-8°C             |
| PP3349       | P090      | 18-22ml | LB Agar (Lennox)   | 10 x 90mm Plates | Cellulose | 70                | 2-8°C             |
| PP3354       | P090      | 18-22ml | LB Agar (Lennox) with Ampicillin (100mg/L)   | 10 x 90mm Plates | Cellulose | 56                | 2-8°C             |
| PP3348       | P090      | 18-22ml | LB Agar (Lennox) with Carbenicillin (100mg/L) & 1% Glucose                                   | 10 x 90mm Plates | Cellulose | 56                | 2-8°C             |
| PP3395       | P090      | 18-22ml | LB Agar (Lennox) with Kanamycin (50mg/L)   | 10 x 90mm Plates | Cellulose | 56                | 2-8°C             |
| PP0201       | P090      | 18-22ml | Legionella BCYE Cysteine Free Agar   | 10 x 90mm Plates | Cellulose | 56                | 2-8°C             |
| PP0070       | P090      | 18-22ml | Legionella BMPA Agar   | 10 x 90mm Plates | Cellulose | 56                | 2-8°C             |
| PP0200       | P090      | 18-22ml | Legionella CYE Agar  | 10 x 90mm Plates | Cellulose | 56                | 2-8°C             |
| PP0860       | P090      | 18-22ml | Legionella MWY Agar  | 10 x 90mm Plates | Cellulose | 56                | 2-8°C             |
| PP0630       | P090      | 18-22ml | Listeria (Oxford) Agar   | 10 x 90mm Plates | Cellulose | 84                | 2-8°C             |
| PP0640       | P090      | 18-22ml | MacConkey Agar   | 10 x 90mm Plates | Cellulose | 56                | 2-8°C             |
| PP1725       | P090      | 18-22ml | MacConkey Agar No.2  | 10 x 90mm Plates | Cellulose | 70                | 2-8°C             |
| PP1720       | P090      | 18-22ml | MacConkey Agar No.3  | 10 x 90mm Plates | Cellulose | 70                | 2-8°C             |
| PP0470       | P090      | 18-22ml | MacConkey Agar without Salt  | 10 x 90mm Plates | Cellulose | 77                | 2-8°C             |
| PP0897       | P090      | 24-26ml | Malt Extract Agar with Chloramphenicol (50mg/L) & Cyclohexamide (Actidione) (400mg/L) (25ml) | 10 x 90mm Plates | Cellulose | 70                | 2-8°C             |

| Product Code | Container  | Volume              | Description   | Pack Size                                     | Packaging | Shelf Life (days) | Storage Temp (°C) |
|--------------|------------|---------------------|---|---|-----------|-------------------|-------------------|
| PP0891       | P090       | 18-22ml             | Malt Extract Agar with Chloramphenicol (80mg/L)   | 10 x 90mm Plates                              | Cellulose | 70                | 2-8°C             |
| PP0660       | P090       | 18-22ml             | Mannitol Salt Agar (USP)  | 10 x 90mm Plates                              | Cellulose | 70                | 2-8°C             |
| PP4085       | P090       | 24-26ml             | Middlebrook 7H10 Selective Agar (25ml)  | 10 x 90mm Plates                              | Cellulose | 70                | 2-8°C             |
| PP4080       | P090       | 24-26ml             | Middlebrook 7H11 Selective Agar (25ml)  | 10 x 90mm Plates                              | Cellulose | 70                | 2-8°C             |
| PP0960       | P140; PLSQ | 60-65ml;<br>55-57ml | Mueller-Hinton Agar   | 5 x 140mm Plates;<br>5 x 120mm Square Plates  | Cellulose | 70                | 2-8°C             |
| PP0963       | P090       | 24-26ml             | Mueller-Hinton Agar (25ml)  | 10 x 90mm Plates                              | Cellulose | 70                | 2-8°C             |
| PP0972       | P090; PLSQ | 24-26ml; 55-57ml    | Mueller-Hinton Agar with 5% Defibrinated Horse Blood & NAD (20mg/L)                         | 10 x 90mm Plates;<br>5 x 120mm Square Plates  | Cellulose | 42                | 2-8°C             |
| PP0973       | P090       | 24-26ml             | Mueller-Hinton Agar with 7% Defibrinated Horse Blood (25ml)                                 | 10 x 90mm Plates                              | Cellulose | 42                | 2-8°C             |
| PP3131       | P090       | 24-26ml             | Mueller-Hinton Agar with 2% Sodium Chloride (25ml)  | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP0965       | P090       | 18-22ml             | Mueller-Hinton GISA   | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP1645       | P090       | 24-26ml             | Mueller-Hinton Agar with 5% Defibrinated Sheep Blood (25ml)                                 | 10 x 90mm Plates                              | Cellulose | 42                | 2-8°C             |
| PP3175       | P090       | 24-26ml             | Mueller-Hinton Agar with 2% Glucose & Methylene Blue (25ml)                                 | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP0235       | P090       | 24-26ml             | Mycological Agar with Chloramphenicol (50mg/L) & Cyclohexamide (Actidione) (400mg/L) (25ml) | 10 x 90mm Plates                              | Cellulose | 70                | 2-8°C             |
| PP0670       | P090       | 18-22ml             | Nagler Medium   | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP2161       | P055       | 17-19ml             | Non-Nutrient Agar   | 10 x 55mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP0690       | P090; PCON | 18-22ml; 15-17ml    | Nutrient Agar   | 10 x 90mm Plates;<br>10 x 55mm Contact Plates | Cellulose | 84                | 2-8°C             |
| PP1740       | P090       | 18-22ml             | Nutrient Agar with 5% Defibrinated Horse Blood  | 10 x 90mm Plates                              | Cellulose | 42                | 2-8°C             |

| Product Code | Container  | Volume           | Description   | Pack Size                                     | Packaging | Shelf Life (days) | Storage Temp (°C) |
|--------------|------------|------------------|---|---|-----------|-------------------|-------------------|
| PP0005       | P090       | 18-22ml          | PALCAM Agar   | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP0004       | P090       | 18-22ml          | PHAGE Agar  | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP2160       | P090       | 18-22ml          | Page's Amoeba Saline & Agar No.1  | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP1020       | P090       | 18-22ml          | Phenylethyl Alcohol (PEA) Agar  | 10 x 90mm Plates                              | Cellulose | 42                | 2-8°C             |
| PP1041       | P090       | 18-22ml          | Potato Dextrose Agar with Chloramphenicol (100mg/L)   | 10 x 90mm Plates                              | Cellulose | 84                | 2-8°C             |
| PP0730       | P090       | 18-22ml          | Proteose Peptone Rice Agar  | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP1060       | P090       | 18-22ml          | Pseudomonas Agar Base with 1% Glycerol & CN Selective Supplement                            | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP1051       | P090       | 18-22ml          | Pseudomonas Agar Base with 1% Glycerol, CFC Selective Supplement                            | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP1870       | P090       | 18-22ml          | RPMI Medium for E-Test  | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP1100       | P090       | 18-22ml          | Salmonella Shigella (SS) Agar (Modified)  | 10 x 90mm Plates                              | Cellulose | 70                | 2-8°C             |
| PP0240       | P090; PCON | 18-22ml; 15-17ml | Sabouraud Dextrose Agar   | 10 x 90mm Plates;<br>10 x 55mm Contact Plates | Cellulose | 70                | 2-8°C             |
| PP0243       | P090       | 24-26ml          | Sabouraud Dextrose Agar (25ml)  | 10 x 90mm Plates                              | Cellulose | 70                | 2-8°C             |
| PP0220       | P090; PCON | 18-22ml; 15-17ml | Sabouraud Dextrose Agar with Chloramphenicol  | 10 x 90mm Plates;<br>10 x 55mm Contact Plates | Cellulose | 84                | 2-8°C             |
| PP0231       | P090       | 18-22ml          | Sabouraud Dextrose Agar with Chloramphenicol (50mg/L) & Cyclohexamide (Actidione) (300mg/L) | 10 x 90mm Plates                              | Cellulose | 70                | 2-8°C             |
| PP2360       | P090       | 24-26ml          | Sabouraud Dextrose Agar with Chloramphenicol (50mg/L) (25ml)                                | 10 x 90mm Plates                              | Cellulose | 84                | 2-8°C             |
| PP0247       | P090       | 24-26ml          | Sabouraud Dextrose Agar with Chloramphenicol (50mg/L) & Cyclohexamide (Actidione) (300mg/L) | 10 x 90mm Plates                              | Cellulose | 70                | 2-8°C             |
| PP1543       | P090       | 18-22ml          | Sabouraud Dextrose Agar with Colistin (20mg/L) & Gentamicin (20mg/L)                        | 10 x 90mm Plates                              | Cellulose | 70                | 2-8°C             |
| PP1660       | P090       | 18-22ml          | Schaedler Anaerobe Agar with 7% Defibrinated Sheep Blood & Vitamin K                        | 10 x 90mm Plates                              | Cellulose | 42                | 2-8°C             |

| Product Code                           | Container  | Volume           | Description   | Pack Size                                     | Packaging | Shelf Life (days) | Storage Temp (°C) |
|--|------------|------------------|---|---|-----------|-------------------|-------------------|
| PP1830                                 | P090       | 18-22ml          | Shadomy Agar  | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP0250                                 | P090       | 18-22ml          | Sorbitol MacConkey Agar                                     | 10 x 90mm Plates                              | Cellulose | 70                | 2-8°C             |
| PP1942                                 | P090       | 18-22ml          | Sorbitol MacConkey Agar with Cefixime Rhamnose              | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP1940                                 | P090       | 18-22ml          | Sorbitol MacConkey Agar with Cefixime & Tellurite (CT-SMAC) | 10 x 90mm Plates                              | Cellulose | 63                | 2-8°C             |
| PP0270                                 | P090       | 18-22ml          | TCBS Cholera Medium   | 10 x 90mm Plates                              | Cellulose | 42                | 2-8°C             |
| PP0003                                 | P090       | 18-22ml          | Thermostable Nuclease Test (TSN) Agar                       | 10 x 90mm Plates                              | Cellulose | 42                | 2-8°C             |
| PP0280                                 | P090; PCON | 18-22ml; 15-17ml | Tryptone Soya Agar  | 10 x 90mm Plates;<br>10 x 55mm Contact Plates | Cellulose | 84                | 2-8°C             |
| PP1651                                 | P090       | 18-22ml          | Tryptone Soya Agar with 5% Defibrinated Sheep Blood         | 10 x 90mm Plates                              | Cellulose | 42                | 2-8°C             |
| PP2000                                 | P090       | 18-22ml          | Urea Agar   | 10 x 90mm Plates                              | Cellulose | 42                | 2-8°C             |
| PP2001                                 | P090       | 18-22ml          | Urea Agar with Kaolin                                       | 10 x 90mm Plates                              | Cellulose | 42                | 2-8°C             |
| PP0320                                 | P090       | 18-22ml          | XLD Agar  | 10 x 90mm Plates                              | Cellulose | 70                | 2-8°C             |
| PP1791                                 | P090       | 24-26ml          | Yeast Sensitivity Agar (25ml)                               | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP0350                                 | P090       | 18-22ml          | Yersinia Selective Agar with CIN                            | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| <b>Food, Water &amp; Environmental</b> |            |                  |   |   |           |                   |                   |
| PP2151                                 | P090       | 18-22ml          | Aeromonas Agar  | 10 x 90mm Plates                              | Cellulose | 77                | 2-8°C             |
| PP0051                                 | P090       | 18-22ml          | Bacillus cereus Selective Agar (PEMBA)                      | 10 x 90mm Plates                              | Cellulose | 70                | 2-8°C             |
| PP0050                                 | P090       | 18-22ml          | Bacillus cereus (Mannitol Egg Yolk Polymyxin – MYP) Agar    | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP1420                                 | P090; PCON | 18-22ml; 15-17ml | Baird Parker Agar with 5% Egg Yolk Tellurite                | 10 x 90mm Plates;<br>10 x 55mm Contact Plates | Cellulose | 84                | 2-8°C             |

| Product Code | Container  | Volume           | Description   | Pack Size                                    | Packaging | Shelf Life (days) | Storage Temp (°C) |
|--------------|------------|------------------|---|--|-----------|-------------------|-------------------|
| PP0030       | P090       | 18-22ml          | Bile Aesculin Agar                                  | 10 x 90mm Plates                             | Cellulose | 56                | 2-8°C             |
| PP0209       | P090       | 18-22ml          | Blood Agar Base with 5% Defibrinated Sheep Blood    | 10 x 90mm Plates                             | Cellulose | 42                | 2-8°C             |
| PP1220       | P090       | 18-22ml          | Blood Agar No.2 with 7% Defibrinated Horse Blood    | 10 x 90mm Plates                             | Cellulose | 42                | 2-8°C             |
| PP0060       | P090       | 18-22ml          | Brilliant Green Modified Agar                       | 10 x 90mm Plates                             | Cellulose | 84                | 2-8°C             |
| PP0130       | P090       | 18-22ml          | Campylobacter Blood Free CCDA                       | 10 x 90mm Plates                             | Barrier   | 56                | 2-8°C             |
| PP1700       | P090       | 18-22ml          | Campylobacter Selective Agar (Preston)              | 10 x 90mm Plates                             | Cellulose | 42                | 2-8°C             |
| PP0110       | P090       | 18-22ml          | Campylobacter Selective Agar (Skirrow)              | 10 x 90mm Plates                             | Cellulose | 42                | 2-8°C             |
| PP2230       | P090       | 18-22ml          | Cetrimide Agar with 1% Glycerol (USP)               | 10 x 90mm Plates                             | Cellulose | 56                | 2-8°C             |
| PP0120       | P090       | 18-22ml          | Columbia Agar Base with 5% Defibrinated Horse Blood | 10 x 90mm Plates                             | Cellulose | 42                | 2-8°C             |
| PP1747       | P090       | 18-22ml          | Columbia Agar Base with 5% Defibrinated Sheep Blood | 10 x 90mm Plates                             | Cellulose | 42                | 2-8°C             |
| PP4253       | P090       | 18-22ml          | Dichloran Rose-Bengal Chloramphenicol (DRBC) Agar   | 10 x 90mm Plates                             | Cellulose | 70                | 2-8°C             |
| PP0400       | P090; PLSQ | 18-22ml; 55-57ml | DNase Agar  | 10 x 90mm Plates;<br>5 x 120mm Square Plates | Cellulose | 56                | 2-8°C             |
| PP0560       | P090       | 18-22ml          | DNase Agar with Methyl Green                        | 10 x 90mm Plates                             | Cellulose | 56                | 2-8°C             |
| PP1630       | P090       | 18-22ml          | Enterococcosel Agar                                 | 10 x 90mm Plates                             | Cellulose | 56                | 2-8°C             |
| PP1520       | P090       | 18-22ml          | Hektoen Enteric Agar                                | 10 x 90mm Plates                             | Cellulose | 70                | 2-8°C             |
| PP0850       | P090       | 18-22ml          | Kanamycin Aesculin Azide Agar                       | 10 x 90mm Plates                             | Cellulose | 84                | 2-8°C             |
| PP0201       | P090       | 18-22ml          | Legionella BCYE Cysteine Free Agar                  | 10 x 90mm Plates                             | Cellulose | 56                | 2-8°C             |
| PP0870       | P090       | 18-22ml          | Legionella BCYE with GVPC                           | 10 x 90mm Plates                             | Cellulose | 56                | 2-8°C             |



| Product Code | Container  | Volume           | Description                                       | Pack Size                                     | Packaging | Shelf Life (days) | Storage Temp (°C) |
|--------------|------------|------------------|---|---|-----------|-------------------|-------------------|
| PP0200       | P090       | 18-22ml          | Legionella CYE Agar                               | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP0860       | P090       | 18-22ml          | Legionella MWY Agar                               | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP0630       | P090       | 18-22ml          | Listeria (Oxford) Agar                            | 10 x 90mm Plates                              | Cellulose | 84                | 2-8°C             |
| PP4150       | P090       | 18-22ml          | Lysine Medium                                     | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP0954       | P090       | 18-22ml          | MRS (de Man, Rogosa, Sharpe) Agar (Low pH)        | 10 x 90mm Plates                              | Cellulose | 70                | 2-8°C             |
| PP0640       | P090       | 18-22ml          | MacConkey Agar                                    | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP1725       | P090       | 18-22ml          | MacConkey Agar No.2                               | 10 x 90mm Plates                              | Cellulose | 70                | 2-8°C             |
| PP1720       | P090       | 18-22ml          | MacConkey Agar No.3                               | 10 x 90mm Plates                              | Cellulose | 70                | 2-8°C             |
| PP0478       | P090       | 18-22ml          | MacConkey Agar No.3 with Cefotaxime (1mg/L)       | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP0890       | P090       | 18-22ml          | Malt Extract Agar                                 | 10 x 90mm Plates                              | Cellulose | 84                | 2-8°C             |
| PP0901       | P090       | 18-22ml          | Malt Extract Agar (pH 3.5)                        | 10 x 90mm Plates                              | Cellulose | 84                | 2-8°C             |
| PP0891       | P090       | 18-22ml          | Malt Extract Agar with Chloramphenicol (80mg/L)   | 10 x 90mm Plates                              | Cellulose | 70                | 2-8°C             |
| PP0935       | P090       | 18-22ml          | Milk Agar with Cetrimide                          | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP4130       | P090       | 18-22ml          | Mineral Modified Glutamate Medium                 | 10 x 90mm Plates                              | Cellulose | 84                | 2-8°C             |
| PP0670       | P090       | 18-22ml          | Nagler Medium                                     | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP0690       | P090; PCON | 18-22ml; 15-17ml | Nutrient Agar                                     | 10 x 90mm Plates;<br>10 x 55mm Contact Plates | Cellulose | 84                | 2-8°C             |
| PP1000       | P090; PCON | 18-22ml; 15-17ml | Oxytetracycline Glucose Yeast Extract Agar (OGYE) | 10 x 90mm Plates;<br>10 x 55mm Contact Plates | Cellulose | 84                | 2-8°C             |
| PP0005       | P090       | 18-22ml          | PALCAM Agar                                       | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |

| Product Code | Container  | Volume           | Description  | Pack Size                                     | Packaging | Shelf Life (days) | Storage Temp (°C) |
|--------------|------------|------------------|--|---|-----------|-------------------|-------------------|
| PP1030       | P090; PCON | 18-22ml; 15-17ml | Plate Count Agar   | 10 x 90mm Plates;<br>10 x 55mm Contact Plates | Cellulose | 84                | 2-8°C             |
| PP1041       | P090       | 18-22ml          | Potato Dextrose Agar with Chloramphenicol (100mg/L)              | 10 x 90mm Plates                              | Cellulose | 84                | 2-8°C             |
| PP1060       | P090       | 18-22ml          | Pseudomonas Agar Base with 1% Glycerol & CN Selective Supplement | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP1051       | P090       | 18-22ml          | Pseudomonas Agar Base with 1% Glycerol, CFC Selective Supplement | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP1050       | P090       | 18-22ml          | Pseudomonas Selective Agar with CFC Supplement                   | 10 x 90mm Plates                              | Cellulose | 42                | 2-8°C             |
| PP1281       | P090       | 18-22ml          | R2A Agar   | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP1285       | P090       | 18-22ml          | Raka Ray Agar with Cyclohexamide & Tween 80                      | 10 x 90mm Plates                              | Cellulose | 70                | 2-8°C             |
| PP0240       | P090; PCON | 18-22ml; 15-17ml | Sabouraud Dextrose Agar  | 10 x 90mm Plates;<br>10 x 55mm Contact Plates | Cellulose | 70                | 2-8°C             |
| PP0243       | P090       | 24-26ml          | Sabouraud Dextrose Agar (25ml)                                   | 10 x 90mm Plates                              | Cellulose | 70                | 2-8°C             |
| PP0260       | P090       | 18-22ml          | Slanetz & Bartley Medium   | 10 x 90mm Plates                              | Cellulose | 70                | 2-8°C             |
| PP0250       | P090       | 18-22ml          | Sorbitol MacConkey Agar  | 10 x 90mm Plates                              | Cellulose | 70                | 2-8°C             |
| PP1940       | P090       | 18-22ml          | Sorbitol MacConkey Agar with Cefixime & Tellurite (CT-SMAC)      | 10 x 90mm Plates                              | Cellulose | 63                | 2-8°C             |
| PP0270       | P090       | 18-22ml          | TCBS Cholera Medium  | 10 x 90mm Plates                              | Cellulose | 42                | 2-8°C             |
| PP0285       | P090       | 18-22ml          | Triphenyltetrazolium Chloride Soya Tryptone Agar                 | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP2190       | P090       | 18-22ml          | Tryptone Glucose Extract Agar                                    | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP0280       | P090; PCON | 18-22ml; 15-17ml | Tryptone Soya Agar   | 10 x 90mm Plates;<br>10 x 55mm Contact Plates | Cellulose | 84                | 2-8°C             |
| PP0300       | P090       | 18-22ml          | Tryptone Soya Yeast Extract (TSYE) Agar                          | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP1950       | P090       | 18-22ml          | Tryptose Sulphite Cycloserine (TSC) Agar                         | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |

| Product Code | Container  | Volume           | Description   | Pack Size                                     | Packaging | Shelf Life (days) | Storage Temp (°C) |
|--------------|------------|------------------|---|---|-----------|-------------------|-------------------|
| PP2000       | P090       | 18-22ml          | Urea Agar   | 10 x 90mm Plates                              | Cellulose | 42                | 2-8°C             |
| PP2001       | P090       | 18-22ml          | Urea Agar with Kaolin                                 | 10 x 90mm Plates                              | Cellulose | 42                | 2-8°C             |
| PP1150       | P090; PCON | 18-22ml; 15-17ml | Violet Red Bile Agar (VRBA)                           | 10 x 90mm Plates;<br>10 x 55mm Contact Plates | Cellulose | 70                | 2-8°C             |
| PP1160       | P090; PCON | 18-22ml; 15-17ml | Violet Red Bile Glucose Agar (VRBGA)                  | 10 x 90mm Plates;<br>10 x 55mm Contact Plates | Cellulose | 70                | 2-8°C             |
| PP1140       | P090       | 18-22ml          | Vogel Johnson Agar                                    | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP1170       | P090       | 18-22ml          | WL (Wallerstein Laboratory) Nutrient Agar             | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP0320       | P090       | 18-22ml          | XLD Agar  | 10 x 90mm Plates                              | Cellulose | 70                | 2-8°C             |
| PP1191       | P090; PCON | 18-22ml; 15-17ml | Yeast Extract Dextrose Chloramphenicol (100mg/L) Agar | 10 x 90mm Plates;<br>10 x 55mm Contact Plates | Cellulose | 100               | 2-8°C             |
| PP0995       | P090       | 18-22ml          | Yeast Malt Agar with Oxytetracycline                  | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |
| PP0350       | P090       | 18-22ml          | Yersinia Selective Agar with CIN                      | 10 x 90mm Plates                              | Cellulose | 56                | 2-8°C             |

| Product Code                                     | Container | Volume  | Description  | Pack Size                | Packaging | Shelf Life (days) | Storage Temp (°C) | USP | Irradiated | Ambient Storage |
|--|-----------|---------|--|--------------------------|-----------|-------------------|-------------------|-----|------------|-----------------|
| <b>Pharmaceutical, Industrial &amp; Cosmetic</b> |           |         |  |                          |           |                   |                   |     |            |                 |
| PP2230   | P090      | 18-22ml | Cetrimide Agar with 1% Glycerol (USP)                      | 10 x 90mm Plates         | Cellulose | 56                | 2-8°C             | ●   |            |                 |
| PP1380   | PCON      | 15-17ml | D/E Neutralising Agar                                      | 10 x 55mm Contact Plates | Cellulose | 42                | 2-8°C             |     |            |                 |
| PP3349   | P090      | 18-22ml | LB Agar (Lennox)   | 10 x 90mm Plates         | Cellulose | 70                | 2-8°C             |     |            |                 |
| PP3354   | P090      | 18-22ml | LB Agar (Lennox) with Ampicillin (100mg/L)                 | 10 x 90mm Plates         | Cellulose | 56                | 2-8°C             |     |            |                 |
| PP3348   | P090      | 18-22ml | LB Agar (Lennox) with Carbenicillin (100mg/L) & 1% Glucose | 10 x 90mm Plates         | Cellulose | 56                | 2-8°C             |     |            |                 |
| PP3395   | P090      | 18-22ml | LB Agar (Lennox) with Kanamycin (50mg/L)                   | 10 x 90mm Plates         | Cellulose | 56                | 2-8°C             |     |            |                 |
| PP1723   | P090      | 18-22ml | MacConkey Agar (USP)                                       | 10 x 90mm Plates         | Cellulose | 70                | 2-8°C             | ●   |            |                 |

| Product Code | Container  | Volume              | Description   | Pack Size                                      | Packaging | Shelf Life (days) | Storage Temp (°C) | USP | Irradiated | Ambient Storage |
|--------------|------------|---------------------|---|--|-----------|-------------------|-------------------|-----|------------|-----------------|
| PP0660       | P090       | 18-22ml             | Mannitol Salt Agar (USP)  | 10 x 90mm Plates                               | Cellulose | 70                | 2-8°C             | ●   |            |                 |
| PP4085       | P090       | 24-26ml             | Middlebrook 7H10 Selective Agar (25ml)  | 10 x 90mm Plates                               | Cellulose | 70                | 2-8°C             |     |            |                 |
| PP4080       | P090       | 24-26ml             | Middlebrook 7H11 Selective Agar (25ml)  | 10 x 90mm Plates                               | Cellulose | 70                | 2-8°C             |     |            |                 |
| PP1030       | P090; PCON | 18-22ml;<br>15-17ml | Plate Count Agar  | 10 x 90mm Plates;<br>10 x 55mm Contact Plates  | Cellulose | 84                | 2-8°C             |     |            |                 |
| PP1040       | P090       | 18-22ml             | Potato Dextrose Agar (USP)  | 10 x 90mm Plates                               | Cellulose | 84                | 2-8°C             | ●   |            |                 |
| PP1281       | P090       | 18-22ml             | R2A Agar  | 10 x 90mm Plates                               | Cellulose | 56                | 2-8°C             |     |            |                 |
| PP8030       | P090; PCON | 24-26ml;<br>15-17ml | Sabouraud Dextrose Agar (VHP) (Irradiated)  | 10 x 90 mm Plates;<br>10 x 55mm Contact Plates | Barrier   | 182               | 2-30°C            | ●   | ●          | ●               |
| PP8031       | P090; PCON | 24-26ml;<br>15-17ml | Sabouraud Dextrose Agar with Lecithin & Tween 80 (VHP) (Irradiated)                                 | 10 x 90mm Plates;<br>10 x 55mm Contact Plates  | Barrier   | 182               | 2-30°C            | ●   | ●          | ●               |
| PP8033       | P090; PCON | 24-26ml;<br>15-17ml | Sabouraud Dextrose Agar with Lecithin, Tween 80, Histidine & Sodium Thiosulphate (VHP) (Irradiated) | 10 x 90mm Plates;<br>10 x 55mm Contact Plates  | Barrier   | 182               | 2-30°C            | ●   | ●          | ●               |
| PP0245       | P140       | 60-65ml             | Sabouraud Dextrose Agar (Irradiated)  | 5 x 140mm Plates                               | Cellulose | 182               | 2-8°C             | ●   | ●          |                 |
| PP0280       | P090; PCON | 18-22ml; 15-17ml    | Tryptone Soya Agar  | 10 x 90mm plates; 10 x 55mm Contact Plates     | Cellulose | 84                | 2-8°C             | ●   |            | ●               |
| PP1350       | P090; PCON | 18-22ml; 15-17ml    | Tryptone Soya Agar with Lecithin & Tween 80   | 10 x 90mm plates; 10 x 55mm Contact Plates     | Cellulose | 84                | 2-8°C             |     |            |                 |
| PP0282       | P140       | 60-65ml             | Tryptone Soya Agar (Irradiated)   | 5 x 140mm Plates                               | Cellulose | 182               | 2-8°C             | ●   | ●          |                 |
| PP8020       | P090; PCON | 24-26ml; 15-17ml    | Tryptone Soya Agar (VHP) (Irradiated)   | 10 x 90mm Plates;<br>10 x 55mm Contact Plates  | Barrier   | 182               | 2-30°C            | ●   | ●          | ●               |
| PP8022       | P090; PCON | 24-26ml; 15-17ml    | Tryptone Soya Agar with Lecithin & Tween 80 (VHP) (Irradiated)                                      | 10 x 90mm Plates; 10 x 55mm Contact Plates     | Barrier   | 182               | 2-30°C            | ●   | ●          | ●               |
| PP8024       | P090; PCON | 24-26ml; 15-17ml    | Tryptone Soya Agar with Lecithin, Tween 80, Histidine & Sodium Thiosulphate (VHP) (Irradiated)      | 10 x 90mm Plates; 10 x 55mm Contact Plates     | Barrier   | 182               | 2-30°C            | ●   | ●          | ●               |
| PP4284       | P090       | 18-22ml             | Tween Modified Lethen Agar  | 10 x 90mm Plates                               | Cellulose | 56                | 2-8°C             |     |            |                 |
| PP1162       | P090       | 18-22ml             | Violet Red Bile Glucose Agar (VRBGA) (USP)  | 10 x 90mm Plates                               | Cellulose | 70                | 2-8°C             | ●   |            |                 |
| PP2035       | P090       | 18-22ml             | Vogel Bonner E  | 10 x 90mm Plates                               | Cellulose | 84                | 2-8°C             |     |            |                 |
| PP1140       | P090       | 18-22ml             | Vogel Johnson Agar  | 10 x 90mm Plates                               | Cellulose | 56                | 2-8°C             |     |            |                 |
| PP0321       | P090       | 18-22ml             | XLD Medium (USP)  | 10 x 90mm Plates                               | Cellulose | 70                | 2-8°C             | ●   |            |                 |

| Product Code       | Container | Volume           | Description                                       | Pack Size           | Packaging | Shelf Life (days) | Storage Temp (°C) |
|--------------------|-----------|------------------|---|---------------------|-----------|-------------------|-------------------|
| <b>Chromogenic</b> |           |                  |   |                     |           |                   |                   |
| PP3007             | PB90      | 10-12ml per side | CLED Agar / Primary UTI Chromogenic Agar (Opaque) | 10 x 90mm Bi-plates | Cellulose | 70                | 2-8°C             |
| PP3032             | P090      | 18-22ml          | Colorex™ Acinetobacter                            | 10 x 90mm Plates    | Cellulose | 42                | 2-8°C             |
| PP3052             | P090      | 18-22ml          | Colorex™ Acinetobacter MDR                        | 10 x 90mm Plates    | Cellulose | 42                | 2-8°C             |
| PP3076             | P090      | 18-22ml          | Colorex™ B.cereus                                 | 10 x 90mm Plates    | Cellulose | 56                | 2-8°C             |
| PP3087             | P090      | 18-22ml          | Colorex™ C. difficile                             | 10 x 90mm Plates    | Cellulose | 42                | 2-8°C             |
| PP3086             | P090      | 18-22ml          | Colorex™ Campylobacter                            | 10 x 90mm Plates    | Cellulose | 42                | 2-8°C             |
| PP3010             | P090      | 18-22ml          | Colorex™ Candida                                  | 10 x 90mm Plates    | Cellulose | 70                | 2-8°C             |
| PP3053             | P090      | 18-22ml          | Colorex™ C3GR                                     | 10 x 90mm Plates    | Cellulose | 42                | 2-8°C             |
| PP3030             | P090      | 18-22ml          | Colorex™ ECC                                      | 10 x 90mm Plates    | Cellulose | 42                | 2-8°C             |
| PP3001             | P090      | 18-22ml          | Colorex™ ESBL                                     | 10 x 90mm Plates    | Cellulose | 42                | 2-8°C             |
| PP7021             | P090      | 18-22ml          | Colorex™ Listeria (ISO)                           | 10 x 90mm Plates    | Cellulose | 56                | 2-8°C             |
| PP3077             | P090      | 18-22ml          | Colorex™ Malassezia                               | 10 x 90mm Plates    | Cellulose | 70                | 2-8°C             |
| PP3027             | PB90      | 10-12ml per side | Colorex™ Mastitis GP/Colorex™ Mastitis GN         | 10 x 90mm Bi-plates | Cellulose | 42                | 2-8°C             |
| PP3056             | P090      | 18-22ml          | Colorex™ MRSA                                     | 10 x 90mm Plates    | Cellulose | 42                | 2-8°C             |
| PP3059             | P090      | 15ml             | Colorex™ MRSA (15ml)                              | 10 x 90mm Plates    | Cellulose | 42                | 2-8°C             |
| PP3029             | PB90      | 10-12ml per side | Colorex™ MRSA/Colorex™ VRE                        | 10 x 90mm Bi-plates | Cellulose | 35                | 2-8°C             |
| PP3043             | PB90      | 10-12ml per side | Colorex™ MRSA/Colorex™ Staph aureus (Opaque)      | 10 x 90mm Bi-plates | Cellulose | 42                | 2-8°C             |
| PP3095             | P090      | 18-22ml          | Colorex™ mSuperCARBA™                             | 10 x 90mm Plates    | Cellulose | 35                | 2-8°C             |
| PP3028             | PB90      | 10-12ml per side | Colorex™ mSuperCARBA™ /Colorex™ C3GR              | 10 x 90mm Bi-plates | Cellulose | 35                | 2-8°C             |

| Product Code | Container  | Volume           | Description   | Pack Size                                    | Packaging | Shelf Life (days) | Storage Temp (°C) |
|--------------|------------|------------------|---|--|-----------|-------------------|-------------------|
| PP3040       | P090       | 18-22ml          | Colorex™ O157   | 10 x 90mm Plates                             | Cellulose | 35                | 2-8°C             |
| PP3041       | P090       | 18-22ml          | Colorex™ O157 with Cefixime & Tellurite   | 10 x 90mm Plates                             | Cellulose | 35                | 2-8°C             |
| PP3000       | P090       | 18-22ml          | Colorex™ Orientation Medium   | 10 x 90mm Plates                             | Cellulose | 42                | 2-8°C             |
| PP3034       | P090       | 18-22ml          | Colorex™ Pseudomonas  | 10 x 90mm Plates                             | Cellulose | 42                | 2-8°C             |
| PP1071       | P090       | 18-22ml          | Colorex™ Salmonella Plus  | 10 x 90mm Plates                             | Cellulose | 42                | 2-8°C             |
| PP2126       | PB90       | 10-12ml per side | Colorex™ Salmonella Plus /XLD Agar  | 10 x 90mm Bi-plates                          | Cellulose | 42                | 2-8°C             |
| PP3044       | P090       | 18-22ml          | Colorex™ Staph aureus   | 10 x 90mm Plates                             | Cellulose | 42                | 2-8°C             |
| PP3038       | P090       | 18-22ml          | Colorex™ STEC   | 10 x 90mm Plates                             | Cellulose | 42                | 2-8°C             |
| PP3065       | P090       | 18-22ml          | Colorex™ Strep B Agar   | 10 x 90mm Plates                             | Cellulose | 42                | 2-8°C             |
| PP3066       | PB90       | 10-12ml per side | Colorex™ Strep B Agar / Columbia Agar Base with 7% Defibrinated Horse Blood & CNA Selective Supplement        | 10 x 90mm Bi-plates                          | Cellulose | 42                | 2-8°C             |
| PP3035       | P090       | 18-22ml          | Colorex™ Vibrio   | 10 x 90mm Plates                             | Cellulose | 42                | 2-8°C             |
| PP3055       | P090       | 18-22ml          | Colorex™ VRE  | 10 x 90mm Plates                             | Cellulose | 35                | 2-8°C             |
| PP3078       | P090       | 18-22ml          | Colorex™ Y.enterocolitica   | 10 x 90mm Plates                             | Cellulose | 42                | 2-8°C             |
| PP2121       | P090       | 18-22ml          | Harlequin™ Salmonella ABC Medium  | 10 x 90mm Plates                             | Cellulose | 70                | 2-8°C             |
| PP6051       | P090       | 18-22ml          | Primary Chromogenic Coliform Agar (CCA)   | 10 x 90mm Plates                             | Cellulose | 42                | 2-8°C             |
| PP6028       | P090       | 18-22ml          | Primary membrane Lactose Glucuronide Agar (mLGA)  | 10 x 90mm Plates                             | Cellulose | 42                | 2-8°C             |
| PP1134       | P090       | 18-22ml          | Primary Tryptone Bile X - Glucuronide (TBX) Agar  | 10 x 90mm Plates                             | Cellulose | 56                | 2-8°C             |
| PP3005       | P090; PLSQ | 18-22ml; 55-57ml | Primary UTI Chromogenic Agar  | 10 x 90mm plates;<br>5 x 120mm Square Plates | Cellulose | 70                | 2-8°C             |
| PP3006       | P090       | 18-22ml          | Primary UTI Chromogenic Agar – Opaque   | 10 x 90mm Plates                             | Cellulose | 70                | 2-8°C             |
| PP3007       | PB90       | 10-12ml per side | CLED Agar / Primary UTI Chromogenic Agar (Opaque)   | 10 x 90mm Bi-plates                          | Cellulose | 70                | 2-8°C             |
| PP3008       | PB90       | 10-12ml per side | Primary UTI Chromogenic Agar / Columbia Agar Base with 7% Defibrinated Horse Blood & CNA Selective Supplement | 10 x 90mm Bi-plates                          | Cellulose | 42                | 2-8°C             |

| Product Code                                  | Container           | Volume                        | Description  | Pack Size  | Shelf Life (days) | Storage Temp (°C) |
|---|---------------------|-------------------------------|--|--|-------------------|-------------------|
| <b>Clinical, Veterinary &amp; Aquaculture</b> |                     |                               |  |  |                   |                   |
| BM0518  | M002                | 5ml                           | 2 x YT Broth with 0.5% Glycine                               | 50 x Thin Wall Universal   | 182               | 15-30°C           |
| BM0023  | M002                | 10ml                          | Alkaline Peptone Water – Double Strength                     | 50 x Thin Wall Universal   | 182               | 15-30°C           |
| BM0020  | M002; R1L           | 10ml, 20ml;<br>1000ml         | Alkaline Peptone Water (pH 8.6)                              | 50 x Thin Wall Universal; 6 x 1000ml Syrup Bottle                          | 182               | 15-30°C           |
| BM5810  | M003                | 4ml                           | Arginine Decarboxylase Medium                                | 50 x Thin Wall Bijou   | 365               | 15-30°C           |
| BM0040  | A125; M003          | 100ml; 3ml                    | Bile Aesculin Agar   | 10 x 125ml Syrup Bottle; 50 x Thin Wall Bijou                              | 365               | 2-8°C             |
| BM0076  | M037                | 15ml                          | Brain Heart Infusion Agar                                    | 50 x Polycarbonate Universal   | 365               | 15-30°C           |
| BM0070  | M002; M003;<br>R125 | 10ml, 9ml; 1ml,<br>3ml; 100ml | Brain Heart Infusion Broth                                   | 50 x Thin Wall Universal; 50 x Thin Wall Bijou;<br>10 x 125ml Syrup Bottle | 365               | 15-30°C           |
| BM0159  | M037                | 20ml                          | Brain Heart Infusion Broth with Ballotini Beads              | 50 x Polycarbonate Universal   | 365               | 15-30°C           |
| BM4295  | M002                | 10ml                          | Brain Heart Infusion Broth with 10% Defibrinated Horse Blood | 50 x Thin Wall Universal   | 90                | 2-8°C             |
| BM4181  | M003                | 2ml                           | Brain Heart Infusion Broth with 10% Glycerol                 | 50 x Thin Wall Bijou   | 182               | 15-30°C           |
| BM4290  | M002                | 10ml                          | Brain Heart Infusion Broth with 10% Sterile Horse Serum      | 50 x Thin Wall Universal   | 182               | 15-30°C           |
| BM0554  | M003; M010          | 3ml; 6ml                      | Brucella Broth   | 50 x Thin Wall Bijou; 50 x (145mm x 17mm) Tube                             | 365               | 15-30°C           |
| BM2022  | R300; R1L           | 225ml; 1000ml                 | Campylobacter Broth Bolton                                   | 10 x 300ml Syrup Bottle; 6 x 1000ml Syrup Bottle                           | 90                | 2-8°C             |
| BM0090  | M003                | 3ml                           | Columbia Agar  | 50 x Thin Wall Bijou   | 365               | 15-30°C           |
| BM0100  | M002; M003          | 10ml; 3ml                     | Columbia Agar Base with 7% Defibrinated Horse Blood Slope    | 50 x Thin Wall Universal; 50 x Thin Wall Bijou                             | 60                | 2-8°C             |

| Product Code | Container        | Volume                  | Description   | Pack Size   | Shelf Life (days) | Storage Temp (°C) |
|--------------|------------------|-------------------------|---|---|-------------------|-------------------|
| BM0560       | M003             | 3ml                     | Columbia Chocolate Slope  | 50 x Thin Wall Bijou  | 120               | 2-8°C             |
| BM0110       | M002             | 10ml, 15ml, 20ml        | Cooked Meat Medium with Brain Heart Infusion Broth Overlay      | 50 x Thin Wall Universal  | 365               | 15-30°C           |
| BM2250       | M002             | 10ml, 15ml, 20ml        | Cooked Meat Medium with Fastidious Anaerobe Broth Overlay       | 50 x Thin Wall Universal  | 365               | 15-30°C           |
| BM0452       | M003             | 3ml                     | CTA Urea Medium with Glucose                                    | 50 x Thin Wall Bijou  | 182               | 2-8°C             |
| BM5200       | M002             | 10ml                    | Czapek Dox Agar   | 50 x Thin Wall Universal  | 182               | 15-30°C           |
| BM0650       | As Requested     | N/A                     | Deionised Water   | N/A   | 730               | 15-30°C           |
| BM0670       | M002             | 10ml                    | Dermatophyte Test Medium with Chloramphenicol (100mg/L) (Slope) | 50 x Thin Wall Universal  | 182               | 2-8°C             |
| BM1351       | R1L              | 1000ml                  | Di-Potassium Hydrogen Phosphate Buffer (pH 7.5)                 | 6 x 1000ml Syrup Bottle   | 365               | 15-30°C           |
| BM0655       | M002; M003; R125 | 10ml; 1ml, 5ml; 100ml   | Distilled Water   | 50 x Thin Wall Universal; 50 x Thin Wall Bijou; 10 x 125ml Syrup Bottle | 730               | 15-30°C           |
| BM1328       | R125             | 100ml                   | Double Strength Acid Buffer                                     | 10 x 125ml Syrup Bottle   | 182               | 15-30°C           |
| BM0140       | N100             | 100ml                   | Egg Yolk Emulsion   | 1 x 100ml Square PETG Bottle  | 182               | 2-8°C             |
| BM0530       | N100             | 100ml                   | Egg Yolk Tellurite  | 1 x 100ml Square PETG Bottle  | 182               | 2-8°C             |
| BM0160       | M002; R125       | 10ml, 20ml, 25ml; 100ml | Fastidious Anaerobe Broth (FAB)                                 | 50 x Thin Wall Universal; 10 x 125ml Syrup Bottle                       | 365               | 15-30°C           |
| BM0164       | M037             | 5ml                     | Fastidious Anaerobe Broth (FAB) with Ballotini Beads            | 50 x Polycarbonate Universal  | 365               | 15-30°C           |
| BM0667       | M003             | 2ml                     | Formal Saline   | 50 x Thin Wall Bijou  | 365               | 15-30°C           |
| BM0666       | AM07             | 10ml                    | Formalin with Ballotini Beads                                   | 50 x Moulded Wall Universal   | 365               | 15-30°C           |



| Product Code | Container  | Volume          | Description                            | Pack Size   | Shelf Life (days) | Storage Temp (°C) |
|--------------|------------|-----------------|--|---|-------------------|-------------------|
| BM0453       | M003       | 3ml             | Group B Strep Broth                    | 50 x Thin Wall Bijou                                  | 90                | 2-8°C             |
| BM0710       | M003       | 5ml             | Iso-Sensitest Broth                    | 50 x Thin Wall Bijou                                  | 365               | 15-30°C           |
| BM2110       | M013       | 10ml, 6ml       | Kohn's Medium No.1                     | 50 x (150mm x 16mm) Tube                              | 365               | 2-8°C             |
| BM2120       | M013       | 10ml, 6ml       | Kohn's Medium No.2                     | 50 x (150mm x 16mm) Tube                              | 182               | 2-8°C             |
| BM1499       | A125       | 100ml           | Lactritmel Agar (Borelli's Medium)     | 10 x 125ml Syrup Bottle                               | 365               | 2-8°C             |
| BM1956       | M003; M010 | 5ml; 7ml        | LIM Broth with 10% Serum               | 50 x Thin Wall Bijou;<br>50 x (145mm x 17mm) Tube     | 90                | 2-8°C             |
| BM5280       | J500       | 300ml           | LB Agar (Miller)                       | 6 x 500ml Syrup Bottle                                | 182               | 15-30°C           |
| BM5300       | R500       | 350ml           | LB Broth (Miller)                      | 6 x 500ml Syrup Bottle                                | 182               | 15-30°C           |
| BM1489       | A300       | 200ml           | Malt Extract Agar                      | 10 x 300ml Syrup Bottle                               | 365               | 15-30°C           |
| BM0370       | M002; M003 | 10ml, 20ml; 5ml | Mannitol Selenite Broth                | 50 x Thin Wall Universal;<br>50 x Thin Wall Bijou     | 182               | 15-30°C           |
| BM1840       | M010; M014 | 6ml; 3ml        | Motility Indole Urea Medium            | 50 x (145mm x 17mm) Tube;<br>50 x (100mm x 16mm) Tube | 182               | 15-30°C           |
| BM2037       | P030       | 5ml             | MRSA Selective Broth                   | 50 x Polystyrene Universal                            | 28                | 2-8°C             |
| BM0800       | M003       | 3ml             | Mueller-Hinton Broth                   | 50 x Thin Wall Bijou                                  | 365               | 15-30°C           |
| BM0805       | M003       | 2ml             | Mueller-Hinton Fastidious Broth (MH-F) | 50 x Thin Wall Bijou                                  | 90                | 2-8°C             |
| BM0806       | M003       | 2ml             | Mueller-Hinton Broth (Cation Adjusted) | 50 x Thin Wall Bijou                                  | 365               | 15-30°C           |
| BM5500       | M003       | 3ml             | Non-Nutrient agar                      | 50 x Thin Wall Bijou                                  | 365               | 2-8°C             |

| Product Code | Container              | Volume                                      | Description                                 | Pack Size  | Shelf Life (days) | Storage Temp (°C) |
|--------------|------------------------|---|---|--|-------------------|-------------------|
| BM0540       | A125; A300; M002; M003 | 100ml; 250ml; 10ml; 3ml, 5ml                | Nutrient Agar                               | 10 x 125ml Syrup Bottle; 10 x 300ml Syrup Bottle; 50 x Thin Wall Universal; 50 x Thin Wall Bijou         | 365               | 15-30°C           |
| BM0550       | M002; M003; R125       | 10ml, 20ml; 3ml, 5ml; 100ml                 | Nutrient Broth                              | 50 x Thin Wall Universal; 50 x Thin Wall Bijou; 10 x 125ml Syrup Bottle                                  | 365               | 15-30°C           |
| BM0556       | M002                   | 10ml  | Nutrient Broth – Double Strength            | 50 x Thin Wall Universal   | 365               | 15-30°C           |
| BM1080       | M003                   | 3ml   | Nutrient Broth with 10% Sterile Horse Serum | 50 x Thin Wall Bijou   | 365               | 2-8°C             |
| BM4250       | M002                   | 10ml  | Nutrient Broth with 2.5% Sodium Chloride    | 50 x Thin Wall Universal   | 365               | 15-30°C           |
| BM0500       | M002                   | 10ml  | Nutrient Broth with 6% Sodium Chloride      | 50 x Thin Wall Universal   | 365               | 15-30°C           |
| BM0510       | M002                   | 10ml  | Nutrient Broth with 6.5% Sodium Chloride    | 50 x Thin Wall Universal   | 365               | 15-30°C           |
| BM5820       | M003                   | 4ml   | Ornithine Decarboxylase Medium              | 50 x Thin Wall Bijou   | 365               | 15-30°C           |
| BM4020       | M002                   | 20ml  | Pages Saline (Pages Amoeba Saline)          | 50 x Thin Wall Universal   | 365               | 15-30°C           |
| BM0301       | M002; R500             | 9ml; 450ml                                  | Peptone Water (0.1%)                        | 50 x Thin Wall Universal; 6 x 500ml Syrup Bottle   | 365               | 15-30°C           |
| BM0298       | M037                   | 5ml   | Peptone Water with Ballotini Beads          | 50 x Polycarbonate Universal   | 182               | 15-30°C           |
| BM2000       | M004                   | 3ml   | Peptone Water with 1% Glucose (Andrade's)   | 50 x Thin Wall Bijou with Durham Tube  | 182               | 15-30°C           |
| BM2031       | M013                   | 10ml  | Phenol Red Mannitol Agar                    | 50 x (150mmx16mm) Tube   | 182               | 2-8°C             |
| BM1358       | R1L                    | 1000ml                                      | Phosphate Buffer (0.01M)                    | 6 x 1000ml Syrup Bottle  | 365               | 15-30°C           |
| BM1359       | M002; M003; M037; P030 | 14ml, 16ml, 18ml; 2ml, 3ml; 20ml; 16ml, 2ml | Phosphate Buffer (0.067M)                   | 50 x Thin Wall Universal; 50 x Thin Wall Bijou; 50 x Polycarbonate Universal; 50 x Polystyrene Universal | 365               | 15-30°C           |
| BM1360       | M002; M003; R125; R300 | 10ml, 16ml, 20ml; 5ml; 100ml; 180ml         | Phosphate Buffered Saline                   | 50 x Thin Wall Universal; 50 x Thin Wall Bijou; 10 x 125ml Syrup Bottle; 10 x 300ml Syrup Bottle         | 365               | 15-30°C           |

| Product Code | Container        | Volume                  | Description   | Pack Size  | Shelf Life (days) | Storage Temp (°C) |
|--------------|------------------|-------------------------|---|--|-------------------|-------------------|
| BM0392       | M003             | 3ml                     | Phosphate Buffered Saline (PBS) with Ballotini Beads  | 50 x Thin Wall Bijou   | 730               | 15-30°C           |
| BM1361       | M002; R125       | 14ml; 100ml             | Potassium Dihydrogen Orthophosphate with Phenol Red   | 50 x Thin Wall Universal; 10 x 125ml Syrup Bottle                          | 182               | 15-30°C           |
| BM5230       | M002             | 10ml                    | Potato Dextrose Agar                                  | 50 x Thin Wall Universal   | 182               | 15-30°C           |
| BM0310       | M002; R500       | 5ml, 9ml; 500ml         | Ringers (¼ Strength)                                  | 50 x Thin Wall Universal; 6 x 500ml Syrup Bottle                           | 365               | 15-30°C           |
| BM0315       | M002; R125       | 10ml, 20ml; 100ml       | Ringers Solution (1:40 Strength)                      | 50 x Thin Wall Universal; 10 x 125ml Syrup Bottle                          | 365               | 15-30°C           |
| BM0313       | M037             | 5ml                     | Ringers with Ballotini Beads (¼ Strength)             | 50 x Polycarbonate Universal   | 365               | 15-30°C           |
| BM0340       | M002             | 10ml                    | Sabouraud Dextrose Agar                               | 50 x Thin Wall Universal   | 365               | 15-30°C           |
| BM0890       | A125; M002       | 100ml; 15ml             | Sabouraud Dextrose Agar with Chloramphenicol (50mg/L) | 10 x 125ml Syrup Bottle; 50 x Thin Wall Universal                          | 365               | 15-30°C           |
| BM0384       | M036             | 3ml                     | Saline (0.45%)  | 50 x (75mm x 12mm) Vitek Tube  | 365               | 15-30°C           |
| BM0380       | As Requested     | N/A                     | Saline (0.85%)  | N/A  | 365               | 15-30°C           |
| BM0385       | M003; M037       | 3ml; 1ml, 5ml           | Saline (0.85%) with Ballotini Beads                   | 50 x Thin Wall Bijou; 50 x Polycarbonate Universal                         | 365               | 15-30°C           |
| BM0381       | M003; M036       | 3ml; 1ml, 3ml           | Saline (0.9%)   | 50 x Thin Wall Bijou; 50 x (75mm x 12mm) Vitek Tube                        | 365               | 15-30°C           |
| BM0360       | M002; M013; R500 | 10ml, 15ml; 10ml; 450ml | Selenite F Broth                                      | 50 x Thin Wall Universal; 50 x (150mm x 16mm) Tube; 6 x 500ml Syrup Bottle | 182               | 15-30°C           |
| BM1940       | M011             | 145ml                   | Shadomy Agar  | 1 x 200ml Medical Flat Bottle  | 182               | 2-8°C             |
| BM0280       | M003; M010; M013 | 3ml; 6ml; 10ml          | Simmons Citrate Agar (Slope)                          | 50 x Thin Wall Bijou; 50 x (145mm x 17mm) Tube; 50 x (150mm x 16mm) Tube   | 182               | 15-30°C           |
| BM1950       | M003             | 5ml                     | Todd-Hewitt Broth                                     | 50 x Thin Wall Bijou   | 365               | 15-30°C           |

| Product Code | Container        | Volume                     | Description  | Pack Size   | Shelf Life (days) | Storage Temp (°C) |
|--------------|------------------|----------------------------|--|---|-------------------|-------------------|
| BM0412       | M003             | 5ml                        | Trichomonas Broth (CPLM)   | 50 x Thin Wall Bijou  | 182               | 2-8°C             |
| BM0413       | M003; M036; R500 | 3ml, 5ml; 2.5ml; 500ml     | Trichomonas Broth (CPLM) with Nystatin   | 50 x Thin Wall Bijou; 50 x (75mm x 12mm) Vitek Tube; 6 x 500ml Syrup Bottle | 182               | 2-8°C             |
| BM0440       | M013             | 6ml                        | Triple Sugar Iron Agar (Slope)   | 50 x (150mm x 16mm) Tube  | 182               | 2-8°C             |
| BM0436       | A500             | 400ml                      | Tryptic Soya Agar Base   | 6 x 500ml Syrup Bottle  | 365               | 15-30°C           |
| BM4701       | R125             | 100ml                      | Tryptone Soya Broth with 20% Glycerol  | 10 x 125ml Syrup Bottle   | 365               | 15-30°C           |
| BM2130       | M002; R1L; R300  | 10ml; 1000ml; 225ml        | Tryptone Soya Broth (Modified) with Novobiocin (20mg/L)                                      | 50 x Thin Wall Universal; 6 x 1000ml Syrup Bottle; 10 x 300ml Syrup Bottle  | 182               | 2-8°C             |
| BM0430       | M002; M003; R125 | 10ml, 20ml, 3ml; 3ml; 90ml | Tryptone Soya Broth  | 50 x Thin Wall Universal; 50 x Thin Wall Bijou; 10 x 125ml Syrup Bottle     | 365               | 15-30°C           |
| BM4704       | M003             | 2ml                        | Tryptone Soya Broth with 10% Glycerol  | 50 x Thin Wall Bijou  | 365               | 15-30°C           |
| BM0439       | J125; R1L        | 100ml; 1000ml              | Tryptone Soya Broth with 0.3% Lecithin & 2% Tween 80   | 10 x 125ml Syrup Bottle; 6 x 1000ml Syrup Bottle                            | 182               | 15-30°C           |
| BM0426       | R125             | 100ml                      | Tryptone Soya Broth with 0.07% Lecithin & 0.5% Tween 80                                      | 10 x 125ml Syrup Bottle   | 182               | 15-30°C           |
| BM0610       | M003             | 3ml, 5ml                   | Tryptone Water   | 50 x Thin Wall Bijou  | 365               | 15-30°C           |
| BM0450       | M003             | 3ml, 5ml                   | Urea Agar (Slope)  | 50 x Thin Wall Bijou  | 182               | 2-8°C             |
| BM0460       | M003             | 3ml                        | Urea Broth   | 50 x Thin Wall Bijou  | 182               | 2-8°C             |
| BM0400       | R125             | 100ml                      | Versene in PBS   | 10 x 125ml Syrup Bottle   | 182               | 2-8°C             |
| BM1671       | P003             | 2.5ml                      | Virus Transport Medium & Penicillin/Streptomycin Solution                                    | 50 x Plastic Bijou  | 120               | 2-8°C             |
| BM1672       | P003             | 4ml                        | Virus Transport Medium with Penicillin, Streptomycin, Polymyxin B, Amphotericin B & Nystatin | 50 x Plastic Bijou  | 120               | 2-8°C             |

| Product Code                           | Container        | Volume                     | Description  | Pack Size   | Shelf Life (days) | Storage Temp (°C) |
|--|------------------|----------------------------|--|---|-------------------|-------------------|
| <b>Food, Water &amp; Environmental</b> |                  |                            |  |   |                   |                   |
| BM0023                                 | M002             | 10ml                       | Alkaline Peptone Water – Double Strength                   | 50 x Thin Wall Universal  | 182               | 15-30°C           |
| BM0020                                 | M002; R1L        | 10ml, 20ml; 1000ml         | Alkaline Peptone Water (pH 8.6)                            | 50 x Thin Wall Universal; 6 x 1000ml Syrup Bottle                       | 182               | 15-30°C           |
| BM5810                                 | M003             | 4ml                        | Arginine Decarboxylase Medium                              | 50 x Thin Wall Bijou  | 365               | 15-30°C           |
| BM0040                                 | A125; M003       | 100ml; 3ml                 | Bile Aesculin Agar   | 10 x 125ml Syrup Bottle; 50 x Thin Wall Bijou                           | 365               | 2-8°C             |
| BM0070                                 | M002; M003; R125 | 10ml, 9ml; 1ml, 3ml; 100ml | Brain Heart Infusion Broth                                 | 50 x Thin Wall Universal; 50 x Thin Wall Bijou; 10 x 125ml Syrup Bottle | 365               | 15-30°C           |
| BM4181                                 | M003             | 2ml                        | Brain Heart Infusion Broth with 10% Glycerol               | 50 x Thin Wall Bijou  | 182               | 15-30°C           |
| BM0060                                 | M008             | 10ml                       | Brilliant Green Bile Broth (with Durham Tube)              | 50 x Thin Wall Universal with Durham Tube                               | 365               | 15-30°C           |
| BM0554                                 | M003; M010       | 3ml; 6ml                   | Brucella Broth   | 50 x Thin Wall Bijou; 50 x (145mm x 17mm) Tube                          | 365               | 15-30°C           |
| BM4015                                 | M002             | 10ml, 20ml                 | Buffered Nitrate Motility Medium                           | 50 x Thin Wall Universal  | 90                | 2-8°C             |
| BM0086                                 | M002             | 10ml                       | Buffered Peptone Water (Double Strength)                   | 50 x Thin Wall Universal  | 182               | 15-30°C           |
| BM0082                                 | As Requested     | N/A                        | Buffered Peptone Water (ISO)                               | N/A   | 365               | 15-30°C           |
| BM2022                                 | R300; R1L        | 225ml; 1000ml              | Campylobacter Broth Bolton                                 | 10 x 300ml Syrup Bottle; 6 x 1000ml Syrup Bottle                        | 90                | 2-8°C             |
| BM0090                                 | M003             | 3ml                        | Columbia Agar  | 50 x Thin Wall Bijou  | 365               | 15-30°C           |
| BM0110                                 | M002             | 10ml, 15ml, 20ml           | Cooked Meat Medium with Brain Heart Infusion Broth Overlay | 50 x Thin Wall Universal  | 365               | 15-30°C           |
| BM5400                                 | M002; R125       | 20ml, 9ml; 100ml           | D/E Neutralising Broth                                     | 50 x Thin Wall Universal; 10 x 125ml Syrup Bottle                       | 90                | 2-8°C             |

| Product Code | Container        | Volume                | Description   | Pack Size   | Shelf Life (days) | Storage Temp (°C) |
|--------------|------------------|-----------------------|---|---|-------------------|-------------------|
| BM0650       | As Requested     | N/A                   | Deionised Water   | N/A   | 730               | 15-30°C           |
| BM4901       | A300             | 250ml                 | de Man, Rogosa & Sharpe (MRS) Agar                        | 10 x 300ml Syrup Bottle   | 182               | 15-30°C           |
| BM4903       | R125             | 90ml                  | de Man, Rogosa & Sharpe (MRS) Broth                       | 10 x 125ml Syrup Bottle   | 182               | 15-30°C           |
| BM1351       | R1L              | 1000ml                | Di-Potassium Hydrogen Phosphate Buffer (pH 7.5)           | 6 x 1000ml Syrup Bottle   | 365               | 15-30°C           |
| BM0655       | M002; M003; R125 | 10ml; 1ml, 5ml; 100ml | Distilled Water   | 50 x Thin Wall Universal; 50 x Thin Wall Bijou; 10 x 125ml Syrup Bottle | 730               | 15-30°C           |
| BM1328       | R125             | 100ml                 | Double Strength Acid Buffer                               | 10 x 125ml Syrup Bottle   | 182               | 15-30°C           |
| BM0140       | N100             | 100ml                 | Egg Yolk Emulsion   | 1 x 100ml Square PETG Bottle  | 182               | 2-8°C             |
| BM0530       | N100             | 100ml                 | Egg Yolk Tellurite  | 1 x 100ml Square PETG Bottle  | 182               | 2-8°C             |
| BM0065       | M002             | 10ml                  | Enterobacteriaceae Enrichment Broth                       | 50 x Thin Wall Universal  | 182               | 15-30°C           |
| BM0640       | M002             | 10ml                  | Fraser Broth  | 50 x Thin Wall Universal  | 182               | 2-8°C             |
| BM0647       | R300; 5LB        | 225ml; 5000ml         | Fraser Broth (½ Strength)                                 | 10 x 300ml Syrup Bottle; 1 x 5L Media Bag                               | 182               | 2-8°C             |
| BM1070       | M002; M013       | 15ml; 10ml            | Glucose Agar with Bromocresol Purple Indicator            | 50 x Thin Wall Universal; 50 x (160mm x 16mm) Tube                      | 182               | 2-8°C             |
| BM1275       | A125             | 100ml                 | Iron Sulphite Agar  | 10 x 125ml Syrup Bottle   | 182               | 2-8°C             |
| BM1870       | M013             | 6ml                   | Kligler Iron Agar (Slope)                                 | 50 x (150mm x 16mm) Tube  | 182               | 15-30°C           |
| BM0775       | M002             | 10ml, 20ml            | Lactose-Gelatin Medium                                    | 50 x Thin Wall Universal  | 182               | 2-8°C             |
| BM0621       | M002             | 10ml                  | MacConkey Broth (Double Strength) with Cefotaxime (2mg/L) | 50 x Thin Wall Universal  | 90                | 2-8°C             |

| Product Code | Container                 | Volume                          | Description  | Pack Size   | Shelf Life (days) | Storage Temp (°C) |
|--------------|---------------------------|---------------------------------|--|---|-------------------|-------------------|
| BM0230       | R125                      | 100ml                           | MacConkey Broth (Purple)   | 10 x 125ml Syrup Bottle   | 365               | 15-30°C           |
| BM0620       | M002                      | 10ml                            | MacConkey Broth (Purple) (Double Strength)                                   | 50 x Thin Wall Universal  | 365               | 15-30°C           |
| BM1489       | A300                      | 200ml                           | Malt Extract Agar  | 10 x 300ml Syrup Bottle   | 365               | 15-30°C           |
| BM1490       | M003                      | 5ml                             | Malt Extract Broth   | 50 x Thin Wall Bijou  | 365               | 2-8°C             |
| BM0770       | R125                      | 100ml                           | Membrane Lauryl Sulphate Broth   | 10 x 125ml Syrup Bottle   | 182               | 2-8°C             |
| BM4000       | A125; A300                | 100ml; 200ml                    | Milk Plate Count Agar  | 10 x 125ml Syrup Bottle; 10 x 300ml Syrup Bottle  | 182               | 15-30°C           |
| BM2201       | M008                      | 10ml                            | Mineral Modified Glutamate Medium (Double Strength)                          | 50 x Thin Wall Universal with Durham Tube   | 182               | 2-8°C             |
| BM2200       | M008                      | 10ml                            | Mineral Modified Glutamate Medium (Single Strength)                          | 50 x Thin Wall Universal with Durham Tube   | 182               | 2-8°C             |
| BM4031       | A500                      | 450ml                           | Modified Semi-Solid Rappaport Vassiliadis (MSRV) Medium                      | 6 x 500ml Syrup Bottle  | 182               | 2-8°C             |
| BM1840       | M010; M014                | 6ml; 3ml                        | Motility Indole Urea Medium  | 50 x (145mm x 17mm) Tube;<br>50 x (100mm x 16mm) Tube   | 182               | 15-30°C           |
| BM0945       | M002                      | 10ml                            | Muller-Kauffmann Tetrathionate Novobiocin (MKTTn) Broth                      | 50 x Thin Wall Universal  | 90                | 2-8°C             |
| BM0946       | AM07                      | 10ml                            | Iodine & Potassium Iodine Complex  | 1 x Moulded Wall Universal  | 182               | 15-30°C           |
| BM0947       | M002                      | 10ml                            | Muller-Kauffmann Tetrathionate Novobiocin (MKTTn) with Iodine Complete Broth | 50 x Thin Wall Universal  | 60                | 2-8°C             |
| BM0540       | A125; A300;<br>M002; M003 | 100ml; 250ml;<br>10ml; 3ml, 5ml | Nutrient Agar  | 10 x 125ml Syrup Bottle; 10 x 300ml Syrup Bottle;<br>50 x Thin Wall Universal; 50 x Thin Wall Bijou | 365               | 15-30°C           |
| BM0550       | M002; M003;<br>R125       | 10ml, 20ml; 3ml,<br>5ml; 100ml  | Nutrient Broth   | 50 x Thin Wall Universal; 50 x Thin Wall Bijou;<br>10 x 125ml Syrup Bottle                          | 365               | 15-30°C           |
| BM0556       | M002                      | 10ml                            | Nutrient Broth – Double Strength   | 50 x Thin Wall Universal  | 365               | 15-30°C           |

| Product Code | Container              | Volume                              | Description  | Pack Size  | Shelf Life (days) | Storage Temp (°C) |
|--------------|------------------------|-------------------------------------|--|--|-------------------|-------------------|
| BM0510       | M002                   | 10ml                                | Nutrient Broth with 6.5% Sodium Chloride                     | 50 x Thin Wall Universal   | 365               | 15-30°C           |
| BM5610       | A300                   | 200ml                               | Orange Serum Agar  | 10 x 300ml Syrup Bottle  | 182               | 2-8°C             |
| BM5820       | M003                   | 4ml                                 | Ornithine Decarboxylase Medium                               | 50 x Thin Wall Bijou   | 365               | 15-30°C           |
| BM0300       | M003; M013             | 3ml; 5ml                            | Peptone Water  | 50 x Thin Wall Bijou; 50 x (150mm x 16mm) Tube   | 365               | 15-30°C           |
| BM2000       | M004                   | 3ml                                 | Peptone Water with 1% Glucose (Andrade's)                    | 50 x Thin Wall Bijou with Durham Tube  | 182               | 15-30°C           |
| BM0860       | M003; M004             | 3ml; 3ml                            | Peptone Water with Lactose (Andrade's)                       | 50 x Thin Wall Bijou; 50 x Thin Wall Bijou with Durham Tube                                      | 182               | 15-30°C           |
| BM2070       | M002; A300             | 15ml; 250ml                         | Perfringens Agar (Base Only)                                 | 50 x Thin Wall Universal; 10 x 300ml Syrup Bottle  | 365               | 15-30°C           |
| BM2031       | M013                   | 10ml                                | Phenol Red Mannitol Agar                                     | 50 x (150mmx16mm) Tube   | 182               | 2-8°C             |
| BM1358       | R1L                    | 1000ml                              | Phosphate Buffer (0.01M)                                     | 6 x 1000ml Syrup Bottle  | 365               | 15-30°C           |
| BM1360       | M002; M003; R125; R300 | 10ml, 16ml, 20ml; 5ml; 100ml; 180ml | Phosphate Buffered Saline                                    | 50 x Thin Wall Universal; 50 x Thin Wall Bijou; 10 x 125ml Syrup Bottle; 10 x 300ml Syrup Bottle | 365               | 15-30°C           |
| BM1363       | M002                   | 10ml                                | Phosphate Buffered Saline with 0.05% Tween 20                | 50 x Thin Wall Universal   | 365               | 2-8°C             |
| BM1986       | M002; P030             | 15ml; 25ml                          | Phosphate Buffered Water                                     | 50 x Thin Wall Universal; 50 x Polystyrene Universal   | 182               | 15-30°C           |
| BM0270       | A300; A500; M002       | 250ml; 400ml; 20ml                  | Plate Count Agar   | 10 x 300ml Syrup Bottle; 6 x 500ml Syrup Bottle; 50 x Thin Wall Universal                        | 365               | 15-30°C           |
| BM5230       | M002                   | 10ml                                | Potato Dextrose Agar   | 50 x Thin Wall Universal   | 182               | 15-30°C           |
| BM0274       | A500; M002             | 450ml; 20ml                         | R2A Agar   | 6 x 500ml Syrup Bottle; 50 x Thin Wall Universal   | 182               | 15-30°C           |
| BM0320       | M002; M010; P030       | 10ml; 6ml; 10ml                     | Rappaport Vassiliadis (RV) Single Component Enrichment Broth | 50 x Thin Wall Universal; 50 x (145mm x 17mm) Tube; 50 x Polystyrene Universal                   | 182               | 2-8°C             |
| BM0310       | M002; R500             | 5ml, 9ml; 500ml                     | Ringers (¼ Strength)   | 50 x Thin Wall Universal; 6 x 500ml Syrup Bottle   | 365               | 15-30°C           |
| BM0315       | M002; R125             | 10ml, 20ml; 100ml                   | Ringers Solution (1:40 Strength)                             | 50 x Thin Wall Universal; 10 x 125ml Syrup Bottle  | 365               | 15-30°C           |
| BM0380       | As Requested           | N/A                                 | Saline (0.85%)   | N/A  | 365               | 15-30°C           |



| Product Code | Container        | Volume                     | Description  | Pack Size  | Shelf Life (days) | Storage Temp (°C) |
|--------------|------------------|----------------------------|--|--|-------------------|-------------------|
| BM0381       | M003; M036       | 3ml; 1ml, 3ml              | Saline (0.9%)  | 50 x Thin Wall Bijou; 50 x (75mm x 12mm) Vitek Tube                        | 365               | 15-30°C           |
| BM0350       | M002             | 10ml                       | Selenite Cystine Broth                                   | 50 x Thin Wall Universal   | 182               | 15-30°C           |
| BM0360       | M002; M013; R500 | 10ml, 15ml; 10ml; 450ml    | Selenite F Broth   | 50 x Thin Wall Universal; 50 x (150mm x 16mm) Tube; 6 x 500ml Syrup Bottle | 182               | 15-30°C           |
| BM0280       | M003; M010; M013 | 3ml; 6ml; 10ml             | Simmons Citrate Agar (Slope)                             | 50 x Thin Wall Bijou; 50 x (145mm x 17mm) Tube; 50 x (150mm x 16mm) Tube   | 182               | 15-30°C           |
| BM5030       | A125             | 100ml                      | Slanetz & Bartley Agar                                   | 10 x 125ml Syrup Bottle  | 90                | 2-8°C             |
| BM0940       | M002             | 10ml                       | Tetrathionate Broth Base                                 | 50 x Thin Wall Universal   | 90                | 2-8°C             |
| BM0463       | A125             | 100ml                      | Tributyryn Agar  | 10 x 125ml Syrup Bottle  | 182               | 2-8°C             |
| BM0440       | M013             | 6ml                        | Triple Sugar Iron Agar (Slope)                           | 50 x (150mmx16mm) Tube   | 182               | 2-8°C             |
| BM0091       | M002             | 10ml                       | Tris Buffered Peptone Water                              | 50 x Thin Wall Universal   | 180               | 15-30°C           |
| BM0436       | A500             | 400ml                      | Tryptic Soya Agar Base                                   | 6 x 500ml Syrup Bottle   | 365               | 15-30°C           |
| BM0535       | A300             | 250ml                      | Tryptone Bile X-Glucuronide (TBX) Agar                   | 10 x 300ml Syrup Bottle  | 365               | 2-8°C             |
| BM1430       | A125; A300; M002 | 100ml; 250ml; 12ml         | Tryptone Soya Agar                                       | 10 x 125ml Syrup Bottle; 10 x 300ml Syrup Bottle; 50 x Thin Wall Universal | 365               | 15-30°C           |
| BM0430       | M002; M003; R125 | 10ml, 20ml, 3ml; 3ml; 90ml | Tryptone Soya Broth                                      | 50 x Thin Wall Universal; 50 x Thin Wall Bijou; 10 x 125ml Syrup Bottle    | 365               | 15-30°C           |
| BM4704       | M003             | 2ml                        | Tryptone Soya Broth with 10% Glycerol                    | 50 x Thin Wall Bijou   | 365               | 15-30°C           |
| BM4701       | R125             | 100ml                      | Tryptone Soya Broth with 20% Glycerol                    | 10 x 125ml Syrup Bottle  | 365               | 15-30°C           |
| BM2130       | M002; R1L; R300  | 10ml; 1000ml; 225ml        | Tryptone Soya Broth (Modified) with Novobiocin (20mg/L ) | 50 x Thin Wall Universal; 6 x 1000ml Syrup Bottle; 10 x 300ml Syrup Bottle | 182               | 2-8°C             |
| BM0610       | M003             | 3ml, 5ml                   | Tryptone Water   | 50 x Thin Wall Bijou   | 365               | 15-30°C           |
| BM0395       | M002             | 9ml                        | Tryptone Water with 0.85% Saline                         | 50 x Thin Wall Universal   | 365               | 15-30°C           |

| Product Code | Container        | Volume                   | Description                          | Pack Size  | Shelf Life (days) | Storage Temp (°C) |
|--------------|------------------|--------------------------|--------------------------------------|--|-------------------|-------------------|
| BM0450       | M003             | 3ml, 5ml                 | Urea Agar (Slope)                    | 50 x Thin Wall Bijou   | 182               | 2-8°C             |
| BM0460       | M003             | 3ml                      | Urea Broth                           | 50 x Thin Wall Bijou   | 182               | 2-8°C             |
| BM4060       | A125; M002       | 100ml; 15ml              | Violet Red Bile Agar (VRBGA)         | 10 x 125ml Syrup Bottle; 50 x Thin Wall Universal                          | 182               | 15-30°C           |
| BM0980       | A125; A300; M002 | 100ml; 250ml; 15ml       | Violet Red Bile Glucose Agar (VRGBA) | 10 x 125ml Syrup Bottle; 10 x 300ml Syrup Bottle; 50 x Thin Wall Universal | 182               | 15-30°C           |
| BM0470       | A500; M002       | 450ml, 500ml; 15ml, 20ml | Yeast Extract Agar                   | 6 x 500ml Syrup Bottle; 50 x Thin Wall Universal                           | 182               | 15-30°C           |


| Product Code                                     | Container        | Volume            | Description  | Pack Size  | Shelf Life (days) | Storage Temp (°C) | USP |
|--|------------------|-------------------|--|--|-------------------|-------------------|-----|
| <b>Pharmaceutical, Industrial &amp; Cosmetic</b> |                  |                   |  |  |                   |                   |     |
| BM5240   | R1L              | 1000ml            | Algal Broth  | 6 x 1000ml Syrup Bottle  | 365               | 15-30°C           |     |
| BM4710   | R300             | 200ml             | Antibiotic Medium No. 1 (USP)                                      | 10 x 300ml Syrup Bottle  | 365               | 15-30°C           | ●   |
| BM4711   | R300             | 200ml             | Antibiotic Medium No. 8 (USP)                                      | 10 x 300ml Syrup Bottle  | 365               | 15-30°C           | ●   |
| BM4712   | R125             | 100ml             | Antibiotic Medium No. 32 (USP)                                     | 10 x 125ml Syrup Bottle  | 365               | 15-30°C           | ●   |
| BM0092   | 5LB; R1L         | 5000ml; 1000ml    | Buffered Sodium Chloride Peptone Solution (USP)                    | 1 x 5 Litre Media Bag; 6 x 1000ml Syrup Bottle   | 365               | 15-30°C           | ●   |
| BM0089   | R125             | 100ml             | Buffered Sodium Chloride Peptone Solution (USP) with 0.1% Tween 80 | 10 x 125ml Syrup Bottle  | 365               | 15-30°C           | ●   |
| BM0093   | R125             | 100ml             | Buffered Sodium Chloride Peptone Solution (USP) with 1% Tween 80   | 10 x 125ml Syrup Bottle  | 365               | 15-30°C           | ●   |
| BM5110   | M002; R125       | 9ml; 50ml         | CEN Neutraliser  | 50 x Thin Wall Universal; 10 x 125ml Syrup Bottle  | 182               | 2-8°C             |     |
| BM5400   | M002; R125       | 20ml, 9ml; 100ml  | D/E Neutralising Broth   | 50 x Thin Wall Universal; 10 x 125ml Syrup Bottle  | 90                | 2-8°C             |     |
| BM0756   | R1L              | 1000ml            | Fluid A (USP)  | 6 x 1000ml Syrup Bottle  | 365               | 15-30°C           | ●   |
| BM0757   | 300l             | 200ml             | Fluid D (USP)  | 10 x 300ml Syrup Bottle  | 365               | 15-30°C           | ●   |
| BM0758   | R125             | 100ml             | Fluid K (USP)  | 10 x 125ml Syrup Bottle  | 365               | 15-30°C           | ●   |
| BM0954   | M002; M010; V100 | 10ml; 15ml; 100ml | Fluid Thioglycollate Medium (USP)                                  | 50 x Thin Wall Universal; 50 x (145mm x 17mm) Tube; 10 x 100ml Pharma Vial with Injectable Stopper | 365               | 15-30°C           | ●   |
| BM1075   | R125             | 100ml             | Lactose Broth  | 10 x 125ml Syrup Bottle  | 182               | 15-30°C           |     |

| Product Code | Container             | Volume                           | Description  | Pack Size   | Shelf Life (days) | Storage Temp (°C) | USP |
|--------------|-----------------------|----------------------------------|--|---|-------------------|-------------------|-----|
| BM0232       | R125                  | 100ml                            | MacConkey Broth (USP)  | 10 x 125ml Syrup Bottle   | 365               | 15-30°C           | ●   |
| BM4284       | R125                  | 90ml                             | Modified Lethen Broth  | 10 x 125ml Syrup Bottle   | 182               | 2-8°C             |     |
| BM0301       | M002; R500            | 9ml; 450ml                       | Peptone Water (0.1%)   | 50 x Thin Wall Universal;<br>6 x 500ml Syrup Bottle   | 365               | 15-30°C           |     |
| BM5230       | M002                  | 10ml                             | Potato Dextrose Agar (USP)   | 50 x Thin Wall Universal  | 182               | 15-30°C           | ●   |
| BM0274       | A500; M002            | 450ml; 20ml                      | R2A Agar   | 6 x 500ml Syrup Bottle;<br>50 x Thin Wall Universal   | 182               | 15-30°C           |     |
| BM0322       | M002                  | 10ml                             | Rappaport Vassiliadis Broth (USP)                                  | 50 x Thin Wall Universal  | 182               | 15-30°C           | ●   |
| BM0995       | A125                  | 100ml                            | Reinforced Clostridial Medium (USP)                                | 10 x 125ml Syrup Bottle   | 182               | 15-30°C           | ●   |
| BM0310       | M002; R500            | 5ml, 9ml; 500ml                  | Ringers (¼ Strength)   | 50 x Thin Wall Universal;<br>6 x 500ml Syrup Bottle   | 365               | 15-30°C           |     |
| BM0342       | A500                  | 500ml                            | Sabouraud Dextrose Agar (USP)                                      | 6 x 500ml Syrup Bottle  | 365               | 15-30°C           | ●   |
| BM0882       | J300                  | 200ml                            | Sabouraud Dextrose Agar with Lecithin (0.7g/L) & Tween 80 (5.0g/L) | 10 x 300ml Syrup Bottle   | 365               | 15-30°C           |     |
| BM0900       | M002                  | 10ml                             | Sabouraud Liquid Medium (USP)                                      | 50 x Thin Wall Universal  | 365               | 15-30°C           | ●   |
| BM0380       | As Requested          | N/A                              | Saline (0.85%)   | N/A   | 365               | 15-30°C           |     |
| BM0381       | M003; M036            | 3ml; 1ml, 3ml                    | Saline (0.9%)  | 50 x Thin Wall Bijou;<br>50 x (75mm x 12mm) Vitek Tube  | 365               | 15-30°C           |     |
| BM0940       | M002                  | 10ml                             | Tetrathionate Broth Base   | 50 x Thin Wall Universal  | 90                | 2-8°C             |     |
| BM0436       | A500                  | 400ml                            | Tryptic Soya Agar Base   | 6 x 500ml Syrup Bottle  | 365               | 15-30°C           |     |
| BM1435       | J300; J500            | 200ml; 400ml                     | Tryptone Soya Agar with Lecithin (0.7g/L) & Tween 80 (5.0g/L)      | 10 x 300ml Syrup Bottle;<br>6 x 500ml Syrup Bottle  | 365               | 15-30°C           |     |
| BM0432       | 5LB; M010; R125; V100 | 5000ml; 15ml; 100ml, 90ml; 100ml | Tryptone Soya Broth (USP)  | 1 x 5 Litre Media Bag;<br>50 x (145mm x 17mm) Tube;<br>10 x 125ml Syrup Bottle;<br>10 x 100ml Pharma Vial with Injectable Stopper | 365               | 15-30°C           | ●   |
| BM1432       | J300; J500            | 250ml; 500ml                     | Tryptone Soya Agar (USP)   | 10 x 300ml Syrup Bottle;<br>5 x 500ml Syrup Bottle  | 365               | 15-30°C           | ●   |
| BM4701       | R125                  | 100ml                            | Tryptone Soya Broth with 20% Glycerol                              | 10 x 125ml Syrup Bottle   | 365               | 15-30°C           |     |
| BM4704       | M003                  | 2ml                              | Tryptone Soya Broth with 10% Glycerol                              | 50 x Thin Wall Bijou  | 365               | 15-30°C           |     |
| BM0439       | J125; R1L             | 100ml; 1000ml                    | Tryptone Soya Broth with 0.3% Lecithin & 2% Tween 80               | 10 x 125ml Syrup Bottle;<br>6 x 1000ml Syrup Bottle   | 182               | 15-30°C           |     |
| BM0426       | R125                  | 100ml                            | Tryptone Soya Broth with 0.07% Lecithin & 0.5% Tween 80            | 10 x 125ml Syrup Bottle   | 182               | 15-30°C           |     |

| Product Code                                 | Container        | Volume        | Description  | Pack Size  | Shelf Life (days) | Storage Temp (°C) |
|--|------------------|---------------|--|--|-------------------|-------------------|
| <b>Media &amp; Reagents for Mycobacteria</b> |                  |               |  |  |                   |                   |
| BM8100                                       | M001; M037       | 4ml; 9ml      | Acidified Glycerol Egg Medium  | 50 x ½oz Bijou; 50 x Polycarbonate Universal                           | 270               | 15-30°C           |
| BM8110                                       | M001; M002; M037 | 4ml; 7ml; 9ml | Acidified Pyruvate Egg Medium  | 50 x ½oz Bijou; 50 x Thin Wall Universal; 50 x Polycarbonate Universal | 270               | 15-30°C           |
| BM0120                                       | M003             | 3ml, 5ml      | Dorset Egg Medium  | 50 x Thin Wall Bijou   | 365               | 2-8°C             |
| BM8400                                       | M001             | 4ml           | Kirchner Medium (Base only)  | 50 x ½oz Bijou   | 365               | 2-8°C             |
| BM8405                                       | M037             | 9ml           | Kirchner Medium with Antibiotics   | 50 x Polycarbonate Universal   | 365               | 2-8°C             |
| BM8410                                       | M001; M037       | 4ml; 9ml      | Kirchner Medium with Antibiotics and Calf Serum                                | 50 x ½oz Bijou; 50 x Polycarbonate Universal                           | 365               | 2-8°C             |
| BM0180                                       | M037             | 10ml          | Kirchner Medium with 2% Glycerol, 10% Sterile Horse Serum & PACT.              | 50 x Polycarbonate Universal   | 365               | 2-8°C             |
| BM8210                                       | M001             | 4ml           | Lowenstein-Jensen Glycerol Medium with Paranitrobenzoic Acid (PNB)             | 50 x ½oz Bijou   | 70                | 15-30°C           |
| BM8220                                       | M002             | 7ml           | Lowenstein-Jensen Glycerol Medium with Thiacetazone                            | 50 x Thin Wall Universal   | 70                | 15-30°C           |
| BM8310                                       | M001             | 4ml           | Lowenstein-Jensen Pyruvate Medium with Thiopenecarboxylic Acid Hydrazide (TCH) | 50 x ½oz Bijou   | 70                | 15-30°C           |
| BM0200                                       | M037             | 7ml           | Lowenstein-Jensen Medium   | 50 x Polycarbonate Universal   | 365               | 2-8°C             |
| BM0210                                       | M037             | 7ml           | Lowenstein-Jensen Medium with Pyruvate   | 50 x Polycarbonate Universal   | 365               | 2-8°C             |
| BM0216                                       | AM07             | 7ml           | Lowenstein-Jensen Pyruvate Medium with Penicillin (100K IUs)                   | 50 x Moulded Wall Universal  | 365               | 2-8°C             |
| BM0786                                       | N100             | 50ml          | Middlebrooks 7H9 Broth with OADC (40ml/L)                                      | 1 x 100ml Square PETG Bottle   | 182               | 2-8°C             |
| BM0780                                       | M037             | 7ml           | Middlebrooks 7H11 Agar with OADC Supplement (Slope)                            | 50 x Polycarbonate Universal   | 365               | 2-8°C             |
| BM0781                                       | AM07             | 7ml           | Middlebrooks 7H11 Selective Agar (No Glycerol)                                 | 50 x Moulded Wall Universal  | 365               | 2-8°C             |
| BM8200                                       | M001; M002; M037 | 4ml; 7ml; 9ml | Non-acidified Glycerol LJ Medium   | 50 x ½oz Bijou; 50 x Thin Wall Universal; 50 x Polycarbonate Universal | 270               | 15-30°C           |
| BM8300                                       | M001; M002; M037 | 4ml; 7ml; 9ml | Non-acidified Pyruvate LJ Medium   | 50 x ½oz Bijou; 50 x Thin Wall Universal; 50 x Polycarbonate Universal | 270               | 15-30°C           |

| Product Code          | Container                       | Volume  | Description                                       | Pack Size   | Shelf Life (days) | Storage Temp (°C) |
|-----------------------|---------------------------------|---|---|---|-------------------|-------------------|
| BM1301                | P030                            | 10ml, 2ml, 20ml                                   | Oxalic Acid (3%)                                  | 50 x Polystyrene Universal  | 182               | 15-30°C           |
| BM1300                | M003                            | 3ml   | Oxalic Acid (5%)                                  | 50 x Thin Wall Bijou  | 182               | 15-30°C           |
| BM1359                | M002;<br>M003;<br>M037;<br>P030 | 14ml, 16ml,<br>18ml; 2ml, 3ml;<br>20ml; 16ml, 2ml | Phosphate Buffer (0.067M)                         | 50 x Thin Wall Universal;<br>50 x Thin Wall Bijou;<br>50 x Polycarbonate Universal;<br>50 x Polystyrene Universal | 365               | 15-30°C           |
| BM1318                | P030                            | 11ml, 3ml, 7ml                                    | Sodium Hydroxide (2%) with Phenol Red             | 50 x Polystyrene Universal  | 182               | 15-30°C           |
| BM1322                | P030; R125                      | 10ml, 2ml, 4ml,<br>7ml; 100ml                     | Sodium Hydroxide (4%)                             | 50 x Polystyrene Universal;<br>10 x 125ml Syrup Bottle  | 182               | 15-30°C           |
| BM1320                | P030                            | 2ml, 7ml  | Sodium Hydroxide (4%) with Phenol Red             | 50 x Polystyrene Universal  | 182               | 15-30°C           |
| BM1325                | P030                            | 10ml  | Sputum Neutralising Buffer                        | 50 x Polystyrene Universal  | 182               | 15-30°C           |
| BM1324                | P030                            | 10ml, 3ml, 6ml                                    | Sputum Neutralising Buffer (No Phenol Red)        | 50 x Polystyrene Universal  | 182               | 15-30°C           |
| BM8350                | M002                            | 7ml   | Stonebrinks Medium                                | 50 x Thin Wall Universal  | 270               | 15-30°C           |
| BM1331                | P030                            | 7ml   | Sulphuric Acid (0.5N)                             | 50 x Polystyrene Universal  | 365               | 15-30°C           |
| BM1323                | P030                            | 4ml   | Sulphuric Acid (4%)                               | 50 x Polystyrene Universal  | 365               | 15-30°C           |
| BM8500                | M037                            | 10ml  | TB Diluent with Ballotini Beads                   | 50 x Polycarbonate Universal  | 182               | 15-30°C           |
| BM1314                | P030                            | 10ml, 14ml  | TB Buffer (pH 4.5)                                | 50 x Polystyrene Universal  | 182               | 15-30°C           |
| <b>Specimen Packs</b> |                                 |   |   |   |                   |                   |
| BM0379                | M037                            | 10ml  | Clean Area Pack (Ballotini Beads in 0.85% Saline) | 1 x Polycarbonate Universal per Pack with<br>50 Packs per Box   | 730               | 15-30°C           |
| BM0378                | S2                              | N/A   | Large Specimen Pack                               | 1 x 2.5 Litre Tub per Pack with 5 Packs per Box   | 730               | 15-30°C           |

Bottles – Tubes – Vials



| Product Code | Description   | Product Code   | Description  |
|--------------|---|----------------|--|
| 10LB         | 10 Litre Blood Bag  | P003           | Polypropylene Bijou with Polypropylene Cap                         |
| 1LB          | 1 Litre Blood Bag   | P030           | 30ml Universal Polystyrene Bottle with Polypropylene Cap           |
| 4LB          | 4 Litre Blood Bag   | V001           | Freeze Dry Vial 5ml  |
| 5LB          | 5 Litre Media Bag   | V002           | Freeze Dry Vial 20ml   |
| H250         | Honey Jar with Metal Cap 250ml                                      | V100           | Pharma Vial (DIN) 100ml with Injectable Stopper and Crimp Closure  |
| M001         | ½oz Moulded Bijou with Metal Cap                                    | V1L            | Pharma Vial (DIN) 1000ml with Injectable Stopper and Crimp Closure |
| M002         | Thin Wall Universal 30ml with 24mm Polypropylene Cap                | V500           | Pharma Vial (DIN) 500ml with Injectable Stopper and Crimp Closure  |
| M003         | Thin Wall Bijou 7ml with 18mm Polypropylene Cap                     | 060            | Alpha (Syrup) Bottle 60ml with Cap Option                          |
| M004         | Thin Wall Bijou with Durham Tube                                    | 125            | Alpha (Syrup) Bottle 100ml with Cap Option                         |
| M008         | Thin Wall Universal with Durham Tube                                | 300            | Alpha (Syrup) Bottle 300ml with Cap Option                         |
| M010         | 145 x 17mm Soda Glass Screw Neck Tube with Gold Cap                 | 500            | Alpha (Syrup) Bottle 500ml with Cap Option                         |
| M011         | Medical Flat 200ml with 28mm Polypropylene Cap                      | 1L             | Alpha (Syrup) Bottle 1000ml with Cap Option                        |
| M013         | 150 x 16mm Soda Glass Test Tube with Red Push Cap                   | Cap Option - I | 28mm Injectable Cap  |
| M014         | 100 x 16mm Borosilicate Screw Neck Tube 15ml with Polypropylene Cap | Cap Option - Y | 28mm Vima Cap  |
| M015         | ¼oz Moulded Bijou with Metal Cap                                    | Cap Option - A | Polypropylene Cap  |
| M036         | 75 x 12mm Vitek Polystyrene Push Cap 5ml                            | Cap Option - J | 28mm Pharma Tamper Evident Cap                                     |
| M037         | Polycarbonate Universal Bottle 30ml                                 | Cap Option - R | 28mm Aluminium Tamper Evident Cap                                  |
| AM07         | Moulded Universal Vial 30ml   | 5000           | Jerry Can 5Lt  |
| N100         | Square graduated PETG Sterile Bottle with Tamper Evident Cap 100ml  | no code        | Durham Tube  |
| N500         | Square graduated PETG Sterile Bottle 500ml with Tamper Evident Cap  | no code        | Aluminium foil for injectable Cap                                  |
| N1L          | Square graduated PETG Sterile Bottle 1000ml with Tamper Evident Cap |                |  |

**Bottles – Tubes – Vials  
Media/Blood Bags**

Aluminium foil for Injectible Cap



Injectable Tamper evident I 28mm



Vima Cap Y 28mm



Polypropylene Cap A 28mm



Pharma Tamper evident J 28mm



Aluminium Tamper evident R 28mm



28mm



1L 1L



500 500ml



300 300ml



125 125ml



060 60ml

Square PETG Bottles (Sterile)



N100



N500



Honey Jar

H250 250ml



Medical Flat

M011 200ml

**Product Code Petri-Dish Description**

|      |   |
|------|---|
| P090 | Triple Vented 90mm Crystal Polystyrene  |
| P055 | Six Vented 55mm Crystal Polystyrene     |
| P100 | Triple Vented 100mm Crystal Polystyrene |
| P140 | Triple Vented 140mm Crystal Polystyrene |

**Product Code Petri-Dish Description**

|      |  |
|------|--|
| PB90 | Triple Vented 90mm Two Compartment Crystal Polystyrene   |
| PT90 | Triple Vented 90mm Three Compartment Crystal Polystyrene |
| PLSQ | Four Vented 120mm Square Crystal Polystyrene             |
| PCON | Four Vented 55mm Grip lid Crystal Polystyrene Contact    |



Media Bag 5Ltr 5LB



Blood Bags 1Ltr 1LB 4Ltr 4LB 10Ltr 10LB




Head Office



Donor Animal Farm



cGMP Manufacturing



Fully Automated Lines

**E&O Laboratories Ltd**  
Burnhouse, Bonnybridge  
FK4 2HH Scotland, UK

Tel: **+44 (0)1324 840404**  
Fax: **+44 (0)1324 841314**  
E-mail: **[info@eolabs.com](mailto:info@eolabs.com)**  
Web: **[www.eolabs.com](http://www.eolabs.com)**