




# LABCOLD™

medical & scientific refrigeration



## RLHE0845 -40°C Sparkfree Chest Freezer

**Warning: Do not store open containers of volatile substances in this freezer.**

Compact chest freezer featuring a spark free interior for laboratory safety.

Fitted with a lid lock as standard, this freezer also comes complete with a digital temperature display and high and low temperature alarms.



Bespoke racking can be provided if required as can baskets at an extra cost. To increase the life of your freezer and to ensure you get optimum performance, preventative maintenance contracts are also available. Why not ask for details when you order?



Digital Temperature Display

Part Number	RLHE0845
Capacity (Litres)	237
Style	Chest Freezer
Exterior Dimensions (HxWxD)	860 x 1050 x 730mm
Interior Dimensions (HxWxD)	640 x 845 x 445mm
Weight	120kg
Operating Temperature	-20°C to -45°C
Interior Construction	Aluminium
Exterior Construction	Painted Steel
Defrost	Manual
Sparkfree Interior	✓
High Temperature Alarm	✓
Low Temperature Alarm	✓
Digital Temperature Display	✓
Door Lock	✓
Feet	4
Energy Consumption	0.61 kWh/24h
Eco Friendly Refrigerant	✓

01256 705 570

[www.labcold.com](http://www.labcold.com)

