

Stronghold (GM006*)

AQL 0.65 · chemical tested · chemotherapy drug tested · purple medical grade · heavy weight · nitrile · examination gloves

Physical Specifications:

Material: Nitrile (latex free)
Colour: Purple
Length: 240 mm
Weight: 6.2 g
Thickness: 0.12 mm
Finish: Powder free, chlorinated
Fit: Ambidextrous
Sterile: No

Order Codes and Size Information:

GM0061 | Extra small
GM0062 | Small
GM0063 | Medium
GM0064 | Large
GM0065 | Extra large

Packaging:

Each dispenser box contains 100 gloves
Each case contains 10 dispenser boxes

Classification:



AQL 0.65 (Pinholes) EN455 (Medical Grade)



For more about the **Advanced** range, please **click on the image**



Features


- Tested to EN374
- Tested to ASTM D6978-05
- AQL 0.65
- Tested to EN1186
- Nitrile formulation conforms to the hand during use
- Beaded cuff
- Manufactured with textured palms and fingertips
- Heavy-weight construction
- Touchscreen compatible

Benefits

- Resistant to a range of chemicals
- Resistant to a range of cytotoxic drugs when tested to ASTM D6978-05
- Ultra-low hole level; exceeds medical standards for pinholes
- Suitable for contact with all food types
- Nitrile formulation conforms to the hand during use for a more comfortable fit
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Textured palms and fingertips provide improved grip, particularly when working in wet or oily conditions
- Heavy-weight construction is tough, durable and long lasting
- Suitable for use with mobile and tablet displays without removing gloves

Sectors

- Aerospace
- Chemical manufacture and use
- Dentistry • Dental hygiene
- Electronics • Engineering
- Food processing
- Healthcare • Laboratory work
- Manufacturing and assembly
- Pharmaceutical work • Printing



UNIGLOVES UK
Registered in England No: 4010200
VAT registration 791817787

3 Ambley Green, Gillingham Business Park,
Gillingham, ME8 0NJ, UK
+44 (0)800 049 6602
enquiries@unigloves.co.uk
unigloves.co.uk

UNIGLOVES GERMANY
Arzt- & Klinikbedarf
Handelsgesellschaft mbH

Camp-Spich-Straße 71
D-53842 Troisdorf-Spich,
+49 (0) 2241-9323-0
unigloves@unigloves.de
unigloves.de

