

Zero (GM005*)

zinc, sulphur and accelerator free · low modulus · anti-static blue · medical grade · medium weight · nitrile · examination gloves

Physical Specifications:

Material: Zinc free, sulphur free, and accelerator free nitrile (latex free)
 Colour: Blue
 Length: 240 mm
 Weight: 3.8 g
 Thickness: 0.08 mm
 Finish: Powder free, chlorinated
 Fit: Ambidextrous
 Sterile: No

Order Codes and Size Information:

GM0051 | Extra small
 GM0052 | Small
 GM0053 | Medium
 GM0054 | Large
 GM0055 | Extra large

Packaging:


Each dispenser box contains 100 gloves
 Each case contains 10 dispenser boxes

Classification:

AQL 1.5 (Pinholes) EN455 (Medical Grade)

Food Safe EN ISO 374-1 EN ISO 374-5: VIRUS

For more about the **Advanced** range, please **click on the image**



Features


- Free from zinc, sulphur and accelerators
- Tested to EN374
- Low modulus formulation
- Anti-static / EN1149 compliant
- AQL 1.5
- Free from powder
- Free from phthalates
- Tested to EN1186
- Nitrile formulation conforms to the hand during use
- Manufactured with textured fingertips
- Medium-weight construction
- Touchscreen compatible

Benefits

- ZAS free formulation reduces risk of Type I and Type IV allergic reactions
- Resistant to a range of chemicals
- Low modulus formulation closely mimics the elasticity and comfort of a latex glove which minimises hand fatigue but retains the other benefits of nitrile
- Anti-static formulation suppresses static charge, preventing sparks
- Meets medical standards for pinholes
- Suitable for those with latex allergies (Type I sensitisation)
- Reduces risks of latex allergies (Type I sensitisation) and minimises particulate contamination
- Suitable for contact with all food types
- Nitrile formulation conforms to the hand during use for a more comfortable fit
- Textured fingertips provide improved grip, particularly when working in wet or oily conditions
- Medium-weight construction provides a balance of strength and sensitivity
- Suitable for use with mobile and tablet displays without removing gloves

Sectors

- Aerospace • Chemical manufacture and use • Dentistry
- Dental hygiene • Electronics • Engineering • Food processing • Healthcare
- Laboratory work • Manufacturing and assembly • Pharmaceutical work • Printing



UNIGLOVES UK
 Registered in England No: 4010200
 VAT registration 791817787

📍 3 Ambley Green, Gillingham Business Park, Gillingham, ME8 0NJ, UK
 ☎ +44 (0)800 049 6602
 ✉ enquiries@unigloves.co.uk
 🌐 unigloves.co.uk

UNIGLOVES GERMANY
 Arzt- & Klinikbedarf
 Handelsgesellschaft mbH

📍 Camp-Spich-Straße 71
 D-53842 Troisdorf-Spich,
 ☎ +49 (0) 2241-9323-0
 ✉ unigloves@unigloves.de
 🌐 unigloves.de



ughealthcarecorporation.com