

# Pearl (GPO\*\*\*)

15 colours · medical grade · medium weight nitrile · examination gloves

**Physical Specifications:**

Material: Nitrile (latex free)  
 Colours: Aqua Blue, Black, Blue, Burgundy, Cobalt, Green, Lime, Magenta, Peach, Pink, Red, Sapphire, Violet, White, Yellow  
 Length: 240 mm  
 Weight: 3.8 g  
 Thickness: 0.08 mm  
 Finish: Powder free, chlorinated  
 Fit: Ambidextrous  
 Sterile: No

**Order Codes and Size Information:**

GPO\*\*1 | Extra small  
 GPO\*\*2 | Small  
 GPO\*\*3 | Medium  
 GPO\*\*4 | Large  
 GPO\*\*5 | Extra large

\*\* 01 Blue / 02 White / 03 Black / 04 Green / 05 Pink / 06 Red / 07 Violet / 08 Sapphire / 09 Aqua Blue / 10 Peach / 11 Yellow / 12 Burgundy / 13 Cobalt / 14 Lime / 15 Magenta

**Note:** Not all colours come in Extra Small

**Packaging:**

Each dispenser box contains 100 gloves  
 Each case contains 10 dispenser boxes

**Classification:**



For more about the Pearl range, please click on the image



**Features**

- Tested to EN374
- Available in fifteen colours
- AQL 1.5 and EN455 parts 1-4
- Free from powder
- Tested to EN1186
- Nitrile formulation conforms to the hand during use
- Beaded cuff
- Manufactured with textured palms and fingertips
- Touchscreen compatible

**Benefits**

- Resistant to a range of chemicals
- Ideal for colour-coding in food, cleaning and healthcare environments
- Meets medical standards including pinholes
- Reduces risks of latex allergies (Type I sensitisation) and minimises particulate contamination
- Suitable for contact with all food types
- Nitrile formulation conforms to the hand during use for a more comfortable fit
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Textured palms and fingertips provide improved grip, particularly when working in wet or oily conditions
- Suitable for use with mobile and tablet displays without removing gloves

**Sectors**

- Care home work
- Chemical manufacture and use
- Cleaning and janitorial
- Dentistry • Dental hygiene • Emergency services
- Floristry and gardening • Food processing
- Hair and beauty • Healthcare
- Hospital and clinical duties
- Janitorial • Painting and decorating
- Piercing • Tattoo

**Colour range:**



**UNIGLOVES UK**  
 Registered in England No: 4010200  
 VAT registration 791817787

3 Ambley Green, Gillingham Business Park,  
 Gillingham, ME8 0NJ, UK  
 +44 (0)800 049 6602  
 enquiries@unigloves.co.uk  
 unigloves.co.uk

**UNIGLOVES GERMANY**  
 Arzt- & Klinikbedarf  
 Handelsgesellschaft mbH

Camp-Spich-Straße 71  
 D-53842 Troisdorf-Spich,  
 +49 (0) 2241-9323-0  
 unigloves@unigloves.de  
 unigloves.de

