

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com Email USA: techsery@sial.com Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Potassium iodide - ReagentPlus*, 99%

Product Number: 207969 CAS Number: 7681-11-0 MFCD00011405 MDL:

Formula: ΙK

Formula Weight: 166.00 g/mol

TEST Specification

Appearance (Color) Appearance (Form) X-Ray Diffraction Titration with AgNO3 % lodine Recommended Retest Period

42 Months

Specification: PRD.3.ZQ5.10000016864

White

Powder or Crystals Conforms to Structure

75.3 - 77.6 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.