Corning Inc Life Sci ISO 900	ences ered Product Description
Catalog Number:	3477
Product Description:	Costar ® Netwell <sup>TM</sup> , 15mm, 74 $\mu$ m mesh size, in a 12 well cluster plate with lid
Component Materials: Plate/Lid Netwell body Mesh bottom	<ul> <li>Virgin Polystyrene, meets USP, Class VI requirements for plastic containers and closures.</li> <li>Virgin Polystyrene, meets USP, Class VI requirements for plastic containers and closures with grey concentrate.</li> <li>Polyester</li> </ul>
<b>Product Dimensions:</b> Length of Plate Width of Plate Tolerances of Dimensions	<ul> <li>5.035 in.</li> <li>3.370 in.</li> <li>+/050 in.</li> <li>Bottom of Netwell to plate060 in.</li> <li>Height with Lid890 in.</li> </ul>

## Sterilization:

The product is irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137 Sterilization of healthcare *products-Requirements for validation and routine control-Radiation sterilization.* Sterility Assurance Level: SAL 10<sup>-3</sup>

**Performance Testing:** Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures. Visual examination of the product. Inspection for seal and barrier integrity, accurate labeling, and correct Visual Attributes: Packaging: product configuration.

## Lot Number Designation:

8 Digit Lot Number: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Rev No: 4