Corning Inc Life Sci ISO 900	ered Product Description
Catalog Number:	3479
Product Description:	Costar ${f B}$ Netwell ${}^{T\!M}$, 24mm, 74 μ m mesh size, in a 6 well cluster plate with lid
Component Materials: Plate/Lid Netwell body Mesh bottom	 Virgin Polystyrene, meets USP, Class VI requirements for plastic containers and closures. Virgin Polystyrene, meets USP, Class VI requirements for plastic containers and closures with grey concentrate. Polyester
Product Dimensions: Length of Plate Width of Plate Tolerances of Dimensions	 5.030 in. 3.365 in. +/050 in. Bottom of Netwell to plate050 in. Height with Lid - 1.100 in.

Sterilization:

The product is irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137 Sterilization of healthcare *products-Requirements for validation and routine control-Radiation sterilization.* Sterility Assurance Level: SAL 10⁻³

Performance Testing: Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures. Visual examination of the product. Inspection for seal and barrier integrity, accurate labeling, and correct Visual Attributes: Packaging: product configuration.

Lot Number Designation:

8 Digit Lot Number: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Rev No: 6